STS-IFU-001-REV1

strauss[®]
S U R G I C A L

Instructions For Use

Cystoscopes/Hysteroscopes
(USA)

TPA766-025-01_C
Version: C
2019 - Dec - 16

	Content	Page
0	Important Information About this Document	4
0.1	Scope of Validity, Identification, Purpose.....	4
0.2	Target Group.....	4
0.3	Using and Storing this Document.....	4
1	General Information About the Device.....	5
1.1	Scope of Delivery	5
1.2	Device Description.....	5
1.2.1	Visual Overview.....	5
1.3	Usage.....	6
1.3.1	Intended Use.....	6
1.3.2	Indications for Use.....	6
1.3.3	Contraindications	7
1.4	Marking	9
1.4.1	Pictograms and Information on the Device and Packaging	9
1.4.2	Pictograms in this Document.....	10
1.5	Service Department Contact Details	10
2	General Safety Information	11
2.1	Safety Messages in this Document.....	11
2.1.1	Safety Messages at the Start of a Chapter.....	11
2.1.2	Safety Messages in the Body of the Text.....	11
2.2	Product Safety	12
2.2.1	Basic Safety Information.....	12
2.2.2	Staff Qualifications.....	12
3	Operation.	13
3.1	Safety Instructions.....	13
3.2	Staff Qualifications.....	14
3.3	Technical Inspection Prior to Use.....	15
3.3.1	Visual Inspection.....	15
3.3.2	Function Check	16
3.4	Connecting a Light Guide.....	16
3.5	Post-Use and Pre-Cleaning Activities	16
4	Reprocessing	17
4.1	Safety Instructions.....	17
4.2	Staff Qualifications.....	18
4.3	Validated Procedures.....	18
4.4	The Process Flow.	19
4.5	Cleaning and Disinfecting Agents.....	19
4.6	Manual Cleaning in Soaking Bath	20
4.6.1	Overview.....	20
4.6.2	Performance.....	20
4.7	Automated Cleaning.....	22
4.7.1	Overview	22
4.7.2	Performance.....	22
4.8	High Level Disinfection (HLD).....	23
4.8.1	Overview	23
4.8.2	Performance.....	23
4.9	Sterilization	24
4.9.1	Steam Sterilization.....	24
5	Maintenance and Repairs.....	25
5.1	Troubleshooting Table.	25
5.2	Repair.....	26
6	Product Data.....	27
6.1	Technical Data.....	27
6.2	Ambient Conditions	29
6.3	Spare Parts and Accessories.....	29
7	Disposal.....	30

0 Important Information About this Document

0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the following rigid endoscopes:

STS-2040-302ASA, STS-2040-302AA,
STS-2040-302BA, STS-2040-302BSA,
STS-2040-302ESA, STS-2040-302CA,
STS-2029-302AA, STS-2029-302BA,
STS-2040-298DSA-RWF, STS-2040-298ASA-RWF,
STS-2040-300GSA-RWF, STS-2040-303CA-RWF,
STS-2040-280ASA-OLY, STS-2040-280ESA-OLY,
STS-2040-283BSA-OLY, STS-2040-283BA-OLY,
STS-2040-283CA-OLY

These instructions for use are an integral component of the device and contain all the information required by users and operators for safe and proper use.

These instructions for use do not describe endoscopic procedures or techniques during surgery.

0.2 Target Group

These instructions for use are intended for physicians, medical assistants, and employees of central sterile supply who are entrusted with the operation, handling and reprocessing of the device.

0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

1 General Information About the Device

1.1 Scope of Delivery

Check scope of delivery! Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

1.2 Device Description

1.2.1 Visual Overview

Cystoscope/Hysteroscope


Figure 1-1: Cystoscope/Hysteroscope.

- 1 Eyepiece cup
- 2 Light guide connector (ACMI)
- 3 Adapter for connecting Storz light guides
- 4 Adapter for connecting Wolf light guides
- 5 Endoscope sheath
- 6 Endoscope tip

1.3 Usage

1.3.1 Intended Use

The device is intended to be used in combination with medical electrical devices that are, as a **minimum**, compliant with the requirements for applied parts classified as type **BF** (that is parts offering increased protection against electric shock), according to IEC 60601-1.

This insulation barrier can be realized by the medical electrical devices themselves or by the connecting cables to the endoscope and must be assured for every connection between the endoscope and connected devices.

Devices or connecting cables that comply with the requirements for applied parts classified as type BF feature the adjacent symbol.

Devices and connecting cables that comply with the requirements for applied parts classified as type CF offer an even higher level of protection. Devices bearing the CF symbol (the adjacent heart symbol) are approved for direct conductive contact with the heart.

Defibrillation-proof devices or connecting cables are also compatible. These will be marked with one of the adjacent symbols depending on the degree of protection against electric shock offered.

1.3.2 Indications for Use

The Cystoscopes/Hysteroscopes are indicated to provide the user with the means for endoscopic diagnostic and therapeutic surgical procedures. Examples for the use of the devices include the visualization and manipulation of anatomy as the surgeon deems appropriate. The Cystoscopes/Hysteroscopes are intended to be used in general urological and gynecological surgery through a minimally invasive approach by utilizing natural orifices to access the surgical site.

Note: For this medical application, you should comply with national statutory regulations and guidelines.

1.3.3 Contraindications

General contraindications

- > Contraindications may exist depending on the patient's general condition or the specific presenting disease.
- > The decision to perform a procedure rests with the responsible surgeon and must be made on the basis of an individual benefit/risk analysis.
- > Acute pelvic inflammatory disease

Hysteroscopy may be contraindicated by the following conditions, depending on their severity or extent:

- > Inability to distend the uterus
- > Cervical stenosis
- > Cervical/vaginal infection
- > Uterine bleeding or menses
- > Known pregnancy
- > Invasive carcinoma of the cervix
- > Recent uterine perforation
- > Medical contraindication or intolerance to anesthesia

Contraindication to Endometrial Ablation

Hysteroscopic endometrial ablation, whether by laser or electrosurgery, should not be undertaken without adequate training, preceptorship, and clinical experience. Additionally, endometrial biopsy should be performed prior to any ablation. The following are clinical conditions that can significantly complicate hysteroscopic endometrial ablation:

- > Adenomatous Endometrial Hyperplasia
- > Uterine Leiomyoma
- > Severe Adenomyosis
- > Pelvic Pain (Subtle PID)
- > Uterine anomalies

Contraindication to Hysteroscopic Myomectomy

Hysteroscopic myomectomy should not be undertaken without adequate training, preceptorship, and clinical experience. The following are clinical conditions that can significantly complicate hysteroscopic myomectomy:

- > Severe anemia
- > Inability to circumnavigate a myoma due to myoma size (e.g. predominantly intramural myomas with small submucous components).

⚠ WARNING

Continuous Flow Hysteroscopy

Risk to the patient

If a liquid distension medium is used, strict fluid intake and output surveillance should be maintained to ensure that fluid deficit is known at all times. Depending on whether non-electrolyte or electrolyte solution is being used, when excessive fluid deficit occurs, consideration should be given to stopping further infusion and concluding the procedure.

Potential Complications of Continuous Flow Hysteroscopy:

- > Hyponatremia
- > Hypothermia
- > Uterine perforation resulting in potential injury to bowel, bladder, major blood vessels, and ureter
- > Pulmonary edema
- > Cerebral edema
- > Death

⚠ WARNING

Additional warnings in conjunction with hysteroscopic or cystoscopic application

Risk to the patient

- > If CO₂ gas is used as a distention medium, operative hysteroscopy is contraindicated due to the risk of gas embolization. CO₂ gas may be used for diagnostic procedures. It is extremely important that a hysteroscopic insufflator is used. Death has been reported when laparoscopic CO₂ insufflators were used during hysteroscopy. Flow of CO₂ should be limited to < 100 mL/min, and intrauterine pressure should not exceed 100 mmHg.
- > For use only by physicians trained in hysteroscopy
- > Suspicion of pregnancy should suggest a pregnancy test before the performance of diagnostic hysteroscopy.
- > The high level disinfectant Cidex OPA is contraindicated for patients with a history of bladder cancer

⚠ CAUTION

Cautions in conjunction with hysteroscopic application

Risk to the patient

- > Vaginal ultrasonography before hysteroscopy may identify clinical conditions that will alter patient management
- > Gravity fed intrauterine fluid distention can usually be accomplished with pressures in the range of 35-75 mmHg. Hanging the fluid distention medium 42" above the patient can generate intrauterine pressure of approximately 80 mmHg. Unless the systemic blood pressure is excessive, it is seldom necessary to use pressure greater than 75-80 mmHg.

1.4 Marking

1.4.1 Pictograms and Information on the Device and Packaging

This section describes the different pictograms that feature on the device or the packaging.

	The device is suitable for autoclave sterilization
	Item number
	Serial number
	Lot number
	Quantity of products of a particular item number
	Direction of view
	Manufacturer
	Date of manufacture
	Caution (IEC 60601-1 3rd edition) / Take note of accompanying documents (IEC 60601-1 2nd edition)
	Non-sterile, sterilize prior to each use
	Not for use in vicinity of MRI equipment
	Permissible storage temperature
	Permissible relative air humidity during storage
	Permissible atmospheric pressure during storage
	Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

1.4.2 Pictograms in this Document

This section describes the pictograms used in this document.

General warning sign

Dangerous electrical voltage warning sign

Biohazard warning sign, risk of infection

Hot surfaces warning sign

Laser beam warning sign

Type BF applied part according to IEC 60601-1

Type CF applied part according to IEC 60601-1

Defibrillation-proof type BF applied part

Defibrillation-proof type CF applied part

1.5 Service Department Contact Details

If you have any questions about our devices, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

2 General Safety Information

2.1 Safety Messages in this Document

2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to material damage:

NOTICE

Nature and origin of the risk of material damage

Potential consequences of non-observance

> Preventive action

2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

NOTICE! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

2.2 Product Safety

2.2.1 Basic Safety Information

Our products are developed and manufactured to the highest quality standards.

Risks – despite the utmost quality

Although this device corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing, or maintenance.

Therefore, it is important that you read through these instructions for use carefully and observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

The original packaging should be retained and reused when returning the device for servicing.

Follow the instructions for use for all devices and instruments that are to be used in conjunction with this device.

WARNING! This device is supplied in a non-sterile condition, which carries a risk of infection. Carry out sterile processing of the device and its accessories prior to first use.

WARNING! Unauthorized modifications to the device carry a risk of serious injury. Do not make any unauthorized modifications.

WARNING! Component failure during a surgical procedure carries a risk to the patient. Keep spares on hand for emergency replacement.

WARNING! High-intensity light source. Risk of injury to eyes. Do not look directly into the light exit at the distal end of the endoscope.

WARNING! Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.

WARNING! Using the device during operation of a defibrillator carries a risk to persons. Before discharging a defibrillator, remove device from the operating field.

CAUTION! Rough handling can damage the device and therefore carries a risk to the patient. Handle the device with appropriate care. If the device is dropped or subjected to high mechanical stress, stop using it and send to the manufacturer for inspection.

WARNING! Improper handling, maintenance, and use carries a risk to the patient and user, or can lead to premature wearing of the device.

2.2.2 Staff Qualifications

Observe qualification requirements

For Operation and Reprocessing special qualifications are required in each instance. The qualifications required are stated in the respective chapters of this document and must be observed.

3 Operation

3.1 Safety Instructions

WARNING

Interactions between devices in simultaneous use

Risk to the patient and user, image distortion, damage to the device

- > Ensure that all the devices in use are, as a minimum, compliant with the requirements for applied parts classified as type BF according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used

WARNING

Displaying a recorded image instead of the live image or an altered image orientation

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera
- > Ensure that the live image is displayed in the correct orientation (not mirrored)

WARNING

Use of non-sterile parts

Infection risk for the patient

- > Use only properly reprocessed endoscopes and endoscopic accessories
- > Accessories supplied in a non-sterile condition must be subject to sterile processing before use
- > Always carry out a visual inspection prior to use
- > Anchor the light guide in the operating field and secure against slipping, allowing some slack for maneuvering

WARNING

Recontamination resulting from improper handling

Infection risk for the patient

- > Adhere to the infection control practices that are in place

**⚠ CAUTION**

High temperatures in combination with light source usage, temperatures on the insertion part up to max. 47 °C

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

- > Use only light guides that are suitable for use with the endoscope
- > Avoid the prolonged use of intense light
- > You should always use the lowest light output setting that will allow you to complete the task at hand
- > Do not touch the light source in the immediate vicinity of the lamp
- > Do not touch light guide connectors
- > Do not allow the distal end of the endoscope to come into contact with patient tissue or with combustible or heat-sensitive materials
- > Use endoscopes only in conjunction with a sheath

NOTICE**Improper handling**

Risk of damage to the device

- > Set down carefully without bending or knocking
- > Avoid scratches to the surface
- > Take care when removing instruments from the operating field

3.2 Staff Qualifications

The device may be used only by physicians and surgical assistants who have received instruction in how to use the device and who have the requisite training or advanced training, knowledge, and practical experience in the endoscopic procedure to be applied as defined by the provisions in place at the site of operation.

3.3 Technical Inspection Prior to Use

3.3.1 Visual Inspection

Performing a visual inspection

*A visual inspection must be performed **prior to each use**.*

Proceed as follows:

- | | |
|---|--|
| Is there any external damage? | <p>1. Check the endoscope and all the components to be used for external damage.</p> <p>Do not use the device if it has any sharp edges or corners, bulges, or rough surfaces that might cause injury to the patient.</p> |
| Are the fiber optics in working order? | <p>2. Check that the fiber optics in the endoscope are in full working order.</p> <p>To do this, hold the distal end toward a bright light (not a cold light source) and the light guide connector toward your eyes. Gently move the endoscope from side to side, taking note of the brightness of the fibers.</p> <p>If more than 20 % of the fibers remain dark, the endoscope must not be used and a repair must be arranged.</p> |
| Are the surfaces clean and smooth? | <p>3. Check that the glass surfaces and fiber optic end faces of the endoscope are clean and smooth.</p> <p>Do not use the endoscope if the surfaces are soiled or scratched. Set aside the endoscope and follow the steps described in <u>Chapter 7</u>.</p> |
| Is the image clear, bright, and complete? | <p>4. Make sure you can see a clear, bright, complete image.</p> <p>Look through the eyepiece and assess the image quality.</p> <p>Do not use the endoscope if the image is yellowed, dark, speckled, or cropped. Set aside the endoscope and follow the steps described in <u>Chapter 7</u>.</p> |
| Is the OR equipment mutually compatible? | <p>5. Ensure that all the equipment to be used for the procedure is mutually compatible.</p> <p>Only use the OR equipment if the isolation barrier for protecting against electric shock, as a minimum, complies with the requirements for applied parts classified as type BF in the IEC 60601-1 (see <u>Intended Use</u>).</p> |
- This concludes the inspection.

3.3.2 Function Check

Performing a function check

*A function check must be performed **before every medical procedure.***

Proceed as follows:

1. Connect a camera.
 - Switch on the system 2. Switch on all the system components that you intend using for the procedure.
 - Focus the image 3. Direct the camera head at an object nearby in the room and focus the image.
Check that you are able to achieve a sharp, bright, good quality image.
Do not use the system if the image contains vertical lines or color variations, flickers, or you are unable to achieve a sharp, bright, good quality image.
- This concludes the function check.

3.4 Connecting a Light Guide

Connecting the light guide

Different adapters are needed depending on the type of light guide being connected to the endoscope.

Proceed as follows:

1. Screw the suitable adapter to the light guide connector on the endoscope and, if applicable, to the light guide.
2. Connect the light guide to the light guide connector.

The light guide is connected.

3.5 Post-Use and Pre-Cleaning Activities

Stopping usage and completing pre-cleaning

*Pre-clean the device **immediately after use.***

Proceed as follows:

1. Remove and discard the sterile drape, if present.
2. Decouple the endoscope from the endocoupler.
3. Detach the light guide and all detachable parts.

NOTICE! The eyepiece cup is nondetachable.

4. Pre-clean all used parts of the fully disassembled product at the operating table at the conclusion of the surgery by wiping with a lint-free wipe wetted (i.e. wetted but not dripping; squeeze out excess solution) with an enzymatic cleaning solution until visibly clean.

Cleaning solution, see [Section 4.5](#)

Before wetting the wipe with the enzymatic cleaning solution, it should be prepared per the manufacturer's instructions.

5. Place all the parts into a dry removal container and seal.
6. Arrange for reprocessing, ensuring that all of the device components are reprocessed within 6 hours.

This concludes use and pre-cleaning.

4 Reprocessing

4.1 Safety Instructions

WARNING

If the device is potentially contaminated with Creutzfeldt-Jakob pathogens, it is not suitable for reprocessing

Risk of cross-contamination during reprocessing

- > Do not reprocess devices that are potentially contaminated
- > Devices that are potentially contaminated must be disposed of

WARNING

High Level Disinfection

Risk to patient

- > High level disinfection is NOT recommended if devices are used for hysteroscopy
- > Cidex OPA is contraindicated for patients with a history of bladder cancer

CAUTION

Improper cleaning

Risk to patient as a result of inadequate cleaning, damage to the device

- > Use an automatic washer that meets the performance requirements of ISO 15883-1
- > Load the washer with care: all contents must be rinsed and cleaned fully (no rinsing 'blind spots')
- > Connect instruments with lumens and channels directly to the dedicated rinsing ports on the cleaning and storage basket
- > Remove all stopcocks on the instruments
- > Ensure that the washer is properly maintained
- > Only use cleaning and disinfecting media that are approved for the device
- > The devices must be pre-cleaned immediately after use and reprocessed within 6 hours
- > When pre-cleaning, do not use fixing temperatures in excess of 45 °C
- > When pre-cleaning, do not use fixing detergents (containing aldehyde or alcohol)

NOTICE

Contact with chloride-based solutions

Corrosion and destruction of the device

- > Avoid contact with chloride-based solutions, which can be found, for example, in surgical residues, lotions, drugs, saline solutions, and cleaning and disinfecting media
- > If the instruments or devices should come into contact with chloride-based solutions, rinse thoroughly with deionized water and dry completely

NOTICE**Use of ultrasonic cleaning bath**

Damage to the endoscope

- > The endoscope must not be exposed to ultrasonic cleaning bath

NOTICE**Unsuitable detergents and process chemicals**

Corrosion damage, premature aging, and visible material deterioration

- > Only use cleaning and disinfecting media that are approved for the device
- > Only use process chemicals that are compatible with the device's materials according to the chemical manufacturer's recommendations
- > Follow all the chemical manufacturer's application specifications regarding temperature, concentration, and contact time
- > Do not use process chemicals that can cause stress cracking or brittleness of plastic materials

4.2 Staff Qualifications

In many countries, the qualifications required by personnel responsible for the reprocessing of medical devices are regulated by law.

In any event, the reprocessing of medical devices must only be carried out by qualified personnel who have the necessary knowledge and expertise.

Professional expertise may be obtained through specialized further education or through training and practical experience, and enhanced through appropriate further training if necessary.

4.3 Validated Procedures

The effectiveness of the following procedures:

- Manual cleaning
- Automated cleaning
- High Level Disinfection (HLD, only if devices are used for cystoscopy)
- Sterilization

as described in this document, has been fully validated.

Operator responsibility

It is the operator's responsibility to introduce, document, implement, and maintain a validated reprocessing procedure.

Make sure that the equipment used for reprocessing is properly maintained.

4.4 The Process Flow

Pre-cleaning The device has to be prepared for reprocessing immediately after use by pre-cleaning.

See [Section 3.5](#) for pre-cleaning instructions.

Steps The reprocessing procedure described in this document consists of the following steps:

- Pre-cleaning immediately after use
- Manual or automated cleaning
- Sterilization or HLD (only if devices are used for cystoscopy)

Protective clothing During reprocessing, staff should wear personal protective clothing.

Careful cleaning The device must be cleaned thoroughly at the beginning of reprocessing. It is essential that the sterilization medium accesses all parts of the endoscope.

Automated cleaning and steam sterilization The optimum and most reliable reprocessing results are achieved using automated cleaning with subsequent steam sterilization using the fractionated vacuum steam procedure.

Legislation and standards Adhere to national statutory regulations, national and international standards and directives, and the infection control practices that are in place at your institution for reprocessing.

Reprocessing prior to each use New and returned devices following repair have to undergo complete reprocessing before first use.

The alternating use of different reprocessing methods can cause premature aging of the product.

Excessive doses of neutralizers and cleaning agents can damage the device and cause laser engravings to fade.

www.a-k-i.org Please see www.a-k-i.org for more detailed information about hygienically safe, material-protecting reprocessing methods that also preserve the instruments' value.

4.5 Cleaning and Disinfecting Agents

The following products can be used for **manual cleaning**:

- Cidezyme/Enzol (Johnson & Johnson)

The following products can be used for **automated cleaning**:

- neodisher MediClean forte 0.5 %
(Chem. Fabrik Dr. Weigert GmbH & Co. KG)

The following product can be used for **high-level disinfection**:

- CIDEX OPA (Johnson & Johnson)

Wherever possible, use the above cleaning and disinfecting agents only.

Follow the manufacturer's instructions! Before use, read the manufacturer's user information carefully and follow the specifications regarding concentration, temperature, usage time, water quantities, and contact time.

4.6 Manual Cleaning in Soaking Bath

4.6.1 Overview

Stage	Work Step	Temperature (°C/°F)	Time (min.)	Water Quality	Cleaning Solution
I	Cleaning	per manufacturer's instructions	2-5	Tap water	Enzymatic cleaning solution
II	Rinse 2x	< 45/113	2x ≥1	Tap water	---
III	Final rinse	< 45/113	≥1	High purity water	---
IV	Drying	---	---	---	---

Table 4-1: Overview of manual cleaning in soaking bath.

4.6.2 Performance

Complete stage I: Cleaning

Prepare the cleaning solution per the manufacturer's instructions.

Proceed as follows:

NOTICE! Damage resulting from rough handling. Handle the device with appropriate care. Do not knock or bend. Set down carefully.

NOTICE! (Endoscope) Scratch-sensitive surfaces. Risk of abrasion. Do not use metal brushes, metal objects, or abrasive cleaners. To remove soiling on optical surfaces, use pads with cleaning solution only.

NOTICE! Scratch-sensitive surfaces. Risk of corrosion due to abrasion. Do not use metal brushes, metal objects, or abrasive cleaners to clean the instruments. To remove soiling use pads with cleaning solution only.

Immerse for
2-5 min

1. Fully immerse all parts of the completely disassembled device for **2-5** minutes in the cleaning solution.

All accessible surfaces must remain immersed in the cleaning solution bath throughout the entire cleaning time.

Wipe or brush
external surfaces

2. After soaking in the cleaning solution, clean all outer surfaces while all parts are fully immersed in the cleaning solution using a soft, clean, lint-free wipe or a soft bristled brush until visibly clean.

Complete stage II: Rinsing

2 rinses in tap
water of at least
1 minute each

1. Fully immerse all parts in a tap water bath (less than (< 45° C / 113° F) and thoroughly rinse all accessible surfaces twice for a minimum of **one** minute each.

Use fresh water for each rinse.

Allow water to drip off for a sufficient length of time.

Complete stage III: Final rinse

For at least 1 min.
in high purity
water

1. Completely immerse the product in high purity water (< 45 °C / 113 °F) and rinse all accessible surfaces completely for at least **one** minute.

Allow water to drip off for a sufficient length of time.

Complete stage IV: Dry

1. Thoroughly dry all parts with a clean, soft, lint-free wipe or operating room towel. Use medical-quality filtered compressed air (max = 0.5 bar) to dry the lumen.
- Visual inspection
2. After drying, visually inspect the components in a well-lit area (see [Section 3.3.1](#)); they should be dry, free of visible residue and undamaged.

Use a magnifying glass. Repeat the cleaning process if necessary.

Immediately set aside any damaged components.
- This concludes the manual cleaning process.

4.7 Automated Cleaning

4.7.1 Overview

Stage	Step	Temperature (°C/°F)	Time (min.)	Cleaning Solution / Water Quality
I	Pre-rinse	Cold	2	Tap water
II	Clean	60/140	10	Alkaline cleaning solution
III	Intermediate rinse	per the equipment manufacturer's standard cycle	1	Tap water
IV	Intermediate rinse	per the equipment manufacturer's standard cycle	1	Tap water
V	Final rinse	90/194	5	Deionized or purified water
VI	Dry	99/210	30	---

Table 4-2: Automated cleaning overview.

Note: Machine reprocessing can cause color anodized or plastic components (e.g. standard rings and the eyepiece cup) to fade.

4.7.2 Performance

Performing automated cleaning

Use a washer that has been validated as effective and meets the performance requirements of ISO 15883-1 (or the respective country-specific version thereof).

Dry removal of the products is recommended prior to automated cleaning.

In the case of wet removal, do not use foaming detergents and rinse the products thoroughly prior to the automated cleaning.

Proceed as follows:

1. Place all parts in a cleaning and storage basket.
2. Put the cleaning and storage basket into the washer.
3. Start the cleaning cycle described in [Section 4.7.1](#) in accordance with the manufacturer's instructions and instructions for use for the washer.

CAUTION! Risk of scalding when unloading the washer. Be sure to wear gloves.

4. Take the cleaning and storage basket out of the washer.

- Visual inspection
5. After drying, visually inspect the device in a well-lit area (see [Section 3.3.1](#)); it should be dry, free of visible residue and undamaged.

Use a magnifying glass. Repeat the cleaning process if necessary.

Immediately set aside any damaged components.

- This concludes the automated cleaning process.

4.8 High Level Disinfection (HLD)

4.8.1 Overview

Stage	Work Step	Temperature (°C/°F)	Time (min.)	Disinfectant solution / Water quality
I	Disinfection	≥ 20/68	12	CIDEX OPA (undiluted)
II	Rinse 3x	< 45/113	3x ≥1	Sterile water
III	Dry	Room temperature	---	---

Table 4-3: Overview of HLD.

4.8.2 Performance

Perform the HLD instead of steam sterilization only if steam sterilization is not possible.

WARNING! High level disinfection is NOT recommended if devices are used for hysteroscopy.

CAUTION! Residues from cleaning agents and organic material can negatively influence the disinfection result. Perform the HLD only on a properly cleaned and dried device.

Complete stage I: Disinfection

Prepare the disinfectant solution per the manufacturer's instructions.

Proceed as follows:

1. Ensure that the disinfectant solution contains the required concentration of active substances. Use an original CIDEX OPA test strip.
2. Fully immerse all parts of the completely disassembled device in the disinfectant solution.

All accessible surfaces must remain immersed in the disinfectant solution bath throughout the entire disinfecting time.

Remove all adherent air bubbles from the component surfaces.

Complete stage II: Rinsing

Proceed as follows:

3 rinses for at least 1 min.

1. Fully immerse all parts in a large amount of sterile water (< 45° C / 113° F) and thoroughly rinse all accessible surfaces three times for a minimum of **one** minute each.

Use fresh water for each rinse.

Allow water to drip off for a sufficient length of time.

Complete stage III: Dry

Proceed as follows:

1. Thoroughly dry all parts with a sterile, soft, lint-free wipe or operating room towel. Use medical-quality filtered compressed air (max = 0.5 bar) to dry the lumen.

Visual inspection

2. After drying, visually inspect the components in a well-lit area (see [Section 3.3.1](#)); they should be dry, free of visible residue and undamaged.

Use a magnifying glass. Repeat the disinfection process if necessary.

Immediately set aside any damaged components.

This concludes the high-level disinfection.

4.9 Sterilization

4.9.1 Steam Sterilization

Performing steam sterilization

AUTOCLAVE

Steam sterilization should be performed only for products that are labeled as autoclavable.

WARNING! The biological indicators, sterilization wraps, pouches, and containers used must be FDA-approved for sterilization under the parameters given.

CAUTION! Residues from cleaning agents and organic material can negatively influence the sterilization result. Perform the sterilization only on a properly cleaned and dried device.

Proceed as follows:

1. Place all parts in a cleaning and storage basket.
2. Double-wrap the basket with two single layers of FDA-cleared polypropylene wrap using a sequential double wrapping technique.

CAUTION! Risk of scalding when loading the sterilizer. Wear suitable gloves.

3. Load the wrapped basket in the sterilizer.

Load the sterilizer as per the validated loading plan.

Follow the manufacturer's instructions!

4. Start the sterilization cycle in accordance with the manufacturer's instructions and sterilizer instructions for use.

Use the following validated parameters for steam sterilization using an autoclave cycle with prevacuum air removal:

Cycle type:	prevacuum
Pulses:	4
Temperature:	270 °F (132 °C)
Holding time:	min. 4 minutes (effective sterilization time)
Minimum dry time:	at least 30 min.

Note: It is the responsibility of the user to establish that the sterilizer used meets the steam sterilization parameters above.

CAUTION! Risk of scalding when unloading the sterilizer. Wear suitable gloves.

NOTICE! Products can be shock sensitive when hot. Avoid knocking or shaking.

NOTICE! Damage resulting from a sudden change in temperature. Allow products to cool to room temperature; do not use additional cooling measures.

5. Remove the sterilized products from the sterilizer.

Ensure that the products remain sterile after reprocessing.

- This concludes sterilization.

5 Maintenance and Repairs

5.1 Troubleshooting Table

Issue	Possible Causes	Fixes
Image cloudy	Glass surfaces soiled	Clean manually then reprocess, check the water quality
	Stubborn deposits on glass surfaces	Remove deposits as per the Cleaning and Disinfection chapter, check the water quality
	Leaky, defective lens system	Send the endoscope for repair
Image too dark, insufficient illumination	Glass surfaces soiled	Clean manually then reprocess, check the water quality
	Stubborn deposits on glass surfaces	Remove deposits as per the Cleaning and Disinfection chapter, check the water quality
	Unsuitable light guide	Use a suitable light guide
	Light guide incorrectly attached to endoscope	Check that the light guide is seated properly
	Defective fiber optics	Check the fiber optics as per Section 3.3
	Defective light guide or light source	Check the light guide and light source
Image is yellowed	Soiled fiber optics	Clean manually then reprocess, check the water quality
	Light guide soiled or defective	Check light guide (e.g. by illuminating white surface)
Corrosion, staining, discoloration	Inadequate cleaning (e.g. remaining protein residue)	Clean manually, rub thoroughly if necessary, then reprocess
	Insufficient rinsing between the reprocessing stages, especially prior to sterilization	Rinse sufficiently between the reprocessing stages
	Excessive chloride concentration	Check water quality
	Excessive concentration of minerals (e.g. limescale) or organic substances	Check the water quality, only use demineralized water if necessary

Issue	Possible Causes	Fixes
	Heavy metal ions and/or silicates, elevated concentration of iron, copper, and manganese in the water or sterilizer steam	Check the water quality, only use demineralized water if necessary
	Contaminated, too frequently used cleaning/disinfecting solution	The cleaning/disinfectant solution should be replaced regularly
	Extraneous rust, e.g. resulting from rust contained in steam or damaged or non-corrosion-resistant instruments reprocessed at the same time	Avoid contact with other components
	Contact corrosion	Avoid contact with other components

Table 5-1: Troubleshooting.

5.2 Repair

Should you need to arrange a repair for the device, please contact our Service Department.

When you return the device, please also enclose a detailed **description of the fault**.

WARNING

Contaminated device

Risk of infection

> Reprocess the device properly prior to shipping

Only ship devices that have been properly reprocessed.

The device should be returned in its original packaging, if possible.

Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return these to sender.

6 Product Data

6.1 Technical Data

Item No.	STS-2040-302ASA
Direction of view	0°
Field of view	standard
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-302AA
Direction of view	0°
Field of view	wide angle
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-302BA
Direction of view	30°
Field of view	wide angle
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-302BSA
Direction of view	30°
Field of view	standard
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-302ESA
Direction of view	12°
Field of view	standard
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-302CA
Direction of view	70°
Field of view	wide angle
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2029-302AA
Direction of view	0°
Field of view	wide angle
Working length	302 mm
Outer diameter of insertion portion	2.9 mm

Item No.	STS-2029-302BA
Direction of view	30°
Field of view	wide angle
Working length	302 mm
Outer diameter of insertion portion	2.9 mm

Item No.	STS-2040-298DSA-RWF
Direction of view	5°
Field of view	standard
Working length	298 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-298ASA-RWF
Direction of view	0°
Field of view	standard
Working length	298 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-300GSA-RWF
Direction of view	25°
Field of view	standard
Working length	300 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-303CA-RWF
Direction of view	70°
Field of view	wide angle
Working length	303mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-280ASA-OLY
Direction of view	0°
Field of view	standard
Working length	280 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-280ESA-OLY
Direction of view	12°
Field of view	standard
Working length	280 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-283BSA-OLY
Direction of view	25°
Field of view	standard
Working length	283 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-283BA-OLY
Direction of view	30°
Field of view	wide angle
Working length	283 mm
Outer diameter of insertion portion	4.0 mm

Item No.	STS-2040-283CA-OLY
Direction of view	70°
Field of view	wide angle
Working length	283 mm
Outer diameter of insertion portion	4.0 mm

6.2 Ambient Conditions

Transport and storage conditions

Temperature	-20 °C to +70 °C
Rel. humidity	5 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

Reprocessed devices must be protected against recontamination, in a dark, dry, well-ventilated, dust-free and light-controlled location with a consistent temperature.

Direct sunlight, high temperatures, high humidity or radiation can damage the device or present a risk of infection.

When placing into storage, make sure that the device cannot be damaged by other instruments. It is therefore best to store the device individually or use containers in which it can be secured.

Operating conditions

Temperature	+15 °C to +40 °C
Rel. humidity	10 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

6.3 Spare Parts and Accessories

Use original spare parts and accessories only.

Photo	Designation	Item Number
	Light guide adapter for Storz systems	STS-ACC-0000-LGUA
	Light guide adapter for Wolf systems	STS-ACC-0000-LGUA-RWF

Table 6-1: Spare parts and accessories.

7 Disposal

WARNING

Contaminated device

Risk of infection

- > The device must be reprocessed prior to disposal

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

STRAUSS SURGICAL

3020 NW 82 Ave.
Miami, FL 33122

Phone: 305 436 0599
Fax: 305 436 0399
Web: www.straussurgical.com

Instructions For Use

Cystoscopes/Hysteroscopes

(Not valid for USA)

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Cistoscopios/Histeroscopios

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Cistoscópios/Histeroscópios

TPA767-025-34_C
Version: C
2020 - Dec - 18

	Content	Page
0	Important Information about this Document	7
0.1	Scope of Validity, Identification, Purpose.....	7
0.2	Target Group.....	7
0.3	Using and Storing this Document.....	7
1	General Information about the Device	8
1.1	Scope of Delivery	8
1.2	Product Description	8
1.2.1	Performance Characteristics and Function	8
1.2.2	Visual Overview	8
1.3	Usage.....	9
1.3.1	Intended Use.....	9
1.3.2	Indications.....	9
1.3.3	Contraindications.	10
1.4	Marking	11
1.4.1	Pictograms and Information on the Device and Packaging	11
1.4.2	Pictograms in this Document.....	12
1.5	Service Department Contact Details	12
2	General Safety Information	13
2.1	Safety Messages in this Document.....	13
2.1.1	Safety Messages at the Start of a Chapter.....	13
2.1.2	Safety Messages in the Body of the Text.....	13
2.2	Product Safety	14
2.2.1	Basic Safety Notices.	14
2.2.2	Staff Qualifications.....	14
3	Operation.	15
3.1	Safety Instructions.....	15
3.2	Application Notes.....	17
3.3	Staff Qualifications.....	18
3.4	Technical Inspection Prior to Use.....	18
3.4.1	Visual Inspection.....	18
3.4.2	Function Check	19
3.5	Connection of a Light Guide	19
3.6	Post-Use and Pre-Cleaning Activities	19
4	Reprocessing	20
4.1	Safety Notices.....	20
4.2	Staff Qualifications.....	21
4.3	Validated Procedures.....	21
4.4	The Process Flow.	22
4.5	Cleaning and Disinfecting Agents.....	22
4.6	Manual Cleaning and Disinfection.....	23
4.6.1	Overview	23
4.6.2	Performance.....	23
4.7	Automated Cleaning and Thermal Disinfection	25
4.7.1	Overview	25
4.7.2	Performance.....	25
4.8	Sterilization	26
4.8.1	Steam Sterilization.....	26
5	Maintenance and Repairs	27
5.1	Troubleshooting Table.	27
5.2	Repair.....	28
6	Product Data	29
6.1	Technical Data.....	29
6.2	Ambient Conditions	31
6.3	Spare Parts and Accessories.....	31
7	Disposal	32

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0 Important Information about this Document

0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the following products:

STS-2040-302ASA, STS-2040-302AA,
STS-2040-302BA, STS-2040-302BSA,
STS-2040-302ESA, STS-2040-302CA,
STS-2029-302AA, STS-2029-302BA,
STS-2040-298DSA-RWF, STS-2040-298ASA-RWF,
STS-2040-300GSA-RWF, STS-2040-303CA-RWF,
STS-2040-280ASA-OLY, STS-2040-280ESA-OLY,
STS-2040-283BSA-OLY, STS-2040-283BA-OLY,
STS-2040-283CA-OLY

These instructions for use are an integral component of the product and contain all the information required by users and operators for safe and proper use.

These instructions for use do not describe endoscopic procedures or techniques during surgery.

0.2 Target Group

These instructions for use are intended for physicians, medical assistants, and employees of central sterile supply who are entrusted with the operation, handling, and reprocessing of the device.

0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

1 General Information about the Device

1.1 Scope of Delivery

Check scope of delivery! Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

1.2 Product Description

1.2.1 Performance Characteristics and Function

The cystoscopes/hysteroscopes described herein are rigid endoscopes for visualizing the urethra, the urinary bladder or uterus, the fallopian tubes and the cervical canal during the performance of endoscopic procedures in urology or gynecology.

Structure and function A rigid endoscope consists of a fiber optic cable and sensitive image transmission system with eyepiece.

The fiber optic cable is used to illuminate the site inside the body. The connector for connecting the light guide to the light source is situated at the proximal end of the endoscope. The adapters required to connect the light guide are included in the scope of delivery.

The distal end of the endoscope features an object lens that captures the image from inside the body. The image is sent through the image transmission system to the eyepiece.

The eyepiece cap connects the eyepiece to an endocoupler that transmits the image to a camera. The camera's controller converts the signal so that the image can be displayed on a monitor.

Note: During this medical application, adhere to national statutory regulations and directives.

1.2.2 Visual Overview

Cystoscope/hysteroscope


Figure 1-1: Cystoscope/hysteroscope.

- 1 Eyepiece cap
- 2 Light guide connector (ACMI)
- 3 Adapter for connecting Storz light guides
- 4 Adapter for connecting Wolf light guides
- 5 Endoscope sheath
- 6 Endoscope tip

1.3 Usage

1.3.1 Intended Use

The cystoscopes/hysteroscopes are indicated to provide the user with the means for endoscopic diagnostic and therapeutic surgical procedures.

Examples of the use of the devices include visualization and tissue manipulation as long as this is considered appropriate by the surgeon.

The cystoscopes/hysteroscopes are intended to be used in general urological and gynecological surgery through a minimally invasive approach by utilizing natural orifices to access the surgical site.

The devices are intended to be used in combination with medical electrical devices that meet at **minimum** the type **BF** classification requirements according to IEC 60601-1 for elevated protection from electric shock.


This insulation barrier can be realized by the medical electrical devices themselves or by the connection cables to the endoscope and must be assured for every connection between the endoscope and connected devices.

Devices or connection cables that meet the type BF classification requirements feature the adjacent symbol.

Devices and connection cables that meet the requirements for applied parts classified as type CF offer an even higher level of protection. They bear the adjacent heart symbol and are approved for direct conductive contact with the heart.

Defibrillation-proof devices or connection cables are also compatible. These will be marked with one of the adjacent symbols depending on the degree of protection against electric shock offered.

1.3.2 Indications

The indications for an endoscopic application depend on the medical condition of the patient and the individual risk/benefit analysis by the surgeon.

Indications for surgical implementation include:

- Examination of the cavities during a diagnostic or therapeutic minimally invasive treatment (e.g. the removal of myomata)
- Anatomical changes (e.g. stenoses) or foreign bodies in the cavities

During this medical application, adhere to national statutory regulations and directives.

1.3.3 Contraindications

Contraindications may exist depending on the patient's general condition or specific symptoms.

The decision to perform a procedure remains with the responsible surgeon and must be made on the basis of an individual risk/benefit analysis.

General contraindications

- Acute pelvic inflammatory disease

Contraindications to Hysteroscopy

Hysteroscopy may be contraindicated by the following clinical conditions, depending on their severity or extent:

- Inability to distend the uterus
- Cervical stenosis
- Cervical/vaginal infection
- Uterine bleeding or menses
- Known pregnancy
- Invasive carcinoma of the cervix
- Recent uterine perforation
- Medical contraindication or intolerance to anesthesia

Contraindication to Endometrial Ablation

The following are clinical conditions that can significantly complicate hysteroscopic endometrial ablation:

- Adenomatous Endometrial Hyperplasia
- Uterine Leiomyoma
- Severe Adenomyosis
- Pelvic Pain (Subtle PID)
- Uterine anomalies

Contraindication to Hysteroscopic Myomectomy

The following are clinical conditions that can significantly complicate hysteroscopic myomectomy:

- Severe anemia
- Inability to circumnavigate a myoma due to myoma size (e.g. predominantly intramural myomas with small submucous components).

1.4 Marking

1.4.1 Pictograms and Information on the Device and Packaging

This section describes the different pictograms that feature on the product or packaging.

	The device is suitable for autoclave sterilization
	Item number
	Serial number
	Lot number
	Quantity of products of a particular item number
	Direction of view
	Manufacturer
	Date of manufacture
	Caution (IEC 60601-1 3rd edition) / Take note of accompanying documents (IEC 60601-1 2nd edition)
	Non-sterile, sterilize prior to each use
	MR unsafe
	Permissible storage and transport temperature
	Permissible relative air humidity during storage and transport
	Permissible atmospheric pressure during storage and transport
	Caution: Federal (US) law restricts this device to sale by or on order of a physician

1.4.2 Pictograms in this Document

This section describes the pictograms used in this document.

General warning sign

Dangerous electrical voltage warning sign

Biohazard warning sign, risk of infection

Hot surfaces warning sign

Laser beam warning sign

Type BF applied part according to IEC 60601-1

Type CF applied part according to IEC 60601-1

Defibrillation-proof type BF applied part

Defibrillation-proof type CF applied part

1.5 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

2 General Safety Information

2.1 Safety Messages in this Document

2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to material damage:

NOTICE

Nature and origin of the risk of material damage

Potential consequences of non-observance

> Preventive action

2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

NOTICE! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

2.2 Product Safety

2.2.1 Basic Safety Notices

Our products are developed and manufactured to the highest quality standards.

Risks- despite highest of quality standards

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

The original packaging should be retained and reused when returning the device for servicing.

Follow the instructions for use for all devices and instruments that are to be used in conjunction with this device.

WARNING! This device is supplied in a non-sterile condition. Risk of infection. Carry out reprocessing of the device and its accessories prior to first use.

WARNING! Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.

WARNING! Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.

WARNING! High-intensity light source. Risk of injury to eyes. Do not look directly into the light exit at the distal end of the endoscope.

WARNING! Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.

WARNING! Using the device during operation of a defibrillator. Risk to persons. Before discharging a defibrillator, remove device from the surgical field.

CAUTION! Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.

WARNING! Improper handling, maintenance, and use carries a risk to the patient and user, or can lead to premature wearing of the device.

2.2.2 Staff Qualifications

Observe qualification requirements

Specific qualifications are required for operation and reprocessing of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

3 Operation

3.1 Safety Instructions

WARNING

Interactions between devices in simultaneous use

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the requirements for applied parts classified as type BF according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used

WARNING

Displaying a recorded image instead of the live image or an altered image orientation

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera
- > Ensure that the live image is displayed in the correct orientation (not mirrored)

WARNING

Use of non-sterile parts

Infection risk for the patient

- > Use only properly reprocessed endoscopes and endoscopic accessories
- > Accessories supplied in a non-sterile condition must be subject to reprocessing before use
- > Always carry out a visual inspection prior to use
- > Anchor the light guide in the surgical field and secure against slipping, allowing some slack for maneuvering

WARNING

Recontamination resulting from improper handling

Infection risk for the patient

- > Adhere to the infection control practices that are in place

**⚠ CAUTION**

High temperatures resulting from light source usage, temperatures on the insertion portion of up to 47 °C

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

- > Use only light guides that are suitable for use with the endoscope
- > Avoid the prolonged use of intense light
- > You should always use the lowest possible light output setting that will allow you to illuminate the target area
- > Do not touch the light source in the vicinity of the lamp
- > Do not touch light guide connectors
- > Do not allow the distal end of the endoscope to come into contact with patient tissue or with combustible or heat-sensitive materials
- > Always use sheath with endoscope

NOTICE**Improper handling**

Risk of damage to the product

- > Set down carefully without bending or knocking
- > Avoid scratching the surface
- > Take care when removing instruments from the surgical field

3.2 Application Notes

**WARNING****Use of CO₂ gas as a distention medium during operative hysteroscopy**

Risk of gas embolization

- > Do not use CO₂ gas as a distention medium in operative hysteroscopy

**WARNING****Use of CO₂ insufflators in diagnostic hysteroscopy**

Risk of gas embolization

- > Do not use a laparoscopic insufflator
- > Use a hysteroscopic insufflator
- > Limit insufflator working pressure to the range of 70-100 mmHg
- > Limit intrauterine pressure to a maximum of 100 mmHg
- > Limit flow of CO₂ to a maximum of 100 mL/min

**WARNING****Resorption of electrolyte-free solutions, high pressure during hysteroscopy**

Hypotonic hyperhydration, hyponatremia, hypothermia, uterine perforation, pulmonary edema, cerebral edema, potential injury to adjacent organs

- > Keep the intrauterine distention pressure as low as possible
- > Keep the duration of the procedure as short as possible
- > Strictly supervise the fluid intake and output during the procedure
- > When excessive fluid deficit occurs quickly complete or interrupt procedure

**WARNING****Contact of bladder cancer patients with Cidex OPA**

Possible anaphylactic reaction

- > Do not bring patients with a history of bladder cancer in contact with devices that were disinfected using Cidex OPA

Vaginal ultrasonography before hysteroscopy may identify clinical conditions that will alter patient management.

Gravity fed intrauterine fluid distention can usually be accomplished with pressures in the range of 35-75 mmHg.

Hanging the fluid distention medium 42" above the patient can generate intrauterine pressure of approximately 80 mmHg.

Unless the systemic blood pressure is excessive, it is seldom necessary to use pressure greater than 75-80 mmHg.

Suspicion of pregnancy should suggest a pregnancy test before the performance of diagnostic hysteroscopy.

3.3 Staff Qualifications

The device may be operated only by physicians and medical assistants who have received instruction in how to use the device and who have the requisite training or advanced training, knowledge, and practical experience in the endoscopic procedure to be applied as defined by the provisions in place at the site of operation.

3.4 Technical Inspection Prior to Use

3.4.1 Visual Inspection

Performing a visual inspection

*A visual inspection must be performed **prior to each use**.*

Proceed as follows:

- | | |
|---|--|
| Is there any external damage? | <p>1. Check the endoscope and all the components to be used for external damage.</p> <p>Do not use the product if it has any sharp edges or corners, bulges, or rough surfaces that might cause injury to the patient.</p> |
| Are the fiber optics in working order? | <p>2. Check that the fiber optics in the endoscope are in full working order.</p> <p>To do this, hold the distal end toward a bright light (not a cold light source) and the light guide connection toward your eyes. Gently move the endoscope from side to side, taking note of the brightness of the fibers.</p> <p>If more than 20 % of the fibers remain dark, the device must not be used. Submit the endoscope to repair.</p> |
| Are the surfaces clean and smooth? | <p>3. Check that the glass surfaces and fiber optic end faces of the endoscope are clean and smooth.</p> <p>Do not use the endoscope if the surfaces are soiled or scratched. Set aside the endoscope and follow the steps described in section 5.1.</p> |
| Is the image clear, bright, and complete? | <p>4. Make sure you can see a clear, bright, complete image.</p> <p>Look through the eyepiece and assess the image quality.</p> <p>Do not use the endoscope if the image is yellowed, dark, speckled, or cropped. Set aside the endoscope and follow the steps described in section 5.1.</p> |
| Is the OR equipment mutually compatible? | <p>5. Check that all the equipment to be used for the procedure is mutually compatible.</p> <p>Only use the OR equipment if the isolation barrier for protecting against electric shock meets at minimum the type BF classification requirements according to IEC 60601-1 (see Intended Use).</p> |
| | <p><input checked="" type="checkbox"/> This concludes the inspection.</p> |

3.4.2 Function Check

Performing a function check

*A function check must be performed **prior to each medical procedure.***

Proceed as follows:

1. Connect a camera.
 - Switch on the system 2. Switch on all the system components that you intend using for the procedure.
 - Focus the image 3. Direct the camera head at an object in the room and focus the image.
 4. Direct the camera head at a nearby object in the room and focus the image.
- Check that you are able to achieve a sharp, bright, good quality image.
- Do not use the system if the image contains vertical lines or color variations, flickers, or you are unable to achieve a sharp, bright, good quality image.
- This concludes the function check.

3.5 Connection of a Light Guide

Connecting the light guide

Different light guides require different adapters to connect to the endoscope.

Proceed as follows:

1. Screw the suitable adapter to the light guide connector on the endoscope and, if applicable, to the light guide.
 2. Connect the light guide to the light guide connector.
- The light guide is connected.

3.6 Post-Use and Pre-Cleaning Activities

Stopping usage and completing pre-cleaning

*Complete the pre-cleaning process **immediately after use.***

Proceed as follows:

1. Remove and discard the sterile drape, if present.
2. Uncouple the endoscope from the endocoupler.
3. Detach the light guide and all detachable parts.

NOTICE! The eyepiece cap is non-detachable.

4. Remove visible surgical residue as completely as possible using a moist (but not dripping), lint-free wipe.

Cleaning agents, see [section 4.5](#)

Before wetting the wipe with an enzymatic cleaning solution, it should be prepared per the manufacturer's instructions.

5. Place all the parts into a dry removal container and seal.
6. Arrange for reprocessing, ensuring that all of the product components are reprocessed within 6 hours.

This concludes use and pre-cleaning.

4 Reprocessing

4.1 Safety Notices

WARNING

Contact of bladder cancer patients with Cidex OPA

Possible anaphylactic reaction

- > Do not bring patients with a history of bladder cancer in contact with devices that were disinfected using Cidex OPA

CAUTION

Improper cleaning and disinfection

Risk to patient as a result of inadequate cleaning and disinfection, damage to the device

- > Use an automatic washer/disinfector that meets the performance requirements of ISO 15883-1
- > Load the washer/disinfector with care: all contents must be rinsed and cleaned fully (no rinsing 'blind spots')
- > Connect instruments with lumens and channels directly to the dedicated ports on the reprocessing tray
- > Disassemble all valves on instruments
- > Ensure that the washer/disinfector is properly maintained
- > Only use cleaning and disinfecting agents that are approved for the product
- > The devices must be pre-cleaned immediately after use and reprocessed within 6 hours
- > When pre-cleaning, do not use fixing temperatures in excess of 45 °C
- > When pre-cleaning, do not use fixing cleaning and disinfecting agents (active ingredients: Aldehyde, alcohol)

NOTICE

Contact with chloride-based solutions

Corrosion and destruction of the device

- > Avoid contact with chloride-based solutions, which can be found, for example, in surgical residue, lotions, drugs, saline solutions, and cleaning and disinfecting agents
- > After contact with chloride-based solutions, rinse products thoroughly with deionized water and dry completely

NOTICE

Use of ultrasonic bath

Damage to the endoscope

- > Do not expose endoscope to ultrasonic bath

NOTICE**Unsuitable cleaning agents and process chemicals**

Corrosion damage, premature aging, and visible material changes

- > Only use cleaning and disinfecting agents that are approved for the product
- > Only use process chemicals that are compatible with the device's materials according to the chemical manufacturer's recommendations
- > Follow all the chemical manufacturer's application specifications regarding temperature, concentration, and contact time
- > Do not use process chemicals that can cause stress cracking or brittleness of plastic materials

4.2 Staff Qualifications

In many countries, the qualifications required by personnel responsible for the reprocessing of medical devices are regulated by law.

In any event, the reprocessing of medical devices must always fall under the responsibility of qualified personnel who have the necessary knowledge and expertise.

This knowledge and expertise can be acquired by completing further training in this field or as a result of completing a dedicated qualification with practical experience, supplemented by appropriate further training measures, where required.

4.3 Validated Procedures

The effectiveness of the following procedures:

- Manual cleaning and disinfection
- Automated cleaning and thermal disinfection
- Sterilization

as described in this document, has been fully validated.

Operator
responsibility

It is the operator's responsibility to introduce, document, implement, and maintain a validated reprocessing procedure.

Make sure that the equipment used for reprocessing is properly maintained.

4.4 The Process Flow

Pre-cleaning The product has to be prepared for reprocessing immediately after use by pre-cleaning.

Refer to [section 3.6](#) for pre-cleaning instructions.

Work steps The reprocessing procedure described in this document consists of the following steps:

- Pre-cleaning immediately after use
- Cleaning and disinfection (manual or automated)
- Sterilization

Protective clothing During reprocessing, staff should wear personal protective clothing.

Careful cleaning The product must be cleaned thoroughly at the beginning of reprocessing. It is essential that the sterilization medium reaches all parts of the product.

Automated cleaning and steam sterilization The optimum and most reliable reprocessing results are achieved using automated cleaning and disinfection with subsequent steam sterilization using the fractionated vacuum steam procedure.

Legislation and standards Adhere to national statutory regulations, national and international standards and directives, and the infection control practices that are in place at your institution for reprocessing.

Reprocessing prior to each use New and returned products following repair have to undergo complete reprocessing before first use.

The alternating use of different reprocessing methods can cause premature aging of the product.

Excessive doses of neutralizers and cleaning agents can damage the product and cause laser engravings to fade.

www.a-k-i.org Please see www.a-k-i.org for more detailed information about hygienically safe, value-preserving reprocessing methods that also protect the materials.

4.5 Cleaning and Disinfecting Agents

The following products can be used for **manual cleaning**:

- Cidezyme/Enzol (Johnson & Johnson)

The following products can be used for **manual disinfection**:

- Cidex OPA (Johnson & Johnson)

The following products can be used for **automated cleaning**:

- neodisher MediClean forte 0.5 %
(Chem. Fabrik Dr. Weigert GmbH & Co. KG)

Wherever possible, use the above cleaning and disinfecting agents only.

Follow the manufacturer's instructions! Before use, read the manufacturer's user information carefully and follow the specifications regarding concentration, temperature, usage time, water quantities, and contact time.

4.6 Manual Cleaning and Disinfection

4.6.1 Overview

Stage	Work step	Temperature (°C/°F)	Time (min)	Water quality	Cleaning/disinfectant solution
I	Cleaning	(as per manufacturer's instructions)	2-5	Tap water	Enzymatic cleaning solution
II	Rinsing 2x	< 45/113	2x ≥ 1	Tap water	---
III	Disinfection	(as per manufacturer's instructions)	12	---	Disinfectant solution (undiluted)
IV	Rinsing 2x	< 45/113	2x ≥ 1	Tap water	---
V	Final rinsing	< 45/113	≥ 1	Deionized water	---
VI	Drying	---	---	---	---

Table 4-1: Overview of manual cleaning and disinfection.

4.6.2 Performance

Complete stage I: Cleaning

Prepare the cleaning solution and disinfectant solution as per the manufacturer's instructions.

Proceed as follows:

NOTICE! Damage resulting from rough handling. Handle the product with care. Do not knock or bend. Set down carefully.

NOTICE! Scratch-sensitive surface. Risk of abrasion. Do not use metal brushes, metal objects, or abrasive cleaners. To remove soiling on optical surfaces, use pads soaked with cleaning solution only.

- | | | |
|---------------------------------|----|---|
| Immerse for 2-5 min | 1. | Fully immerse all parts of the fully disassembled product that were used in the cleaning solution bath for two to five minutes.

Make sure that all accessible surfaces remain immersed in the cleaning solution throughout the entire cleaning time. |
| Wipe or brush external surfaces | 2. | While the components soak in the solution, use a soft wipe or soft brush to remove all visible residue from all exterior surfaces of the components. |
| | 3. | While the product is soaking in the cleaning solution, clean it for a minimum of one minute using a soft, lint-free wipe or a soft bristled brush. |

Complete stage II: Rinsing*Proceed as follows:*

2 rinses of at least 1 minute each

1. Fully immerse all parts in tap water (< 45 °C / 113 °F) and thoroughly rinse all accessible surfaces in two rinses of at least **one** minute each.

Use fresh water for each rinse.

Allow residual water to drip off for a sufficient length of time.

Complete stage III: Disinfection*Proceed as follows:*

Disinfect for at least 12 min

1. Fully immerse all parts in disinfectant solution for at least **twelve** minutes.
2. Remove all adherent air bubbles from the component surfaces.

All accessible surfaces must remain immersed in the disinfectant solution bath throughout the entire disinfecting time.

Complete stage IV: Rinsing*Proceed as follows:*

2 rinses of at least 1 min each

1. Fully immerse the product in a tap water bath (less than 45 °C / 113 °F) and thoroughly rinse all accessible surfaces twice for a minimum of **one** minute each.

Use fresh water for each rinse.

Allow residual water to drip off for a sufficient length of time.

Complete stage V: Final rinsing*Proceed as follows:*

At least 1 min in deionized water

1. Fully immerse all parts in deionized water (< 45 °C / 113 °F) and thoroughly rinse all accessible surfaces for a minimum of **one** minute.

Allow residual water to drip off for a sufficient length of time.

Complete stage VI: Drying*Proceed as follows:*

Visual inspection

1. Dry all parts with a clean, soft, lint-free wipe/ operating room towel or using medical-quality filtered compressed air (pmax = 0.5 bar).
2. After drying, visually inspect the product (see [section 3.4.1](#)) in a well-lit area and make sure it is dry, undamaged and free of visible residue.

Use a magnifying glass. Repeat the cleaning and disinfection process if necessary.

Immediately set aside any damaged products.

- This concludes the manual cleaning and disinfection process.

4.7 Automated Cleaning and Thermal Disinfection

4.7.1 Overview

Stage	Step	Temperature (°C/°F)	Time (min)	Cleaning solution / water quality
I	Pre-rinsing	Cold	2	Tap water
II	Cleaning	60/140	10	Alkaline cleaning solution
III	Intermediate rinsing	per the equipment manufacturer's standard cycle	1	Tap water
IV	Intermediate rinsing	per the equipment manufacturer's standard cycle	1	Tap water
V	Thermal disinfection	90/194	5	Deionized or purified water A ₀ value: > 3000
VI	Drying	99/210	30	---

Table 4-2: Overview of automated cleaning and thermal disinfection.

Note: Automated reprocessing can cause color-anodized or plastic components (e.g. serial rings and the eyepiece cap) to fade.

4.7.2 Performance

Performing automated cleaning and thermal disinfection

AUTOCLAVE

Thermal disinfection should be performed only for products that are labeled as autoclavable.

Use a washer/disinfector that has been validated as effective and meets the performance requirements of ISO 15883-1 (or the respective country-specific version thereof).

Dry removal of the products is recommended for automated cleaning.

In the case of wet removal, do not use foaming detergents, and rinse the products thoroughly prior to the automated cleaning.

Proceed as follows:

1. Place all parts in a reprocessing tray.
2. Place the reprocessing tray into the cleaning device.
3. Start the cleaning cycle described in [section 4.7.1](#) in accordance with the manufacturer's instructions and instruction for use for the washer/disinfector.

CAUTION! Risk of scalding when unloading the washer/disinfector. Be sure to wear gloves.

4. Remove the reprocessing tray from the washer.

- Visual inspection
5. After drying, visually inspect the product ([section 3.4.1](#)) in a well-lit area and make sure it is dry, undamaged and free of visible residue.

Use a magnifying glass. Repeat the cleaning and disinfection process if necessary.

Immediately set aside any damaged products.

- This concludes the automated cleaning and disinfection process.

4.8 Sterilization

4.8.1 Steam Sterilization

Performing steam sterilization

AUTOCLAVE

Steam sterilization should be performed only for products that are labeled as autoclavable.

CAUTION! Perform the sterilization only on a properly cleaned and dried product. Residue from cleaning agents and organic material can negatively influence the sterilization result.

Proceed as follows:

1. Place all parts in a reprocessing tray.
2. Pack the tray with two single layers of polypropylene wrap.

CAUTION! Risk of scalding when loading the sterilizer. Wear suitable gloves.

3. Load the wrapped tray in the sterilizer.

Load the sterilizer as per the validated loading plan.

Follow the manufacturer's instructions!

4. Start the sterilization cycle in accordance with the manufacturer's instructions and instruction for use for the sterilizer.

For steam sterilization in a fractionated vacuum autoclave cycle use the following validated parameters:

Cycle type	Fractionated vacuum processes
Pulses	4
Temperature	134 °C (273 °F)
Holding time	at least 3 min (effective sterilization time)

Note: It is the responsibility of the operator to establish that the sterilizer used meets the parameters above.

CAUTION! Risk of scalding when unloading the sterilizer. Wear suitable gloves.

NOTICE! Products can be shock sensitive when hot. Avoid knocking or shaking.

NOTICE! Damage resulting from a sudden change in temperature. Allow products to cool to room temperature; do not use additional cooling measures.

5. Remove the sterilized products from the sterilizer.

Ensure that the products remain sterile after reprocessing.

- This concludes sterilization.

5 Maintenance and Repairs

5.1 Troubleshooting Table

Issue	Possible causes	Fixes
Image cloudy	Glass surfaces soiled	Manual cleaning as per section 4.6 , then re-process, check the water quality
	Leaky, defective lens system	Send the endoscope in for repair
Image too dark, too little illumination	Glass surfaces soiled	Manual cleaning as per section 4.6 , then re-process, check the water quality
	Unsuitable light guide	Use a suitable light guide
	Light guide incorrectly attached to endoscope	Check that the light guide is seated properly
	Defective fiber optics	Check fiber optics as per section 3.4
	Defective light guide or light source	Check the light guide and light source
Image is yellowed	Soiled fiber optics	Manual cleaning as per section 4.6 , then re-process, check the water quality
	Light guide soiled or defective	Check the light guide (e.g. by illuminating a white surface)
Corrosion, staining, discoloration	Inadequate cleaning (e.g. remaining protein residue)	Manual cleaning as per section 4.6 , if necessary, rub thoroughly, then re-process
	Insufficient rinsing between the reprocessing stages, especially prior to sterilization	Rinse sufficiently between the reprocessing stages
	Excessive chloride concentration	Check water quality
	Excessive concentration of minerals (e.g. limescale) or organic substances	Check the water quality, only use demineralized water if necessary
	Heavy metal ions and/or silicates, elevated concentration of iron, copper, and manganese in the water or sterilizer steam	Check the water quality, only use demineralized water if necessary

Issue	Possible causes	Fixes
	Contaminated, too frequently reused cleaning and disinfecting solution	The cleaning and disinfecting solution should be replaced regularly
	Extraneous rust, e.g. resulting from rust contained in steam or damaged or non-corrosion-resistant instruments reprocessed at the same time	Check the supply systems, if multiple items are being reprocessed together, pay attention to material compatibility and signs of prior damage; prevent the different items from touching one another
	Contact corrosion	Avoid contact with other devices

Table 5-1: Troubleshooting.

5.2 Repair

Should you need to arrange a repair for the device, please contact our Service Department.

When you return the device, please also enclose a detailed **description of the fault**.

WARNING

Contaminated device

Risk of infection

- > Reprocess the product properly prior to shipping

Only ship devices that have been properly reprocessed.

The device should be returned in its original packaging, if possible.

Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

6 Product Data

6.1 Technical Data

Item no.	STS-2040-302ASA
Direction of view	0°
Field of view	Standard
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-302AA
Direction of view	0°
Field of view	Wide angle
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-302BA
Direction of view	30°
Field of view	Wide angle
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-302BSA
Direction of view	30°
Field of view	Standard
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-302ESA
Direction of view	12°
Field of view	Standard
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-302CA
Direction of view	70°
Field of view	Wide angle
Working length	302 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2029-302AA
Direction of view	0°
Field of view	Wide angle
Working length	302 mm
Outer diameter of insertion portion	2.9 mm

Item no.	STS-2029-302BA
Direction of view	30°
Field of view	Wide angle
Working length	302 mm
Outer diameter of insertion portion	2.9 mm

Item no.	STS-2040-298DSA-RWF
Direction of view	5°
Field of view	Standard
Working length	298 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-298ASA-RWF
Direction of view	0°
Field of view	Standard
Working length	298 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-300GSA-RWF
Direction of view	25°
Field of view	Standard
Working length	300 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-303CA-RWF
Direction of view	70°
Field of view	Wide angle
Working length	303 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-280ASA-OLY
Direction of view	0°
Field of view	Standard
Working length	280 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-280ESA-OLY
Direction of view	12°
Field of view	Standard
Working length	280 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-283BSA-OLY
Direction of view	30°
Field of view	Standard
Working length	283 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-283BA-OLY
Direction of view	30°
Field of view	Wide angle
Working length	283 mm
Outer diameter of insertion portion	4.0 mm

Item no.	STS-2040-283CA-OLY
Direction of view	70°
Field of view	Wide angle
Working length	283 mm
Outer diameter of insertion portion	4.0 mm

6.2 Ambient Conditions

Transport and storage conditions

Temperature	-20 °C to +70 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

Reprocessed devices must be protected against recontamination, in a dark, dry, well-ventilated, dust-free and light-controlled location with a consistent temperature.

Direct sunlight, high temperatures, high humidity or radiation can damage the device or present a risk of infection.

When placing into storage, make sure that the device cannot be damaged by other instruments. It is therefore best to store the device individually or use containers in which it can be secured.

Operating conditions

Temperature	+15 °C to +40 °C
Relative air humidity	10 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

6.3 Spare Parts and Accessories

Use original spare parts and accessories only.

Image	Designation	Item number
	Light guide adapter for Storz systems	STS-ACC-0000-LGUA
	Light guide adapter for Wolf systems	STS-ACC-0000-LGUA-RWF

Table 6-1: Spare parts and accessories.

7 Disposal

WARNING

Contaminated device

Risk of infection

- > The device must be reprocessed prior to disposal

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

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