strauss[®]
S U R G I C A L

Instructions For Use

Ureterorenoscopes
Endoscopes for Ureterorenoscopy

USA

TPA702-025-01_B
Version: B
2019 - Dec - 18

	Content	Page
0	Important Information about this Document.....	6
0.1	Scope of Validity, Identification, Purpose.....	6
0.2	Target Group.....	6
0.3	Using and Storing this Document.....	6
1	General Information about the Device.....	7
1.1	Scope of Delivery.....	7
1.2	Device Description.....	7
1.2.1	Performance Characteristics and Function.....	7
1.2.2	Visual Overview.....	8
1.2.3	Maximum Compatible Instruments Diameter.....	9
1.2.4	Recommended Light Guides.....	9
1.3	Usage.....	10
1.3.1	Intended Use.....	10
1.3.2	Indications.....	10
1.3.3	Contraindications.....	10
1.4	Conformity.....	11
1.4.1	Standards and Directives.....	11
1.4.2	Medical Device Classification.....	11
1.5	Marking.....	12
1.5.1	Pictograms and Information on the Device and Packaging.....	12
1.5.2	Pictograms in this Document.....	13
1.6	Service Department Contact Details.....	13
2	General Safety Information.....	14
2.1	Safety Messages in this Document.....	14
2.1.1	Safety Messages at the Start of a Chapter.....	14
2.1.2	Safety Messages in the Body of the Text.....	14
2.2	Product Safety.....	15
2.2.1	Basic Safety Instructions.....	15
2.2.2	Staff Qualifications.....	15
3	Operation.....	16
3.1	Safety Instructions.....	16
3.2	Staff Qualifications.....	17
3.3	Technical Inspection Prior to Use.....	18
3.3.1	Visual Inspection.....	18
3.3.2	Function Check.....	19
3.4	Assembly.....	19
3.5	Connecting a Light Guide.....	21
3.6	Post-Use and Pre-Cleaning Activities.....	21

4	Reprocessing	22
4.1	Safety Instructions	22
4.2	Staff Qualifications	23
4.3	Validated Procedures	23
4.4	The Process Flow	24
4.5	Cleaning Agents	24
4.6	Manual Cleaning with Immersion	25
4.6.1	Overview	25
4.6.2	Performance	25
4.7	Automated Cleaning	27
4.7.1	Manual Pre-Cleaning	27
4.7.2	Overview	28
4.7.3	Performance	29
4.8	Steam Sterilization	30
5	Maintenance and Repairs	31
5.1	Troubleshooting Table	31
5.2	Repair	32
6	Product Data	33
6.1	Technical Data	33
6.1.1	Ureterorenoscopes	33
6.2	Ambient Conditions	33
6.3	Spare Parts and Accessories	34
7	Disposal	35

0 Important Information about this Document

0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the following devices:

- Ureterorenoscopes with the item numbers STS-6075-425HA, STS-6075-310HA, STS-6060-425HA, STS-6060-310HA
- Adapters with the item numbers STS-ACC-6000-EMSA-LC1, STS-ACC-6000-EMSA-LC2
- Instrument bridge with the item number STS-URO-6000-BRGA-002
- Three-slot membranes with the item number STS-ACC-6000-MEMD-000
- Sealing caps with the item number STS-ACC-5000-CAPD-000
- Reprocessing tray with the item number STS-URO-6000-BRGA-002

These instructions for use are an integral component of the device set and contain all the information required by users and operators for safe and proper use.

These instructions for use do not describe endoscopic procedures or techniques during surgery.

0.2 Target Group

These instructions for use are intended for physicians, medical assistants, and employees of central sterile supply who are entrusted with the operation, handling, and reprocessing of the device.

0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

1 General Information about the Device

1.1 Scope of Delivery

Check scope of delivery! Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

1.2 Device Description

1.2.1 Performance Characteristics and Function

The ureterorenoscopes described herein are semiflexible endoscopes for visualizing the surgical field during surgical procedures within the urinary tract and kidneys.

Structure and function A semiflexible endoscope consists of a fiber optic cable and sensitive image transmission system with eyepiece.

The fiber optic cable is used to illuminate the site inside the body. The connector for connecting the light guide to the light source is situated at the proximal end of the endoscope. The adapters required to connect the light guide are included in the scope of delivery.

The distal end of the endoscope features an object lens that captures the image from inside the body. The image is sent through the image transmission system to the eyepiece.

The eyepiece cup connects the eyepiece to an endocoupler that transmits the image to a camera. The camera's controller converts the signal so that the image can be displayed on a monitor.

Use with an instrument bridge The endoscope is intended to be used in conjunction with an instrument bridge that allows the insertion of instruments, probes and irrigation fluid during the procedure.

Additional instruments Additional instruments are required for endoscopic procedures.

It is the surgeon's responsibility to make sure the instruments required for a particular endoscopic procedure are assembled and available for use.

Note: During this medical application, adhere to national statutory regulations and directives. For the treatment of urolithiasis, for example, you can refer to the European Association of Urology (EAU) **Guidelines on Urolithiasis**.

1.2.2 Visual Overview

Ureterorenoscope


Figure 1-1: Ureterorenoscope and instrument bridge.

- 1** Ureterorenoscope
- 2** Instrument bridge
- 3** Adapter for connecting Storz light guides
- 4** Adapter for connecting Wolf light guides
- 5** Light guide connector (ACMI)
- 6** Eyepiece

Instrument bridge


Figure 1-2: Instrument bridge.

- 1** Instrument bridge coupling piece
- 2** Stopcock nuts
- 3** Stopcocks for irrigation channels and angled instrument channel
- 4** Release lever
- 5** Membrane port
- 6** Three-slot membrane
- 7** Instrument receptacle (standard)
- 8** Sealing cap

Adapters for EMS Swiss Lithoclast


Figure 1-3: Adapter for EMS Swiss Lithoclast 1.

- 1** Instrument receptacle (adapter)
- 2** Adapter for EMS Swiss Lithoclast 1


Figure 1-4: Adapter for EMS Swiss Lithoclast 2.

- 1** Locking screw
- 2** Adapter for EMS Swiss Lithoclast 2

1.2.3 Maximum Compatible Instruments Diameter

The endoscopes are available with working channels of various diameters. The following table specifies the maximum compatible instrument diameter for each of the endoscopes.

Item number	Maximum instrument diameter
STS-6075-425HA	1 x 5 Fr or 2 x 3 Fr
STS-6075-310HA	1 x 5 Fr or 2 x 3 Fr
STS-6060-425HA	1 x 4 Fr or 2 x 2.4 Fr
STS-6060-310HA	1 x 4 Fr or 2 x 2.4 Fr

1.2.4 Recommended Light Guides

The specifications for compatible light guides are:

- Fiber bundles of 3.5 to 4.8 mm Ø
- Length up to 300 cm

1.3 Usage

1.3.1 Intended Use

The URS are intended to be used for the visualization of anatomic structures in the urethra, bladder, ureter and renal pelvis during minimally invasive interventions by access through the natural body orifice of the urethra. The URS is intended to allow the insertion of working instruments through the integrated working channel.

The instrument bridge is intended to allow the insertion of working instruments into the integrated working channel of the URS and to allow the connection of irrigation equipment and the channeling of irrigation liquids. Furthermore adapters for EMS lithotriptors can be mounted on the instrument bridge.

The membrane and the sealing cap are intended to block the outflow of liquids from inside the body to maintain the pressure inside the anatomic structures.

The device is intended to be used in combination with medical electrical devices that are, **as a minimum**, compliant with the requirements for applied parts classified as type **BF** (that is parts offering increased protection against electric shock), according to IEC 60601-1.


This isolation barrier can be realized by the medical electrical devices themselves or by the connecting cables to the endoscope and must be assured for every connection between the endoscope and connected devices.

Devices or connecting cables that meet the requirements for applied parts classified as type BF feature the adjacent symbol.

Type CF is the most stringent classification for devices and connecting cables considered applied parts within the meaning of IEC 60601-1. They bear the adjacent heart symbol and are approved for direct conductive contact with the heart.

Defibrillation-proof devices or connecting cables are also compatible. These will be marked with one of the adjacent symbols depending on the degree of protection against electric shock offered.

1.3.2 Indications

The ureterorenoscopes are indicated for the transurethral examination and treatment of the upper urinary tract, including the ureter and kidney, and, using additional accessoires, to perform various diagnostic and therapeutic procedures.

The indications for an endoscopic application depend on the medical condition of the patient and the individual risk/benefit analysis by the surgeon.

1.3.3 Contraindications

Contraindications may exist depending on the patient's general condition or specific symptoms.

The decision to perform an endoscopic procedure remains with the responsible surgeon and must be made on the basis of an individual risk/benefit analysis.

1.4 Conformity

1.4.1 Standards and Directives

The device is compliant with the requirements of the following applicable standard:

- **IEC 60601-2-18** Medical electrical equipment - Part 2-18:
Particular requirements for basic safety and essential performance of endoscopic equipment

The device meets the requirements of the following applicable EU directive:

- Medical Devices Directive **93/42/EEC**

1.4.2 Medical Device Classification

The device is a Class IIa medical device within the meaning of Directive 93/42/EEC.

1.5 Marking

1.5.1 Pictograms and Information on the Device and Packaging

This section describes the different pictograms that feature on the device or packaging.

	The device is suitable for autoclave sterilization
	Item number
	Serial number
	Lot number
	Quantity of products of a particular item number
	Manufacturer
	Date of manufacture
	Caution (IEC 60601-1 3rd edition) / Take note of accompanying documents (IEC 60601-1 2nd edition)
	Non-sterile, sterilize prior to each use
	MR unsafe
	Permissible storage temperature
	Permissible relative air humidity during storage
	Permissible atmospheric pressure during storage
	Do not reuse
	Caution: Federal (US) law restricts this device to sale by or on order of a physician

1.5.2 Pictograms in this Document

This section describes the pictograms used in this document.

General warning sign

Biohazard warning sign, risk of infection

Hot surfaces warning sign

Laser beam warning sign

Type BF applied part according to IEC 60601-1

Type CF applied part according to IEC 60601-1

Defibrillation-proof type BF applied part

Defibrillation-proof type CF applied part

1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

2 General Safety Information

2.1 Safety Messages in this Document

2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to material damage:

NOTICE

Nature and origin of the risk of material damage

Potential consequences of non-observance

> Preventive action

2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

NOTICE! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

2.2 Product Safety

2.2.1 Basic Safety Instructions

Our products are developed and manufactured to the highest quality standards.

Risks – despite the high quality

Although this device corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully and observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

The original packaging should be retained and reused when returning the device for servicing.

Follow the instructions for use for all devices and instruments that are to be used in conjunction with this device.

WARNING! This device is supplied in a non-sterile condition. Risk of infection. Carry out reprocessing of the device and its accessories prior to first use.

WARNING! Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.

WARNING! Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.

WARNING! High-intensity light source. Risk of injury to eyes. Do not look directly into the open end of the light guide or light exit of the endoscope.

WARNING! Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.

CAUTION! Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.

WARNING! Improper handling, maintenance, and use carries a risk to the patient and user, or can lead to premature wearing of the device.

CAUTION! When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. This is particularly important if a type CF applied part endoscope is used, in which case a type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

2.2.2 Staff Qualifications

Observe qualification requirements

Specific qualifications are required for operation and reprocessing of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

3 Operation

3.1 Safety Instructions

WARNING

Interactions between devices in simultaneous use (e.g. lasers, lithotriptors)

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the type BF classification requirements according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used

WARNING

Displaying a recorded image instead of the live image or an altered image orientation

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera
- > Ensure that the live image is displayed in the correct orientation (not mirrored)

WARNING

Use of non-sterile parts

Infection risk for the patient

- > Use only properly reprocessed endoscopes and endoscopic accessories
- > Accessories supplied in a non-sterile condition (e.g. three slot membrane and sealing cap) must be subject to reprocessing before use
- > Always carry out a visual inspection prior to use
- > Anchor the light guide in the operating field and secure against slipping, allowing some slack for maneuvering

WARNING

Recontamination resulting from improper handling

Infection risk for the patient

- > Adhere to the infection control practices that are in place

⚠ CAUTION

High temperatures resulting from light source usage

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

- > Use only light guides that are suitable for use with the endoscope
- > Avoid the prolonged use of intense light
- > You should always use the lowest possible light output setting that will allow you to illuminate the target area
- > Do not touch the light source in the vicinity of the lamp
- > Do not touch light guide connectors
- > Do not allow the distal end of the endoscope to come into contact with patient tissue or with combustible or heat-sensitive materials
- > Do not exceed the maximum permitted ambient temperature for the light source

⚠ CAUTION

Coupling laser beams in the endoscope

Damage to the eyes caused by looking directly into the eyepiece

- > Wear laser safety glasses

NOTICE

Improper handling

Risk of damage to the product

- > Set down carefully and do not bend or knock
- > Avoid lever forces
- > Avoid scratching the surface
- > Do not activate the laser until the tip of the laser fiber can be seen through the endoscope
- > Take care when removing instruments from the surgical field

3.2 Staff Qualifications

The device may be operated only by physicians and medical assistants who have received instruction in how to use the device and who have the requisite training or advanced training, knowledge, and practical experience in the endoscopic procedure to be applied as defined by the provisions in place at the site of operation.

3.3 Technical Inspection Prior to Use

3.3.1 Visual Inspection

Performing a visual inspection

*A visual inspection must be performed **prior to each use**.*

Proceed as follows:

- | | |
|---|---|
| Is there any external damage? | <p>1. Check the endoscope and all the components to be used for external damage.</p> <p>Do not use the device if it has any sharp edges or corners, bulges, or rough surfaces that might cause injury to the patient.</p> |
| Are the fiber optics in working order? | <p>2. Check that the fiber optics in the endoscope are in full working order.</p> <p>To do this, hold the distal end toward a bright light (not a cold light source) and the light guide connection toward your eyes. Gently move the endoscope from side to side, taking note of the brightness of the fibers.</p> <p>If more than 20 % of the fibers remain dark, the endoscope must not be used. Submit the endoscope to repair.</p> |
| Are the surfaces clean and smooth? | <p>3. Check that the glass surfaces and fiber optic end faces of the endoscope are clean and smooth.</p> <p>Do not use the endoscope if the surfaces are soiled or scratched. Set aside the endoscope and follow the steps described in <u>section 5.1</u>.</p> |
| Is the image clear, bright, and complete? | <p>4. Make sure you can see a clear, bright, complete image.</p> <p>Look through the eyepiece and assess the image quality.</p> <p>Do not use the endoscope if the image is yellowed, dark, speckled, or cropped. Set aside the endoscope and follow the steps described in <u>section 5.1</u>.</p> |
| Is the OR equipment mutually compatible? | <p>5. Check that all the equipment to be used for the procedure is mutually compatible.</p> <p>Only use the OR equipment if the isolation barrier for protecting against electric shock meets at minimum the type BF classification requirements according to IEC 60601-1 (see <u>Intended Use</u>).</p> |
| | <p><input checked="" type="checkbox"/> This concludes the inspection.</p> |

3.3.2 Function Check

Performing a function check

A function check must be performed **prior to each medical procedure**.

Proceed as follows:

1. Connect a camera and a monitor.
2. Switch on all the system components that you intend using for the procedure.

Switching on the system
3. Direct the camera head at an object in the room and focus the image.

Focusing the image

Check that you are able to achieve a sharp, bright, good quality image.

Do not use the system if the image contains vertical lines or color variations, flickers, or you are unable to achieve a sharp, bright, good quality image.

This concludes the function check.

3.4 Assembly

Assembling the instrument bridge

Proceed as follows:

CAUTION! Non-sterile, single-use three-slot membranes and sealing caps. Infection risk for the patient. Sterilize the three-slot membranes and sealing caps prior to use and do not reuse.

1. Prepare the instrument bridge.

Insert the stopcock plugs and fasten them using the stopcock nuts. Check the stopcocks for smooth movement and bring them into an open position.

Insert a three-slot membrane into the membrane port.

Attach the instrument receptacle (standard) to the instrument bridge.

Plug a sealing cap onto the instrument receptacle (standard).


Figure 3-1: Preparation of the instrument bridge.

2. Lock the instrument bridge on the ureterorenoscope.

Note: To hermetically seal, the coupling piece must be clean and free of grease - do not grease.

Align the instrument bridge coupling piece and the endoscope until the arrows are facing each other.

Press the coupling piece into the receptacle on the endoscope body, until you feel it click into place.

Make sure that the instrument bridge is firmly connected to the endoscope.

This concludes the assembly of the instrument bridge.

Assembling the adapter for EMS Swiss Lithoclast 1*Proceed as follows:*Screw on only
lightly

1. Replacing the instrument receptacle (standard), attach the instrument receptacle (adapter) to the instrument bridge.
2. Lightly screw on the adapter to the instrument receptacle (adapter).

Note: Do not over-tighten the adapter, otherwise it may be difficult to remove later.


Figure 3-2: Screw the adapter to the instrument bridge.

- This concludes the assembly of the adapter.

Assembling the adapter for EMS Swiss Lithoclast 2*Proceed as follows:*

1. Completely open the locking screw on the adapter.
2. Fully slide the adapter onto the instrument bridge so the cut-outs on the sides (underneath the locking screw) fit around the stopcock shafts on the instrument bridge.

Lock the adapter in place using the locking screw, which needs to fit into the indentation on the instrument bridge.


Figure 3-3: Indentation on the instrument bridge and fitted adapter.

- This concludes the assembly of the adapter.

3.5 Connecting a Light Guide

Connecting the light guide

Different light guides require different adapters to connect to the endoscope.

Proceed as follows:

1. Screw the suitable adapter to the light guide connector on the endoscope and, if applicable, to the light guide.
 2. Connect the light guide to the light guide connector.
- The light guide is connected.

3.6 Post-Use and Pre-Cleaning Activities

Post-use cleaning

*Complete the pre-cleaning process **immediately after use**.*

Proceed as follows:

1. Uncouple the endoscope from the endocoupler and detach the light guide.
2. Remove the instrument bridge from the endoscope and detach the adapters from the light guide connection.
3. Break down the individual components of the instrument bridge and dispose of the single-use items sealing cap and three slot membrane. Also remove the instrument holder and, if applicable, lithoclast adapter from the instrument bridge.

NOTICE! The eyepiece cap is non-detachable.

4. Remove visible surgical residue as completely as possible using a lint-free wipe that has been moistened with an enzymatic cleaning solution.
 5. Rinse the lumens of all components with tap water. Use a disposable 20 ml syringe.
 6. Dry the device using a soft, lint-free wipe.
 7. Place all the parts into a dry removal container and seal.
 8. Arrange for reprocessing, ensuring that all of the device components are reprocessed within 6 hours.
- This concludes the post-use cleaning.

4 Reprocessing

4.1 Safety Instructions

⚠ WARNING

If the product is potentially contaminated with Creutzfeldt-Jakob pathogens, it is not suitable for reprocessing

Risk of cross contamination during use and reprocessing

- > Do not reprocess products that are potentially contaminated
- > Products that are potentially contaminated must be disposed of

⚠ CAUTION

Improper cleaning

Risk to patient as a result of inadequate cleaning, damage to the device

- > Use an automatic washer that meets the performance requirements of ISO 15883-1
- > Load the washer with care: all contents must be rinsed and cleaned fully (no rinsing 'blind spots')
- > Connect instruments with lumens and channels directly to the dedicated ports on the reprocessing tray
- > Remove all stopcocks and the instrument receptacle from the instrument bridge
- > Ensure that the washer is properly maintained
- > Only use cleaning agents that are approved for the device
- > The devices must be pre-cleaned immediately after use and reprocessed within 6 hours
- > When pre-cleaning, do not use fixing temperatures in excess of 45 °C
- > When pre-cleaning, do not use fixing cleaning agents (containing aldehyde or alcohol)

NOTICE

Contact with chloride-based solutions

Corrosion and destruction of the device

- > Avoid contact with chloride-based solutions, which can be found, for example, in surgical residue, lotions, drugs, saline solutions, and cleaning agents
- > If the products should come into contact with chloride-based solutions, rinse thoroughly with deionized water and dry completely

NOTICE

Use of ultrasonic bath

Damage to the endoscope

- > Do not expose endoscope to ultrasonic bath

NOTICE**Unsuitable cleaning agents and process chemicals**

Corrosion damage, premature aging, and visible material deterioration

- > Only use cleaning agents that are approved for the device
- > Only use process chemicals that are compatible with the device's materials according to the chemical manufacturer's recommendations
- > Follow all the chemical manufacturer's application specifications regarding temperature, concentration, and contact time
- > Do not use process chemicals that can cause stress cracking or brittleness of plastic materials

4.2 Staff Qualifications

In many countries, the qualifications required by personnel responsible for the reprocessing of medical devices are regulated by law.

In any event, the reprocessing of medical devices must only be carried out by qualified personnel who have the necessary knowledge and expertise.

Professional expertise may be obtained through specialized further education or through training and practical experience, and enhanced through appropriate further training if necessary.

4.3 Validated Procedures

The effectiveness of the following procedures:

- Manual cleaning
- Automated cleaning
- Sterilization

as described in this document, has been fully validated.

Operator responsibility It is the operator's responsibility to introduce, document, implement, and maintain a validated reprocessing procedure.

Make sure that the equipment used for reprocessing is properly maintained.

4.4 The Process Flow

Pre-cleaning The product has to be prepared for reprocessing immediately after use by pre-cleaning.

Refer to [section 3.6](#) for pre-cleaning instructions.

Work steps The reprocessing procedure described in this document consists of the following steps:

- Pre-cleaning immediately after use
- Manual or automated cleaning
- Sterilization

Protective clothing During reprocessing, staff should wear personal protective clothing.

Careful cleaning The product must be cleaned thoroughly at the beginning of reprocessing. It is essential that the sterilization medium reaches all parts of the product.

Automated cleaning and steam sterilization The optimum and most reliable reprocessing results are achieved using automated cleaning with subsequent steam sterilization using the fractionated vacuum steam procedure.

Legislation and standards Adhere to national statutory regulations, national and international standards and directives, and the infection control practices that are in place at your institution for reprocessing.

Reprocessing prior to each use New and returned products following repair have to undergo complete reprocessing before first use.

The alternating use of different reprocessing methods can cause premature aging of the product.

Excessive doses of cleaning agents can damage the product and cause laser engravings to fade.

www.a-k-i.org Please see www.a-k-i.org for more detailed information about hygienically safe, value-preserving reprocessing methods that also protect the materials.

4.5 Cleaning Agents

The following products can be used for **manual cleaning**:

- Cidezyme/Enzol (Johnson & Johnson)

The following products can be used for **automated cleaning**:

- neodisher MediClean forte 0.5 %
(Chem. Fabrik Dr. Weigert GmbH & Co. KG)

Wherever possible, use the above cleaning agents only.

Follow the manufacturer's instructions! Before use, read the manufacturer's user information carefully and follow the specifications regarding concentration, temperature, usage time, water quantities, and contact time.

4.6 Manual Cleaning with Immersion

4.6.1 Overview

Stage	Work Step	Temperature (°C/°F)	Time (min)	Water Quality	Cleaning Solution
I	Cleaning	35-45 / 95-113	5	Tap water	Enzymatic cleaning solution
II	Rinsing	< 45 / 113	2x ≥1	Tap water	---
III	Final rinsing	< 45 / 113	1x ≥1	High purity water	---
IV	Drying	---	---	---	---

Table 4-1: Overview of manual cleaning with immersion.

4.6.2 Performance

Complete stage I: Cleaning

Prepare the cleaning solution according to the manufacturer's instructions.

Proceed as follows:

NOTICE! Damage resulting from rough handling. Handle the product with care. Do not knock or bend. Set down carefully.

NOTICE! Scratch-sensitive surfaces. Risk of corrosion due to abrasion. Do not use metal brushes, metal objects, or abrasive cleaners to clean the instruments. To remove soiling, use pads and cleaning solution only.

- | | |
|---------------------------------|---|
| Immerse for 5 minutes | 1. Completely immerse all parts of the endoscope and the instrument bridge (and, if applicable, the adapter for lithoclast) in the cleaning solution for five minutes.

Make sure that all accessible surfaces remain immersed in the cleaning solution throughout the entire cleaning time.

Make sure that cavities are free of air bubbles. |
| Wipe or brush external surfaces | 2. Then, while they remain completely immersed in the cleaning solution, clean the exterior surfaces of all parts using a soft, clean, lint-free wipe or a soft brush until no more residue can be seen. |
| Brush cavities | 3. While they are completely immersed in the cleaning solution, thoroughly clean the cavities of all parts. Brush until no more residue can be seen, but for at least one minute.

During the cleaning, move mobile components in all directions as far as they will go several times. |
| | 4. After brushing, thoroughly rinse the exterior surfaces and lumens with cleaning solution at least five times. Use a disposable 20 ml syringe. |

Complete stage II: Rinsing

Proceed as follows:

- | | |
|--|---|
| 2 rinses in tap water of at least 1 min each | 1. Completely immerse all parts in a water bath (< 45 °C / 113 °F) and thoroughly rinse all accessible surfaces twice for a minimum of one minute each.

Use fresh water for each rinse. |
|--|---|

- Rinse cavities 2. Rinse the lumens of all parts at least **five** times with 50 ml tap water each time. Use a disposable 50 ml syringe.

Allow residual water to drip off for a sufficient length of time.

Complete stage III: Final rinsing

Proceed as follows:

- 1 minute in high purity water 1. Completely immerse the product in a high purity water bath (< 45 °C / 113 °F) and thoroughly rinse all accessible surfaces for a minimum of **one** minute.

- Rinse cavities 2. Rinse the lumens of all parts at least **five** times with 50 ml high purity water. Use a disposable 50 ml syringe.

During the rinsing, maneuver mobile parts in all directions as far as they will go for at least **three** times.

Allow residual water to drip off for a sufficient length of time.

Complete stage IV: Drying

Proceed as follows:

1. Dry all parts with a clean, soft, lint-free wipe and use medical-quality filtered compressed air (max = 0,5 bar) to dry cavities.

- Visual inspection 2. After drying, visually inspect the components in a well-lit area (see [section 3.3.1](#)) and make sure they are dry, undamaged and free of visible residue.

Use a magnifying glass. Repeat the cleaning process if necessary. Immediately set aside any damaged products.

3. As needed, apply lubricant to stopcocks for better performance (e.g. PARALIQ GB 363, Klüber Lubrication).

Lubricant requirements: steam sterilization according to ISO 17665-1 and ANSI/AAMI ST8, microbial efficiency according to ISO 17664 and AAMI TIR 12, biocompatibility testing according to ISO 10993-1.

- This concludes the manual cleaning process.

4.7 Automated Cleaning

4.7.1 Manual Pre-Cleaning

Performing manual pre-cleaning

Perform manual pre-cleaning immediately before automated cleaning.

Prepare the cleaning solution according to the manufacturer's instructions.

Proceed as follows:

NOTICE! Damage resulting from rough handling. Handle the device with care. Do not knock or bend. Set down carefully.

NOTICE! Scratch-sensitive surfaces. Risk of corrosion due to abrasion. Do not use metal brushes, metal objects, or abrasive cleaners to clean the instruments. To remove soiling, use pads and cleaning solution only.

Immerse for
5 minutes

1. Completely immerse all parts of the endoscope and the instrument bridge (and, if applicable, the lithoclast adapter) in the cleaning solution for **five** minutes.

Make sure that all accessible surfaces remain immersed in the cleaning solution throughout the entire cleaning time.

Make sure that cavities are free of air bubbles.

2. Then brush all parts of the product until they are free of visible soiling (at least **one** minute). In doing so, pay particular attention to the contact points between the endoscope and the instrument bridge.
3. Brush the lumens of the instrument bridge with a suitable brush until they are free of visible soiling, but for at least **one** minute.

This concludes the manual pre-cleaning process.

4.7.2 Overview

Stage	Work step	Temperature (°C/°F)	Time (min)	Water quality	Notes
	Manual Pre-Cleaning				as per <u>section 4.7.1</u>
I	Pre-rinsing	< 25/77	2	Tap water	---
II	Cleaning	55/131	10	Tap water	Alkaline detergent pH > 10, e.g. neodisher MediClean forte 2 ml - 10 ml per liter
III	Rinsing I	> 10/50	1	Tap water	---
IV	Rinsing II	> 10/50	1	Deionized water	---
V	Final rinsing	90/194	5	Deionized water	A ₀ value: > 3000
VI	Drying	---	---	---	standard program of the washer

Table 4-2: Overview automated cleaning.

Note: Color anodized parts or plastic components (e.g. standard rings and the eyepiece cup) may bleach out during automated reprocessing.

4.7.3 Performance

Performing automated cleaning

Use a washer that has been validated as effective and meets the performance requirements of ISO 15883-1 (or the respective country-specific version thereof).

Dry removal of the products is recommended prior to automated cleaning.

In the case of wet removal, do not use foaming detergents and rinse the products thoroughly prior to the automated cleaning.

Proceed as follows:

1. Make sure that manual pre-cleaning was completed.
2. Place all components in the reprocessing tray and the small parts in the integrated small parts tray. Make sure to include the light guide adapters and, where applicable, the lithoclast adapter as well.

NOTICE! The reprocessing tray can hold only **one** ureterorenoscope including instrument bridge.

3. Connect the parts with lumens to the rinsing ports of the reprocessing tray.
4. Place the reprocessing tray into the washer.
5. Attach the reprocessing tray to the washer.
6. Start the cleaning cycle described in section 4.7.2 in accordance with the manufacturer's instructions and instruction for use for the washer.

CAUTION! Risk of scalding when unloading the washer. Be sure to wear gloves.

7. Remove the reprocessing tray from the washer.

Visual inspection

8. After drying, visually inspect the product in a well-lit area (see section 3.3.1) and make sure it is undamaged and free of visible residue.
Use a magnifying glass for the visual inspection. Repeat the cleaning process if necessary.
Immediately set aside any damaged products.
9. If necessary, apply lubricant to stopcocks for better performance (e.g. PARALIQ GB 363, Klüber Lubrication).

Lubricant requirements: steam sterilization according to ISO 17665-1 and ANSI/AAMI ST8, microbial efficiency according to ISO 17664 and AAMI TIR 12, biocompatibility testing according to ISO 10993-1.

- This concludes the automated cleaning process.

4.8 Steam Sterilization

Performing steam sterilization

AUTOCLAVE

Steam sterilization should be performed only for products that are labeled as autoclavable.

CAUTION! Residue from cleaning agents and organic material can negatively influence the sterilization result. Perform the sterilization only on a properly cleaned and dried product.

CAUTION! The biological indicators, sterilization wraps, pouches, and containers used must be FDA-cleared for sterilization under the parameters given.

Proceed as follows:

1. Load all parts into the reprocessing tray and small parts like the light guide adapter, instrument receptacle, three-slot membrane or sealing cap into the integrated small component tray.

CAUTION! Non-sterile, single-use three-slot membrane and sealing cap. Infection risk for the patient. Sterilize the three-slot membrane and sealing cap prior to use and do not reuse.

NOTICE! The reprocessing tray is designed to accommodate one ureterorenoscope including instrument bridge.

Make sure that the product is fully disconnected from the tube connectors on the reprocessing tray.

2. Double wrap the tray with two single layers of FDA-cleared sterilization wrap using a sequential double wrapping technique.

CAUTION! Risk of scalding when loading the sterilizer. Wear suitable gloves.

3. Load the wrapped tray in the sterilizer.
4. Start the sterilization cycle in accordance with the manufacturer's instructions and instruction for use for the sterilizer.

For steam sterilization in a fractionated vacuum autoclave cycle use the following validated parameters:

Cycle type	Fractionated vacuum process
Pulses	4
Temperature	270° F (132° C)
Holding time	at least 4 minutes (effective sterilization time)
Drying time	at least 30 minutes

Note: It is the responsibility of the operator to establish that the sterilizer used meets the parameters above.

CAUTION! Risk of scalding when unloading the sterilizer. Wear suitable gloves.

NOTICE! Products can be shock sensitive when hot. Avoid knocking or shaking.

NOTICE! Damage resulting from a sudden change in temperature. Allow products to cool to room temperature; do not use additional cooling measures.

5. Remove the sterilized products from the sterilizer.

Ensure that the products remain sterile after reprocessing.

This concludes steam sterilization.

5 Maintenance and Repairs

5.1 Troubleshooting Table

Issue	Possible Causes	Fixes
Image cloudy	Glass surfaces soiled	Manual cleaning (section 4.6), then reprocess, check the water quality
	Leaky, defective lens system	Send the endoscope in for repair
Image too dark, too little illumination	Glass surfaces soiled	Manual cleaning (section 4.6), then reprocess, check the water quality
	Unsuitable light guide	Use a suitable light guide
	Light guide incorrectly attached to endoscope	Check that the light guide is seated properly
	Defective fiber optics	Check fiber optics (section 3.3)
	Defective light guide or light source	Check the light guide and light source
Image is yellowed	Soiled fiber optics	Manual cleaning (section 4.6), then reprocess, check the water quality
	Light guide soiled or defective	Check the light guide (e.g. by illuminating a white surface)
Corrosion, staining, discoloration	Inadequate cleaning (e.g. remaining protein residue)	Manual cleaning (section 4.6), rub thoroughly if necessary, then reprocess
	Insufficient rinsing between the reprocessing stages, especially prior to sterilization	Rinse sufficiently between the reprocessing stages
	Excessive chloride concentration	Check water quality
	Excessive concentration of minerals (e.g. limescale) or organic substances	Check the water quality, only use demineralized water if necessary
	Heavy metal ions and/or silicates, elevated concentration of iron, copper, and manganese in the water or sterilizer steam	Check the water quality, only use demineralized water if necessary
	Contaminated, too frequently reused cleaning solution	The cleaning solution should be replaced regularly

Issue	Possible Causes	Fixes
	Extraneous rust, e.g. resulting from rust contained in steam or damaged or non-corrosion-resistant instruments reprocessed at the same time	If multiple items are being reprocessed together, check the supply systems for material compatibility and signs of prior damage; prevent the different items from touching one another
	Contact corrosion	Avoid contact with other components

Table 5-1: Troubleshooting.

5.2 Repair

Should you need to arrange a repair for the device, please contact our Service Department.

When you return the device, please also enclose a detailed **description of the fault**.

WARNING

Contaminated device

Risk of infection

> Reprocess the device properly prior to shipping

Only ship devices that have been properly reprocessed.

The device should be returned in its original packaging, if possible.

Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return these to sender.

6 Product Data

6.1 Technical Data

6.1.1 Ureterorenoscopes

Item no.	STS-6075-425HA
Direction of view	8°
Field of view	85°
Working length	425 mm
Outer diameter of insertion portion (distal end / maximum diameter)	7.5 / 9.8 Fr

Item no.	STS-6075-310HA
Direction of view	8°
Field of view	85°
Working length	310 mm
Outer diameter of insertion portion (distal end / maximum diameter)	7.5 / 9.8 Fr

Item no.	STS-6060-425HA
Direction of view	8°
Field of view	100°
Working length	425 mm
Outer diameter of insertion portion (distal end / maximum diameter)	6.0 / 7.5 Fr

Item no.	STS-6060-310HA
Direction of view	8°
Field of view	100°
Working length	310 mm
Outer diameter of insertion portion (distal end / maximum diameter)	6.0 / 7.5 Fr

6.2 Ambient Conditions

Transport and storage conditions

Temperature	-20 °C to +70 °C
Rel. air humidity	5 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

Reprocessed devices must be protected against recontamination, in a dark, dry, well-ventilated, dust-free and light-controlled location with a consistent temperature.

Direct sunlight, high temperatures, high humidity or radiation can damage the device or present a risk of infection.

When placing into storage, make sure that the device cannot be damaged by other instruments. It is therefore best to store the device individually or use containers in which it can be secured.

Operating conditions

Temperature	+15 °C to +32 °C
Rel. air humidity	30 % to 85 %
Atmospheric pressure	70 kPa to 106 kPa

6.3 Spare Parts and Accessories

Use original spare parts and accessories only.

Photo	Designation	Item Number
	Sealing caps, unsterile single-use item, pack of 10	STS-ACC-5000-CAPD-000
	Stopcock plugs, steel, pack of 12 (4x clockwise, 8x counterclockwise)	STS-ACC-6000-SCPA-000
	Lubricant PARALIQ GB 363, 40 g	RFP00200
	Three-slot membrane, unsterile single-use item, pack of 10	STS-ACC-6000-MEMD-000
	Light guide adapter for Storz systems	STS-ACC-0000-LUGA
	Light guide adapter for Wolf systems	STS-ACC-0000-LUGA-RWF
	Adapter for EMS Swiss Lithoclast 1	STS-ACC-6000-EMSA-LC1
	Adapter for EMS Swiss Lithoclast 2	STS-ACC-6000-EMSA-LC2
	Instrument bridge	STS-URO-6000-BRGA-002
	Reprocessing tray for URS	STS-ACC-5600-TRYA-000
	Cleaning brush for working channel (URS with an outer diameter of insertion portion on the distal end of 7.5 Fr.)	STS-ACC-6075-520N-BRU
	Cleaning brush for working channel (URS with an outer diameter of insertion portion on the distal end of 6 Fr.)	STS-ACC-6060-520N-BRU
	Cleaning brush for instrument bridge irrigation channels	STS-ACC-5030-100N-BRU
	Cleaning brush for instrument bridge working channel	STS-ACC-5043-100N-BRU

Table 6-1: Spare parts and accessories.

7 Disposal

WARNING

Contaminated device

Risk of infection

- > The device must be reprocessed prior to disposal

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

STRAUSS SURGICAL

3020 NW 82 Ave.
Miami, FL 33122

Phone: 305 436 0599

Fax: 305 436 0399

Web: www.straussurgical.com

strauss[®]
S U R G I C A L

Instructions For Use

Ureterorenoscopes
Endoscopes for Ureterorenoscopy
(Not valid for USA)

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Ureterorenoscopios

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Ureterorenoscópios

TPA768-025-34_A
Version: A
2019 - Dec - 18

	Content	Page
0	Important Information about this Document	10
0.1	Scope of Validity, Identification, Purpose	10
0.2	Target Group	10
0.3	Using and Storing this Document	10
1	General Information about the Device	11
1.1	Scope of Delivery	11
1.2	Product Description	11
1.2.1	Performance Characteristics and Function	11
1.2.2	Visual Overview	12
1.2.3	Maximum Compatible Instruments Diameter	13
1.2.4	Recommended Light Guides	13
1.3	Usage	14
1.3.1	Intended Use	14
1.3.2	Indications	14
1.3.3	Contraindications	14
1.4	Conformity	15
1.4.1	Standards and Directives	15
1.4.2	Medical Device Classification	15
1.5	Marking	16
1.5.1	Pictograms and Information on the Device and Packaging	16
1.5.2	Pictograms in this Document	17
1.6	Service Department Contact Details	17
2	General Safety Information	18
2.1	Safety Messages in this Document	18
2.1.1	Safety Messages at the Start of a Chapter	18
2.1.2	Safety Messages in the Body of the Text	18
2.2	Product Safety	19
2.2.1	Basic Safety Notices	19
2.2.2	Staff Qualifications	19
3	Operation	20
3.1	Safety Notices	20
3.2	Staff Qualifications	21
3.3	Technical Inspection Prior to Use	22
3.3.1	Visual Inspection	22
3.3.2	Function Check	23
3.4	Assembly	23
3.5	Connection of a Light Guide	25
3.6	Post-Use and Pre-Cleaning Activities	25
4	Reprocessing	26
4.1	Safety Notices	26
4.2	Staff Qualifications	27
4.3	Validated Procedures	27
4.4	The Process Flow	28
4.5	Cleaning Agents	28
4.6	Manual Cleaning and Disinfection	29
4.6.1	Overview	29
4.6.2	Performance	29
4.7	Automated Cleaning and Thermal Disinfection	31
4.7.1	Manual Pre-Cleaning	31
4.7.2	Overview	32
4.7.3	Performance	33
4.8	Steam Sterilization	34
5	Maintenance and Repairs	35
5.1	Troubleshooting Table	35
5.2	Repair	36

6	Product Data	37
6.1	Technical Data.....	37
6.1.1	Ureterorenoscopes.....	37
6.2	Ambient Conditions	37
6.3	Spare Parts and Accessories.....	38
7	Disposal.....	39

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6.1.1 Ureterorenoscópios 102

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7 Descarte 104

0 Important Information about this Document

0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the following products:

- Ureterorenoscopes with the item numbers STS-6075-425HA, STS-6075-310HA, STS-6060-425HA, STS-6060-310HA
- Adapters with the item numbers STS-ACC-6000-EMSA-LC1, STS-ACC-6000-EMSA-LC2
- Instrument bridge with the item number STS-URO-6000-BRGA-002
- Three-slot membranes with the item number STS-ACC-6000-MEMD-000
- Sealing caps with the item number STS-ACC-5000-CAPD-000
- Reprocessing tray with the item number STS-URO-6000-BRGA-002

These instructions for use are an integral component of the product set and contain all the information required by users and operators for safe and proper use.

These instructions for use do not describe endoscopic procedures or techniques during surgery.

0.2 Target Group

These instructions for use are intended for physicians, medical assistants, and employees of central sterile supply who are entrusted with the operation, handling, and reprocessing of the device.

0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

1 General Information about the Device

1.1 Scope of Delivery

Check scope of delivery! Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

1.2 Product Description

1.2.1 Performance Characteristics and Function

The ureterorenoscopes described herein are semiflexible endoscopes for visualizing the surgical field during surgical procedures within the urinary tract and kidneys.

Structure and function A semiflexible endoscope consists of a fiber optic cable and sensitive image transmission system with eyepiece.

The fiber optic cable is used to illuminate the site inside the body. The connector for connecting the light guide to the light source is situated at the proximal end of the endoscope. The adapters required to connect the light guide are included in the scope of delivery.

The distal end of the endoscope features an object lens that captures the image from inside the body. The image is sent through the image transmission system to the eyepiece.

The eyepiece cap connects the eyepiece to an endocoupler that transmits the image to a camera. The camera's controller converts the signal so that the image can be displayed on a monitor.

Use with an instrument bridge The endoscope is intended to be used in conjunction with an instrument bridge that allows the insertion of instruments, probes and irrigation fluid during the procedure.

Additional instruments Additional instruments are required for endoscopic procedures.

It is the surgeon's responsibility to make sure the instruments required for a particular endoscopic procedure are assembled and available for use.

Note: During this medical application, adhere to national statutory regulations and directives. For the treatment of urolithiasis, for example, you can refer to the European Association of Urology (EAU) **Guidelines on Urolithiasis**.

1.2.2 Visual Overview

Ureterorenoscope


Figure 1-1: Ureterorenoscope and instrument bridge.

- 1** Ureterorenoscope
- 2** Instrument bridge
- 3** Adapter for connecting Storz light guides
- 4** Adapter for connecting Wolf light guides
- 5** Light guide connector (ACMI)
- 6** Eyepiece

Instrument bridge


Figure 1-2: Instrument bridge.

- 1** Instrument bridge coupling piece
- 2** Stopcock nuts
- 3** Stopcocks for irrigation channels and angled instrument channel
- 4** Release lever
- 5** Membrane port
- 6** Three-slot membrane
- 7** Instrument receptacle (standard)
- 8** Sealing cap

Adapters for EMS Swiss Lithoclast


Figure 1-3: Adapter for EMS Swiss Lithoclast 1.

- 1** Instrument receptacle (adapter)
- 2** Adapter for EMS Swiss Lithoclast 1


Figure 1-4: Adapter for EMS Swiss Lithoclast 2.

- 1** Locking screw
- 2** Adapter for EMS Swiss Lithoclast 2

1.2.3 Maximum Compatible Instruments Diameter

The endoscopes are available with working channels of various diameters. The following table specifies the maximum compatible instrument diameter for each of the endoscopes.

Item number	Maximum instrument diameter
STS-6075-425HA	1 x 5 Fr or 2 x 3 Fr
STS-6075-310HA	1 x 5 Fr or 2 x 3 Fr
STS-6060-425HA	1 x 4 Fr or 2 x 2.4 Fr
STS-6060-310HA	1 x 4 Fr or 2 x 2.4 Fr

1.2.4 Recommended Light Guides

The specifications for compatible light guides are:

- Active diameter of 3.5 to 4.8 mm
- Length up to 300 cm

1.3 Usage

1.3.1 Intended Use

The URS are intended to be used for the visualization of anatomic structures in the urethra, bladder, ureter and renal pelvis during minimally invasive interventions by access through the natural body orifice of the urethra. The URS is intended to allow the insertion of working instruments through the integrated working channel.

The instrument bridge is intended to allow the insertion of working instruments into the integrated working channel of the URS and to allow the connection of irrigation equipment and the channeling of irrigation liquids. Furthermore adapters for EMS lithotriptors can be mounted on the instrument bridge.

The membrane and the sealing cap are intended to block the outflow of liquids from inside the body to maintain the pressure inside the anatomic structures.

The device is intended to be used in combination with medical electrical devices that meet at **minimum** the type **BF** classification requirements for applied parts according to IEC 60601-1 for elevated protection from electric shock.


This insulation barrier can be realized by the medical electrical devices themselves or by the connection cables to the endoscope and must be assured for every connection between the endoscope and connected devices.

Devices or connection cables that meet the type **BF** classification requirements feature the adjacent symbol.

Devices and connection cables that meet the requirements for applied parts classified as type **CF** offer an even higher level of protection. They bear the adjacent heart symbol and are approved for direct conductive contact with the heart.

Defibrillation-proof devices or connection cables are also compatible. These will be marked with one of the adjacent symbols depending on the degree of protection against electric shock offered.

1.3.2 Indications

The indications for an endoscopic application depend on the medical condition of the patient and the individual risk/benefit analysis by the surgeon.

1.3.3 Contraindications

Contraindications may exist depending on the patient's general condition or specific symptoms.

The decision to perform an endoscopic procedure remains with the responsible surgeon and must be made on the basis of an individual risk/benefit analysis.

1.4 Conformity

1.4.1 Standards and Directives

The device meets the requirements of the following applicable standard:

- **IEC 60601-2-18** Medical electrical equipment - Part 2-18:
Particular requirements for basic safety and essential performance of endoscopic equipment

The device meets the requirements of the following applicable EU directive:

- Medical Devices Directive **93/42/EEC**

1.4.2 Medical Device Classification

The device is a Class IIa medical device within the meaning of Directive 93/42/EEC.

1.5 Marking

1.5.1 Pictograms and Information on the Device and Packaging

This section describes the different pictograms that feature on the product or packaging.

	The device is suitable for autoclave sterilization
	Item number
	Serial number
	Lot number
	Quantity of products of a particular item number
	Manufacturer
	Date of manufacture
	Caution (IEC 60601-1 3rd edition) / Take note of accompanying documents (IEC 60601-1 2nd edition)
	Non-sterile, sterilize prior to each use
	MR unsafe
	Permissible storage and transport temperature
	Permissible relative air humidity during storage and transport
	Permissible atmospheric pressure during storage and transport
	Do not reuse
	Caution: Federal (US) law restricts this device to sale by or on order of a physician

1.5.2 Pictograms in this Document

This section describes the pictograms used in this document.

General warning sign

Biohazard warning sign, risk of infection

Hot surfaces warning sign

Laser beam warning sign

Type BF applied part according to IEC 60601-1

Type CF applied part according to IEC 60601-1

Defibrillation-proof type BF applied part

Defibrillation-proof type CF applied part

1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

2 General Safety Information

2.1 Safety Messages in this Document

2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION

Nature and origin of the injury risk

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to material damage:

NOTICE

Nature and origin of the risk of material damage

Potential consequences of non-observance

> Preventive action

2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

WARNING! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

CAUTION! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

NOTICE! Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

2.2 Product Safety

2.2.1 Basic Safety Notices

Our products are developed and manufactured to the highest quality standards.

Risks despite highest of quality standards

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

The original packaging should be retained and reused when returning the device for servicing.

Follow the instructions for use for all devices and instruments that are to be used in conjunction with this device.

WARNING! This device is supplied in a non-sterile condition. Risk of infection. Carry out reprocessing of the device and its accessories prior to first use.

WARNING! Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.

WARNING! Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.

WARNING! High-intensity light source. Risk of injury to eyes. Do not look directly into the open end of the light guide or light exit of the endoscope.

WARNING! Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.

CAUTION! Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.

WARNING! Improper handling, maintenance, and use carries a risk to the patient and user, or can lead to premature wearing of the device.

WARNING! When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. This is particularly important if a type CF applied part endoscope is used, in which case a type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

2.2.2 Staff Qualifications

Observe qualification requirements

Specific qualifications are required for operation and reprocessing of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

3 Operation

3.1 Safety Notices

WARNING

Interactions between devices in simultaneous use (e.g. lasers, lithotriptors)

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the requirements classified as type BF according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used

WARNING

Displaying a recorded image instead of the live image or an altered image orientation

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera
- > Ensure that the live image is displayed in the correct orientation (not mirrored)

WARNING

Use of non-sterile parts

Infection risk for the patient

- > Use only properly reprocessed endoscopes and endoscopic accessories
- > Accessories supplied in a non-sterile condition (e.g. three-slot membrane and sealing cap) must be subject to reprocessing before use
- > Always carry out a visual inspection prior to use
- > Anchor the light guide in the surgical field and secure against slipping, allowing some slack for maneuvering

WARNING

Recontamination resulting from improper handling

Infection risk for the patient

- > Adhere to the infection control practices that are in place

**⚠ CAUTION****High temperatures resulting from light source usage**

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

- > Use only light guides that are suitable for use with the endoscope
- > Avoid the prolonged use of intense light
- > You should always use the lowest possible light output setting that will allow you to illuminate the target area
- > Do not touch the light source in the vicinity of the lamp
- > Do not touch light guide connections
- > Do not allow the distal end of the endoscope to come into contact with patient tissue or with combustible or heat-sensitive materials
- > Do not exceed the maximum permitted ambient temperature for the light source

**⚠ CAUTION****Coupling laser beams in the endoscope**

Damage to the eyes caused by looking directly into the eyepiece

- > Wear laser safety glasses

NOTICE**Improper handling**

Risk of damage to the product

- > Set down carefully without bending or knocking
- > Avoid lever forces
- > Avoid scratching the surface
- > Do not activate the laser until the tip of the laser fiber can be seen through the endoscope
- > Take care when removing instruments from the surgical field

3.2 Staff Qualifications

The device may be operated only by physicians and medical assistants who have received instruction in how to use the device and who have the requisite training or advanced training, knowledge, and practical experience in the endoscopic procedure to be applied as defined by the provisions in place at the site of operation.

3.3 Technical Inspection Prior to Use

3.3.1 Visual Inspection

Performing a visual inspection

*A visual inspection must be performed **prior to each use**.*

Proceed as follows:

- | | |
|---|---|
| Is there any external damage? | 1. Check the endoscope and all the components to be used for external damage.

Do not use the product if it has any sharp edges or corners, bulges, or rough surfaces that might cause injury to the patient. |
| Are the fiber optics in working order? | 2. Check that the fiber optics in the endoscope are in full working order.

To do this, hold the distal end toward a bright light (not a cold light source) and the light guide connection toward your eyes. Gently move the endoscope from side to side, taking note of the brightness of the fibers.

If more than 20 % of the fibers remain dark, the device must not be used. Submit the endoscope to repair. |
| Are the surfaces clean and smooth? | 3. Check that the glass surfaces and fiber optic end faces of the endoscope are clean and smooth.

Do not use the endoscope if the surfaces are soiled or scratched. Set aside the endoscope and follow the steps described in <u>section 5.1</u> . |
| Is the image clear, bright, and complete? | 4. Make sure you can see a clear, bright, complete image.

Look through the eyepiece and assess the image quality.

Do not use the endoscope if the image is yellowed, dark, speckled, or cropped. Set aside the endoscope and follow the steps described in <u>section 5.1</u> . |
| Is the OR equipment mutually compatible? | 5. Check that all the equipment to be used for the procedure is mutually compatible.

Only use the OR equipment if the isolation barrier for protecting against electric shock meets at minimum the type BF classification requirements according to IEC 60601-1 (see <u>Intended Use</u>). |
| | <input checked="" type="checkbox"/> This concludes the inspection. |

3.3.2 Function Check

Performing a function check

A function check must be performed **prior to each medical procedure**.

Proceed as follows:

1. Connect a camera and a monitor.
 2. Switch on all the system components that you intend using for the procedure.
 3. Direct the camera head at an object in the room and focus the image.
- Check that you are able to achieve a sharp, bright, good quality image.
- Do not use the system if the image contains vertical lines or color variations, flickers, or you are unable to achieve a sharp, bright, good quality image.
- This concludes the function check.

3.4 Assembly

Assembling the instrument bridge

Proceed as follows:

CAUTION! Non-sterile, single-use three-slot membranes and sealing caps. Infection risk for the patient. Sterilize the three-slot membranes and sealing caps prior to use and do not reuse.

1. Prepare the instrument bridge.

Insert the stopcock plugs and fasten them using the stopcock nuts. Check the stopcocks for smooth movement and bring them into an open position.

Insert a three-slot membrane into the membrane port.

Attach the instrument receptacle (standard) to the instrument bridge.

Plug a sealing cap onto the instrument receptacle (standard).


Figure 3-1: Preparation of the instrument bridge.

2. Lock the instrument bridge on the ureterorenoscope.

Note: To hermetically seal, the coupling piece must be clean and free of grease - do not grease.

Align the instrument bridge coupling piece and the endoscope until the arrows are facing each other.

Press the coupling piece into the receptacle on the endoscope body, until you feel it click into place.

Make sure that the instrument bridge is firmly connected to the endoscope.

- This concludes the assembly of the instrument bridge.

Assembling the adapter for EMS Swiss Lithoclast 1

Proceed as follows:

Screw on only
lightly

1. Replacing the instrument receptacle (standard), attach the instrument receptacle (adapter) to the instrument bridge.
2. Lightly screw on the adapter to the instrument receptacle (adapter).

Note: Do not over-tighten the adapter, otherwise it may be difficult to remove later.


Figure 3-2: Screw the adapter to the instrument bridge.

- This concludes the assembly of the adapter.

Assembling the adapter for EMS Swiss Lithoclast 2

Proceed as follows:

1. Completely open the locking screw on the adapter.
2. Fully slide the adapter onto the instrument bridge so the cut-outs on the sides (underneath the locking screw) fit around the stopcock shafts on the instrument bridge.

Lock the adapter in place using the locking screw, which needs to fit into the indentation on the instrument bridge.


Figure 3-3: Indentation on the instrument bridge and fitted adapter.

- This concludes the assembly of the adapter.

3.5 Connection of a Light Guide

Connecting the light guide

Different light guides require different adapters to connect to the endoscope.

Proceed as follows:

1. Screw the suitable adapter to the light guide connector on the endoscope and, if applicable, to the light guide.
 2. Connect the light guide to the light guide connector.
- The light guide is connected.

3.6 Post-Use and Pre-Cleaning Activities

Post-use cleaning

*Complete the pre-cleaning process **immediately after use**.*

Proceed as follows:

1. Uncouple the endoscope from the endocoupler and detach the light guide.

NOTICE! The eyepiece cap is non-detachable.

2. Remove the instrument bridge from the endoscope and detach the adapters from the light guide connection.
3. Break down the individual components of the instrument bridge and dispose of the single-use items sealing cap and three-slot membrane. Remove the instrument receptacle and, if applicable, adapter for lithoclast from the instrument bridge.
4. Remove visible surgical residue as completely as possible using a lint-free wipe that has been moistened with an enzymatic cleaning solution.
5. Rinse the lumens of all components with tap water. Use a disposable 20 ml syringe.
6. Dry the product using a soft, lint-free wipe.
7. Place all the parts into a dry removal container and seal.
8. Arrange for reprocessing, ensuring that all of the product components are reprocessed within 6 hours.

This concludes the post-use cleaning.

4 Reprocessing

4.1 Safety Notices

WARNING

If the product is potentially contaminated with Creutzfeldt-Jakob pathogens, it is not suitable for reprocessing

Risk of cross contamination during use and reprocessing

- > Do not reprocess products that are potentially contaminated
- > Products that are potentially contaminated must be disposed of

WARNING

Contact of bladder cancer patients with Cidex OPA

Possible anaphylactic reaction

- > Do not bring patients with a history of bladder cancer in contact with devices that were disinfected using Cidex OPA
- > Use automated cleaning and thermal disinfection, if applicable

CAUTION

Improper cleaning and disinfection

Risk to patient as a result of inadequate cleaning and disinfection, damage to the device

- > Use an automatic washer/disinfector that meets the performance requirements of ISO 15883-1
- > Load the washer/disinfector with care: all contents must be rinsed and cleaned fully (no rinsing 'blind spots')
- > Connect instruments with lumens and channels directly to the dedicated ports on the reprocessing tray
- > Remove all stopcocks and the instrument receptacle from the instrument bridge
- > Ensure that the washer/disinfector is properly maintained
- > Only use cleaning and disinfecting agents that are approved for the product
- > The devices must be pre-cleaned immediately after use and reprocessed within 6 hours
- > When pre-cleaning, do not use fixing temperatures in excess of 45 °C
- > When pre-cleaning, do not use fixing cleaning and disinfecting agents (active ingredients: Aldehyde, alcohol)

NOTICE

Contact with chloride-based solutions

Corrosion and destruction of the device

- > Avoid contact with chloride-based solutions, which can be found, for example, in surgical residue, lotions, drugs, saline solutions, and cleaning and disinfecting agents
- > If the products should come into contact with chloride-based solutions, rinse thoroughly with DI water and dry completely

NOTICE**Use of ultrasonic bath**

Damage to the endoscope

- > Do not expose endoscope to ultrasonic bath

NOTICE**Unsuitable cleaning agents and process chemicals**

Corrosion damage, premature aging, and visible material deterioration

- > Only use cleaning and disinfecting agents that are approved for the product
- > Only use process chemicals that are compatible with the device's materials according to the chemical manufacturer's recommendations
- > Follow all the chemical manufacturer's application specifications regarding temperature, concentration, and contact time
- > Do not use process chemicals that can cause stress cracking or brittleness of plastic materials

4.2 Staff Qualifications

In many countries, the qualifications required by personnel responsible for the reprocessing of medical devices are regulated by law.

In any event, the reprocessing of medical devices must always fall under the responsibility of qualified personnel who have the necessary knowledge and expertise.

This knowledge and expertise can be acquired by completing further training in this field or as a result of completing a dedicated qualification with practical experience, supplemented by appropriate further training measures, where required.

4.3 Validated Procedures

The effectiveness of the following procedures:

- Manual cleaning and disinfection
- Automated cleaning and thermal disinfection
- Sterilization

as described in this document, has been fully validated.

Operator responsibility It is the operator's responsibility to introduce, document, implement, and maintain a validated reprocessing procedure.

Make sure that the equipment used for reprocessing is properly maintained.

4.4 The Process Flow

Pre-cleaning The product has to be prepared for reprocessing immediately after use by pre-cleaning.

Refer to [section 3.6](#) for pre-cleaning instructions.

Work steps The reprocessing procedure described in this document consists of the following steps:

- Pre-cleaning immediately after use
- Cleaning and disinfection (manual or automated)
- Sterilization

Protective clothing During reprocessing, staff should wear personal protective clothing.

Careful cleaning The product must be cleaned thoroughly at the beginning of reprocessing. It is essential that the sterilization medium reaches all parts of the product.

Automated cleaning and steam sterilization The optimum and most reliable reprocessing results are achieved using automated cleaning and disinfection with subsequent steam sterilization using the pre-vacuum steam procedure.

Legislation and standards Adhere to national statutory regulations, national and international standards and directives, and the infection control practices that are in place at your institution for reprocessing.

Reprocessing prior to each use New and returned products following repair have to undergo complete reprocessing before first use.

The alternating use of different reprocessing methods can cause premature aging of the product.

Excessive doses of cleaning agents can damage the product and cause laser engravings to fade.

www.a-k-i.org Please see www.a-k-i.org for more detailed information about hygienically safe, value-preserving reprocessing methods that also protect the materials.

4.5 Cleaning Agents

The following products can be used for **manual cleaning and disinfection**:

- Cidezyme/Enzol (Johnson & Johnson)
- Cidex OPA (Johnson & Johnson)

The following products can be used for **automated cleaning**:

- neodisher MediClean forte 0.5 %
(Chem. Fabrik Dr. Weigert GmbH & Co. KG)

Wherever possible, use the above cleaning and disinfecting agents only.

Follow the manufacturer's instructions! Before use, read the manufacturer's user information carefully and follow the specifications regarding concentration, temperature, usage time, water quantities, and contact time.

4.6 Manual Cleaning and Disinfection

4.6.1 Overview

Stage	Work step	Temperature (°C/°F)	Time (min)	Water quality	Cleaning solution
I	Cleaning	35-45 / 95-113	5	Tap water	Enzymatic cleaning solution
II	Rinsing	< 45 / 113	3x ≥1	Deionized water	---
III	Disinfection	as per manufacturer's instructions	12	---	Disinfectant solution (undiluted)
IV	Final rinsing	< 45 / 113	3x ≥2	Deionized water	---
V	Drying	---	---	---	---

Table 4-1: Overview of manual cleaning and disinfection.

4.6.2 Performance

Complete stage I: Cleaning

Prepare the cleaning solution according to the manufacturer's instructions.

Proceed as follows:

NOTICE! Damage resulting from rough handling. Handle the product with care. Do not knock or bend. Set down carefully.

NOTICE! Scratch-sensitive surfaces. Risk of corrosion due to abrasion. Do not use metal brushes, metal objects, or abrasive cleaners to clean the instruments. To remove soiling, use pads and cleaning solution only.

- | | | |
|---------------------------------|----|--|
| Immerse for 5 minutes | 1. | Completely immerse all parts of the endoscope and the instrument bridge (and, if applicable, the adapter for lithoclast) in the cleaning solution for five minutes.

Make sure that all accessible surfaces remain immersed in the cleaning solution throughout the entire cleaning time.

Make sure that cavities are free of air bubbles. |
| Wipe or brush external surfaces | 2. | Then, while they remain completely immersed in the cleaning solution, clean the exterior surfaces of all parts using a soft, clean, lint-free wipe or a soft brush until no more residue can be seen. |
| Brush cavities | 3. | While they are completely immersed in the cleaning solution, thoroughly clean the cavities of all parts. Brush until no more residue can be seen, but for at least one minute.

During the cleaning, move mobile components in all directions as far as they will go several times. |
| | 4. | After brushing, thoroughly rinse the exterior surfaces and lumens with cleaning solution at least five times. Use a disposable 20 ml syringe. |

Complete stage II: Rinsing*Proceed as follows:*

- 3 rinses in tap water of at least 1 min each
1. Completely immerse all parts in a water bath (< 45 °C / 113 °F) and thoroughly rinse all accessible surfaces three times for a minimum of **one** minute each.
Use fresh, deionized water for each rinse.
- Rinse cavities
2. Rinse the lumens of all parts at least **five** times with 50 ml deionized water each time. Use a disposable 50 ml syringe.
Allow residual water to drip off for a sufficient length of time.

Complete stage III: Disinfection*Proceed as follows:*

- Disinfect for at least 12 minutes
1. Fully immerse all parts in disinfectant solution for at least **twelve** minutes.
 2. Remove all adherent air bubbles from the component surfaces.
 3. Make sure that the disinfectant solution reaches all mobile parts by maneuvering them in all directions as far as they will go for at least **three** times.
 4. Rinse the lumens of all parts at least **five** times thoroughly with disinfectant solution. Use a disposable 50 ml syringe.
All accessible surfaces must remain immersed in the disinfectant solution bath throughout the entire disinfecting time.

Complete stage IV: Final rinsing*Proceed as follows:*

- 2 minutes in deionized water
1. Completely immerse the product in a deionized water bath (< 45 °C / 113 °F) and thoroughly rinse all accessible surfaces for a minimum of **two** minutes. Repeat this work step **three** times.
- Rinse cavities
2. Rinse the lumens of all parts at least **five** times with 50 ml deionized water. Use a disposable 50 ml syringe.
During the rinsing, maneuver mobile parts in all directions as far as they will go for at least **three** times.
Allow residual water to drip off for a sufficient length of time.

Complete stage V: Drying*Proceed as follows:*

1. Dry all parts with a clean, soft, lint-free wipe and use medical-quality filtered compressed air (max = 0,5 bar) to dry cavities.
- Visual inspection
2. After drying, visually inspect the components in a well-lit area (see [section 3.3.1](#)) and make sure they are dry, undamaged and free of visible residue.
Use a magnifying glass. Repeat the cleaning process if necessary. Immediately set aside any damaged products.
 3. As needed, apply lubricant to stopcocks for better performance (e.g. PARALIQ GB 363, Klüber Lubrication).

Lubricant requirements: steam sterilization according to ISO 17665-1 and ANSI/AAMI ST8, microbial efficiency according to ISO 17664 and AAMI TIR 12, biocompatibility testing according to ISO 10993-1.

- This concludes the manual cleaning process.

4.7 Automated Cleaning and Thermal Disinfection

4.7.1 Manual Pre-Cleaning

Performing manual pre-cleaning

Perform manual pre-cleaning immediately before automated cleaning.

Prepare the cleaning solution according to the manufacturer's instructions.

Proceed as follows:

NOTICE! Damage resulting from rough handling. Handle the product with care. Do not knock or bend. Set down carefully.

NOTICE! Scratch-sensitive surfaces. Risk of corrosion due to abrasion. Do not use metal brushes, metal objects, or abrasive cleaners to clean the instruments. To remove soiling, use pads and cleaning solution only.

Immerse for
5 minutes

1. Completely immerse all parts of the endoscope and the instrument bridge (and, if applicable, the adapter for lithoclast) in the cleaning solution for **five** minutes.

Make sure that all accessible surfaces remain immersed in the cleaning solution throughout the entire cleaning time.

Make sure that cavities are free of air bubbles.

2. Then brush all parts of the product until they are free of visible soiling (at least **one** minute). In doing so, pay particular attention to the contact points between the endoscope and the instrument bridge.
3. Brush the lumens of the instrument bridge with a suitable brush until they are free of visible soiling, but for at least **one** minute.

- This concludes the manual pre-cleaning process.

4.7.2 Overview

Stage	Work step	Temperature (°C/°F)	Time (min)	Water quality	Notes
	Manual Pre-Cleaning				as per <u>section 4.7.1</u>
I	Pre-rinsing	< 25/77	2	Tap water	---
II	Cleaning	55/131	10	Tap water	Alkaline detergent pH > 10, e.g. neodisher MediClean forte 2 ml - 10 ml per liter
III	Rinsing I	> 10/50	1	Tap water	---
IV	Rinsing II	> 10/50	1	Deionized water	---
V	Thermal disinfection	90/194	5	Deionized water	A ₀ value > 3000
VI	Drying	---	---	---	standard program of the washer

Table 4-2: Overview automated cleaning.

Note: Automated reprocessing can cause color-anodized or plastic components (e.g. serial rings and the eyepiece cap) to fade.

4.7.3 Performance

Performing automated cleaning

Use a washer/disinfector that has been validated as effective and meets the performance requirements of ISO 15883-1 (or the respective country-specific version thereof).

Dry removal of the products is recommended prior to automated cleaning.

In the case of wet removal, do not use foaming detergents and rinse the products thoroughly prior to the automated cleaning.

Proceed as follows:

1. Make sure that manual pre-cleaning was completed.
2. Place all components in the reprocessing tray and the small parts in the integrated small parts tray. Make sure to include the light guide adapters and, where applicable, the lithoclast adapter as well.

NOTICE! The reprocessing tray can hold only **one** ureterorenoscope including instrument bridge.

3. Connect the parts with lumens to the rinsing ports of the reprocessing tray.
4. Place the reprocessing tray into the washer/disinfector.
5. Attach the reprocessing tray to the washer/disinfector.
6. Start the cleaning cycle described in [section 4.7.2](#) in accordance with the manufacturer's instructions and instruction for use for the washer/disinfector.

CAUTION! Risk of scalding when unloading the washer/disinfector. Be sure to wear gloves.

7. Remove the reprocessing tray from the washer/disinfector.

Visual inspection

8. After drying, visually inspect the product in a well-lit area (see [section 3.3.1](#)) and make sure it is undamaged and free of visible residue.

Use a magnifying glass for the visual inspection. Repeat the cleaning process if necessary.

Immediately set aside any damaged products.

9. If necessary, apply lubricant to stopcocks for better performance (e.g. PARALIQ GB 363, Klüber Lubrication).

Lubricant requirements: steam sterilization according to ISO 17665-1 and ANSI/AAMI ST8, microbial efficiency according to ISO 17664 and AAMI TIR 12, biocompatibility testing according to ISO 10993-1.

- This concludes the automated cleaning process.

4.8 Steam Sterilization

Performing steam sterilization

AUTOCLAVE

Steam sterilization should be performed only for products that are labeled as autoclavable.

CAUTION! Residue from cleaning agents and organic material can negatively influence the sterilization result. Perform the sterilization only on a properly cleaned and dried product.

Proceed as follows:

1. Load all parts into the reprocessing tray and small parts like the light guide adapter, instrument receptacle, three-slot membrane or sealing cap into the integrated small component tray.

CAUTION! Non-sterile, single-use three-slot membrane and sealing cap. Infection risk for the patient. Sterilize the three-slot membrane and sealing cap prior to use and do not reuse.

NOTICE! The reprocessing tray is designed to accommodate one ureterorenoscope including instrument bridge.

Make sure that the product is fully disconnected from the tube connectors on the reprocessing tray.

2. Pack the tray with two single layers of polypropylene wrap.

CAUTION! Risk of scalding when loading the sterilizer. Wear suitable gloves.

3. Load the wrapped tray in the sterilizer.
4. Start the sterilization cycle in accordance with the manufacturer's instructions and instructions for use for the sterilizer.

For steam sterilization in a pre-vacuum autoclave cycle use the following validated parameters:

Cycle type	Pre-vacuum process
Pulses	4
Temperature	134 °C (273 °F)
Holding time	at least 5 min (effective sterilization time)
Drying time	at least 30 min

Note: It is the responsibility of the operator to establish that the sterilizer used meets the parameters above.

CAUTION! Risk of scalding when unloading the sterilizer. Wear suitable gloves.

NOTICE! Products can be shock sensitive when hot. Avoid knocking or shaking.

NOTICE! Damage resulting from a sudden change in temperature. Allow products to cool to room temperature; do not use additional cooling measures.

5. Remove the sterilized products from the sterilizer.
Ensure that the products remain sterile after reprocessing.
- This concludes steam sterilization.

5 Maintenance and Repairs

5.1 Troubleshooting Table

Issue	Possible causes	Fixes
Image cloudy	Glass surfaces soiled	Manual cleaning (section 4.6), then re-process, check the water quality
	Leaky, defective lens system	Send the endoscope in for repair
Image too dark, too little illumination	Glass surfaces soiled	Manual cleaning (section 4.6), then re-process, check the water quality
	Unsuitable light guide	Use a suitable light guide
	Light guide incorrectly attached to endoscope	Check that the light guide is seated properly
	Defective fiber optics	Check fiber optics (section 3.3)
	Defective light guide or light source	Check the light guide and light source
Image is yellowed	Soiled fiber optics	Manual cleaning (section 4.6), then re-process, check the water quality
	Light guide soiled or defective	Check the light guide (e.g. by illuminating a white surface)
Corrosion, staining, discoloration	Inadequate cleaning (e.g. remaining protein residue)	Manual cleaning (section 4.6), rub thoroughly if necessary, then reprocess
	Insufficient rinsing between the reprocessing stages, especially prior to sterilization	Rinse sufficiently between the reprocessing stages
	Excessive chloride concentration	Check water quality
	Excessive concentration of minerals (e.g. limescale) or organic substances	Check the water quality, only use demineralized water if necessary
	Heavy metal ions and/or silicates, elevated concentration of iron, copper, and manganese in the water or sterilizer steam	Check the water quality, only use demineralized water if necessary
	Contaminated, too frequently reused cleaning and disinfecting solution	The cleaning and disinfecting solution should be replaced regularly

Issue	Possible causes	Fixes
	Extraneous rust, e.g. resulting from rust contained in steam or damaged or non-corrosion-resistant instruments reprocessed at the same time	Check the supply systems, if multiple items are being reprocessed together, pay attention to material compatibility and signs of prior damage; prevent the different items from touching one another
	Contact corrosion	Avoid contact with other products

Table 5-1: Troubleshooting.

5.2 Repair

Should you need to arrange a repair for the device, please contact our Service Department.

When you return the device, please also enclose a detailed **description of the fault**.

WARNING

Contaminated device

Risk of infection

> Reprocess the product properly prior to shipping

Only ship devices that have been properly reprocessed.

The device should be returned in its original packaging, if possible.

Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

6 Product Data

6.1 Technical Data

6.1.1 Ureterorenoscopes

Item no.	STS-6075-425HA
Direction of view	8°
Field of view	85°
Working length	425 mm
Outer diameter of insertion portion (distal end / maximum diameter)	7.5 / 9.8 Fr

Item no.	STS-6075-310HA
Direction of view	8°
Field of view	85°
Working length	310 mm
Outer diameter of insertion portion (distal end / maximum diameter)	7.5 / 9.8 Fr

Item no.	STS-6060-425HA
Direction of view	8°
Field of view	100°
Working length	425 mm
Outer diameter of insertion portion (distal end / maximum diameter)	6.0 / 7.5 Fr

Item no.	STS-6060-310HA
Direction of view	8°
Field of view	100°
Working length	310 mm
Outer diameter of insertion portion (distal end / maximum diameter)	6.0 / 7.5 Fr

6.2 Ambient Conditions

Transport and storage conditions

Temperature	-20 °C to +70 °C
Relative air humidity	5 % to 95 %
Atmospheric pressure	70 kPa to 106 kPa

Reprocessed devices must be protected against recontamination, in a dark, dry, well-ventilated, dust-free and light-controlled location with a consistent temperature.

Direct sunlight, high temperatures, high humidity or radiation can damage the device or present a risk of infection.

When placing into storage, make sure that the device cannot be damaged by other instruments. It is therefore best to store the device individually or use containers in which it can be secured.

Operating conditions

Temperature	+15 °C to +32 °C
Relative air humidity	30 % to 85 %
Atmospheric pressure	70 kPa to 106 kPa

6.3 Spare Parts and Accessories

Use original spare parts and accessories only.

Photo	Designation	Item Number
	Sealing caps, unsterile single-use item, pack of 10	STS-ACC-5000-CAPD-000
	Stopcock plugs, steel, pack of 12 (4x clockwise, 8x counterclockwise)	STS-ACC-6000-SCPA-000
	Lubricant PARALIQ GB 363, 40 g	RFP00200
	Three-slot membrane, unsterile single-use item, pack of 10	STS-ACC-6000-MEMD-000
	Light guide adapter for Storz systems	STS-ACC-0000-LUGA
	Light guide adapter for Wolf systems	STS-ACC-0000-LUGA-RWF
	Adapter for EMS Swiss Lithoclast 1	STS-ACC-6000-EMSA-LC1
	Adapter for EMS Swiss Lithoclast 2	STS-ACC-6000-EMSA-LC2
	Instrument bridge	STS-URO-6000-BRGA-002
	Reprocessing tray for URS	STS-ACC-5600-TRYA-000
	Cleaning brush for working channel (URS with an outer diameter of insertion portion on the distal end of 7.5 Fr.)	STS-ACC-6075-520N-BRU
	Cleaning brush for working channel (URS with an outer diameter of insertion portion on the distal end of 6 Fr.)	STS-ACC-6060-520N-BRU
	Cleaning brush for instrument bridge irrigation channels	STS-ACC-5030-100N-BRU
	Cleaning brush for instrument bridge working channel	STS-ACC-5043-100N-BRU

Table 6-1: Spare parts and accessories.

7 Disposal

WARNING

Contaminated device

Risk of infection

- > The device must be reprocessed prior to disposal

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

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