STS-IFU-023-REV1

**strauss**<sup>®</sup>  
S U R G I C A L

## Instructions For Use

Sinuscopes  
USA

TPA764-025-01\_B

Version: B

2019 - Dec - 16

|          | <b>Content</b>  | <b>Page</b> |
|----------|---|-------------|
| <b>0</b> | <b>Important Information About this Document .....</b>          | <b>5</b>    |
| 0.1      | Scope of Validity, Identification, Purpose.....                 | 5           |
| 0.2      | Target Group.....   | 5           |
| 0.3      | Using and Storing this Document.....                            | 5           |
| <b>1</b> | <b>General Information About the Device.....</b>                | <b>6</b>    |
| 1.1      | Device Description.....   | 6           |
| 1.1.1    | Performance Characteristics and Function .....                  | 6           |
| 1.1.2    | Visual Overview .....   | 6           |
| 1.1.3    | Required Accessories .....                                      | 6           |
| 1.2      | Usage.....  | 7           |
| 1.2.1    | Intended Use.....   | 7           |
| 1.2.2    | Indications for Use.....  | 7           |
| 1.2.3    | Contraindications .....   | 7           |
| 1.3      | Marking .....   | 8           |
| 1.3.1    | Pictograms and Information on the Device and<br>Packaging ..... | 8           |
| 1.3.2    | Pictograms in this Document.....                                | 9           |
| 1.3.3    | Direction of View Color Coding.....                             | 9           |
| 1.4      | Service Department Contact Details .....                        | 9           |
| <b>2</b> | <b>General Safety Information .....</b>                         | <b>10</b>   |
| 2.1      | Safety Messages in this Document.....                           | 10          |
| 2.1.1    | Safety Messages at the Start of a Chapter .....                 | 10          |
| 2.1.2    | Safety Messages in the Body of the Text.....                    | 10          |
| 2.2      | Product Safety .....  | 11          |
| 2.2.1    | Basic Safety Information .....                                  | 11          |
| 2.2.2    | Staff Qualifications.....                                       | 11          |
| 2.2.3    | Compatibility with Light Guides and Light Sources .....         | 12          |
| <b>3</b> | <b>Operation .....</b>  | <b>13</b>   |
| 3.1      | Safety Instructions.....  | 13          |
| 3.2      | Staff Qualifications.....                                       | 14          |
| 3.3      | Technical Inspection Prior to Use.....                          | 15          |
| 3.3.1    | Visual inspection.....  | 15          |
| 3.3.2    | Function Check .....  | 16          |
| 3.4      | Connecting a Light Guide.....                                   | 16          |
| 3.5      | Post-Use and Pre-Cleaning Activities .....                      | 16          |
| <b>4</b> | <b>Reprocessing .....</b>                                       | <b>17</b>   |
| 4.1      | Safety Instructions.....  | 17          |
| 4.2      | Staff Qualifications.....                                       | 18          |
| 4.3      | Validated Procedures.....                                       | 19          |
| 4.4      | The Process Flow .....  | 19          |
| 4.5      | Cleaning Agents .....   | 20          |
| 4.6      | Manual Cleaning with Immersion .....                            | 20          |
| 4.6.1    | Overview .....  | 20          |
| 4.6.2    | Performance .....   | 21          |
| 4.7      | Automated Cleaning.....   | 22          |
| 4.7.1    | Overview .....  | 22          |
| 4.7.2    | Performance .....   | 22          |
| 4.8      | Sterilization.....  | 23          |
| 4.8.1    | Steam Sterilization.....  | 23          |
| 4.8.2    | STERRAD® Sterilization Procedure (100S, NX,<br>100NX).....      | 24          |
| <b>5</b> | <b>Maintenance and Repairs.....</b>                             | <b>25</b>   |
| 5.1      | Troubleshooting Table .....                                     | 25          |
| 5.2      | Repair.....   | 26          |
| <b>6</b> | <b>Product Data .....</b>                                       | <b>27</b>   |
| 6.1      | Technical Data.....   | 27          |
| 6.2      | Ambient Conditions .....  | 29          |
| <b>7</b> | <b>Disposal.....</b>  | <b>30</b>   |

## 0 Important Information About this Document

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These instructions for use apply to the following rigid endoscopes:

- STS-7027-110ASA
- STS-7027-110BSA
- STS-7027-110CSA
- STS-7040-175AA
- STS-7040-175BA
- STS-7040-175CA
- STS-7040-175BSA
- STS-7040-175FA
- STS-7027-187ASA
- STS-7027-187BSA
- STS-7027-187CSA
- STS-7020-110AA
- STS-7020-110BA

These instructions for use are an integral component of the device and contain all the information required by users and operators for safe and proper use.

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## 1 General Information About the Device

### 1.1 Device Description

#### 1.1.1 Performance Characteristics and Function

Rigid endoscopes permit the visualization of the inside of the body and body cavities. The approach to the body cavity may be through a natural orifice or a surgically created opening.

**Structure and function** A rigid endoscope consists of a fiber optic cable and sensitive image transmission system with eyepiece.

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The distal end of the endoscope features an object lens that captures the image from inside the body. The image is sent through the image transmission system to the eyepiece.

The eyepiece cup connects the eyepiece to an optocoupler that transmits the image to a camera. The camera's controller converts the signal so that the image can be displayed on a monitor.

#### 1.1.2 Visual Overview


Figure 1-1: Overview.

- 1 Lens
- 2 Eyepiece cup
- 3 Light guide connector
- 4 Adapter for connecting Wolf light guides
- 5 Adapter for connecting Storz light guides

#### 1.1.3 Required Accessories

**Additional instruments** Additional instruments are required for surgical procedures within the ENT area. It is the surgeon's responsibility to make sure the instruments required are assembled and available for use.

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Sinuscopes are used to visualize the frontal, maxillary and nasal sinuses.

The devices are intended to be used in combination with medical electrical devices that are, as a **minimum**, compliant with the requirements for applied parts classified as type **BF** (that is parts offering increased protection against electric shock), according to IEC 60601-1.

This insulation barrier can be realized by the devices themselves or by the connecting cables to the endoscope and must be assured for every connection between the endoscope and connected devices.

Devices or connecting cables that comply with the requirements for applied parts classified as type BF feature the adjacent symbol.

Type CF is the most stringent classification for devices and connecting cables considered applied parts within the meaning of IEC 60601-1. Devices bearing the CF symbol (the adjacent heart symbol) are approved for direct conductive contact with the heart.

Defibrillation-proof devices or connecting cables are also compatible. These will be marked with one of the adjacent symbols depending on the degree of protection against electric shock offered.

### 1.2.2 Indications for Use

The sinuscopes are intended for use in otolaryngology and head and neck procedures, including rhinology, and endoscopic plastic and reconstructive surgery.

### 1.2.3 Contraindications

Depending on the condition of the patient, contraindications may exist in respect of the patient's general condition or the specific symptoms.

The decision to perform an endoscopic procedure rests with the responsible surgeon and should be made on the basis of an individual benefit/risk analysis.

Contraindications for surgical implementation include:

- General condition of the patient
- Pregnancy
- Anesthesia intolerance
- The taking of medication

## 1.3 Marking

### 1.3.1 Pictograms and Information on the Device and Packaging

This section describes the different pictograms that feature on the device or the packaging.

|   |   |
|---|---|
|    | The device is suitable for autoclave sterilization  |
|    | Item number   |
|    | Serial number   |
|    | Quantity of products of a particular item number  |
|    | Direction of view   |
|    | Manufacturer  |
|  | Date of manufacture   |
|  | Caution (IEC 60601-1 3rd edition) /<br>Take note of accompanying documents<br>(IEC 60601-1 2nd edition) |
|  | Non-sterile, sterilize prior to each use  |
|  | Not for use in vicinity of MRI equipment  |
|  | Permissible storage temperature   |
|  | Permissible relative air humidity during storage  |
|  | Permissible atmospheric pressure during storage   |
|  | Caution: US federal law restricts this device to sale by or<br>on order of a physician.                 |

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## 1.4 Service Department Contact Details

If you have any questions about our products, installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

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The warnings described in this section will be listed, where applicable, within the body of the instructions directly prior to any steps that carry a particular risk.

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Risks – despite the utmost quality

Although this device meets the current state of the art, risks could still arise during initial operation, use, reprocessing, or maintenance.

Therefore, it is important that you read through these instructions for use carefully and observe the warnings indicated.

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**WARNING!** This device is supplied in a non-sterile condition, which carries a risk of infection. Carry out sterile processing of the device and its accessories prior to first use.

**WARNING!** Unauthorized modifications to the device carry a risk of serious injury. Do not make any unauthorized modifications.

**WARNING!** Component failure during a surgical procedure carries a risk to the patient. Keep spares on hand for emergency replacement.

**WARNING!** High-intensity light source. Risk of injury to eyes. Do not look directly into the distal end of the light guide and into the light exit of the endoscope.

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**WARNING!** Using the device during operation of a defibrillator carries a risk to persons. Before discharging a defibrillator, remove device from the operating field.

**CAUTION!** Rough handling can damage the device and therefore carries a risk to the patient. Handle the device with appropriate care. If the device is dropped or subjected to high mechanical stress, stop using it and send to the manufacturer for inspection.

**WARNING!** Do not point the light beam into the patient and/or other medically credentialed healthcare provider's eyes. Looking into the light beam may cause eye damage.

**WARNING!** When in operation, the external light source heats up the endoscope. Switch off the light source before putting the endoscope down and do not place it in the proximity of the patient (as per IEC 60601-1-1).

### 2.2.2 Staff Qualifications

Observe the qualification requirements

Specific qualifications are required for installation, operation and reprocessing of the device. The qualifications required are stated in the respective chapters and must be observed.

### 2.2.3 Compatibility with Light Guides and Light Sources

The use of different adapters allows the endoscopes to be used with light sources and light guides manufactured by all standard manufacturers. Use legally marketed light sources and light guides only!

To avoid excessive heat development on the coupling points, light guides and endoscopes that match each other must be selected appropriately.

The active diameter of the light guide must be identical to or smaller than the active diameter of the endoscope light post.

For endoscopes with a light post with an active diameter of less than 4 mm, use a light guide with active diameter of 2.6 mm - 3.5 mm. For endoscopes with a light post with an active diameter of more than 4 mm, there is no restriction. See [section 6.1](#) for technical data.

**Note:** High temperatures may be present at the coupling points of the light guide if used in conjunction with external endoscopic xenon light sources. This may lead to a transmission decrease and therefore increasing heating at the coupling points on the light guides. As a result, when combined with xenon light sources, only high-temperature light guides should be used.

The light output of the used light source is not restricted. However, use the minimum level of illumination necessary to satisfactorily view the target area.

**Note:** Only use the endoscopes in combination with light sources that offer a luminous flux between 300lm and 1500lm.

Only use the endoscopes in combination with light sources that offer a wavelength between 400nm and 700nm.

## 3 Operation

### 3.1 Safety Instructions

#### **WARNING**

##### **Interactions between devices in simultaneous use (e.g. lasers, HF surgery)**

Risk to the patient and user, image distortion, damage to the device

- > Ensure that all the devices in use are, as a minimum, compliant with the requirements for applied parts classified as type BF according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used

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##### **Use of non-sterile parts**

Infection risk for the patient

- > Use only properly reprocessed endoscopes and endoscopic accessories
- > Always carry out a visual inspection prior to use
- > Anchor the light guide in the operating field and secure against slipping, allowing some slack for maneuvering

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**⚠ CAUTION**

**High temperatures in combination with light source usage**

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- > Use only light guides that are suitable for use with the endoscope
- > Avoid the prolonged use of intense light
- > You should always use the lowest light output setting that will allow you to complete the task at hand
- > Do not touch the light source in the immediate vicinity of the lamp
- > Do not touch light guide connectors
- > Do not allow the distal end or the outer tube of the endoscope to come into contact with patient tissue or with combustible or heat-sensitive materials.
- > Do not exceed the maximum permitted ambient temperature for the light source

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Risk of damage to the device

- > Set down carefully without bending or knocking
- > Avoid scratching the surface
- > Take care when removing instruments from the operating field

### 3.2 Staff Qualifications

The device may be used only by physicians, surgeons, and surgical assistants who have received instruction in how to use the device and who have the requisite training or advanced training, knowledge, and practical experience in the endoscopic procedure to be applied as defined by the provisions in place at the site of operation.

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*Proceed as follows:*

- |   |  |
|---|--|
| Is there any external damage?             | <p>1. Check the endoscope and all the components to be used for external damage.</p> <p>Do not use the device if it has any sharp edges or corners, bulges, or rough surfaces that might cause injury to the patient.</p>  |
| Are the fiber optics in working order?    | <p>2. Check that the fiber optics in the endoscope are in full working order.</p> <p>To do this, hold the distal end toward a bright light (not a cold light source) and the light guide connector toward your eyes. Gently move the endoscope from side to side, taking note of the brightness of the fibers.</p> <p>If more than 20 % of the fibers remain dark, the endoscope must not be used and a repair must be arranged.</p> |
| Are the surfaces clean and smooth?        | <p>3. Check that the glass surfaces and fiber optic end faces of the endoscope are clean and smooth.</p> <p>Do not use the endoscope if the surfaces are soiled or scratched. Set aside the endoscope and follow the steps described in <u>section 5.1</u>.</p>  |
| Is the image clear, bright, and complete? | <p>4. Make sure you can see a clear, bright, complete image.</p> <p>Look through the eyepiece and assess the image quality.</p> <p>Do not use the endoscope if the image is yellowed, dark, speckled, or cropped. Set aside the endoscope and follow the steps described in <u>section 5.1</u>.</p>  |
| Is the OR equipment mutually compatible?  | <p>5. Ensure that all the equipment to be used for the procedure is mutually compatible.</p> <p>Only use the OR equipment if the isolation barrier for protecting against electric shock, <b>as a minimum</b>, complies with the requirements for applied parts classified as type <b>BF</b> according to IEC 60601-1 (see <u>Intended Use</u>).</p>   |
- This concludes the inspection.

### 3.3.2 Function Check

#### Performing a function check

*A function check must be performed **before every medical procedure.***

*Proceed as follows:*

1. Connect a camera.
  - Switch on the system 2. Switch on all the system components that you intend using for the procedure.
  - Focus the image 3. Direct the camera head at an object in the room and focus the image.  
Check that you are able to achieve a sharp, bright, good quality image.  
Do not use the system if the image contains vertical lines or color variations, flickers, or you are unable to achieve a sharp, bright, good quality image.
- This concludes the function check.

### 3.4 Connecting a Light Guide

#### Connecting the light guide

*Different adapters are needed depending on the type of light guide being connected to the endoscope.*

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*Perform pre-cleaning **immediately after use.***

*Proceed as follows:*

1. Remove and discard the sterile drape, if present.
2. Decouple the endoscope from the optocoupler.
3. Detach the light guide and all detachable parts.

**NOTICE!** The eyepiece cup is non-detachable.

4. Pre-clean the endoscope at the operating table at the conclusion of the surgery by wiping with a lint-free wipe wetted (i.e., wetted but not dripping; squeeze out excess solution) with an enzymatic cleaning solution for a minimum of **one (1) minute** and until visibly clean.

Recommended cleaning agents see [section 4.5](#)

Before wetting the wipe with the enzymatic cleaning solution, it should be prepared per the manufacturer's instructions.

5. Place all the parts into a dry removal container and seal.
6. Arrange for reprocessing, ensuring that all of the device components are reprocessed within 6 hours.

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#### **CAUTION**

##### **Improper cleaning**

Risk to patient as a result of inadequate cleaning, damage to the device

- > Use an automatic washer that meets the performance requirements of ISO 15883-1
- > Load the automatic washer with care: all contents must be rinsed and cleaned fully (no rinsing 'blind spots')
- > Ensure that the automatic washer is properly maintained
- > Only use cleaning agents that are approved for the device
- > The devices must be pre-cleaned immediately after use and reprocessed within 6 hours
- > When pre-cleaning, do not use fixing temperatures in excess of 45 °C
- > When pre-cleaning, do not use fixing detergents (containing aldehyde or alcohol)

#### **NOTICE**

##### **Contact with chloride-based solutions**

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- > Avoid contact with chloride-based solutions, which can be found, for example, in surgical residues, lotions, drugs, saline solutions, and cleaning agents
- > After contact with chloride-based solutions, rinse devices thoroughly with deionized water and dry completely

#### **NOTICE**

##### **Use of ultrasonic cleaning bath**

Damage to the endoscope

- > Do not expose endoscope to ultrasonic cleaning bath

**NOTICE****Unsuitable detergents and process chemicals**

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- > Only use cleaning agents that are approved for the device
- > Only use process chemicals that are compatible with the device's materials according to the chemical manufacturer's recommendations
- > Follow all the chemical manufacturer's application specifications regarding temperature, concentration, and contact time
- > Do not use process chemicals that can cause stress cracking or brittleness of plastic materials

**4.2 Staff Qualifications**

In many countries, the qualifications required by personnel responsible for the reprocessing and sterilization of medical devices are regulated by law.

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Professional expertise may be obtained through specialized further education or through training and practical experience, and enhanced through appropriate further training if necessary.

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- Manual cleaning
- Automated cleaning
- Sterilization

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### 4.4 The Process Flow

**Pre-cleaning** The device has to be prepared for reprocessing immediately after use by pre-cleaning.

Refer to [section 3.5](#) for information on how to perform pre-cleaning.

**Steps** The reprocessing procedure described in this document consists of the following steps:

- Pre-cleaning immediately after use
- Manual or automated cleaning
- Sterilization

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**Cleaning prior to sterilization** The endoscope must be cleaned thoroughly prior to sterilization. It is essential that the sterilization medium accesses all parts of the endoscope.

**Automated cleaning and steam sterilization** The optimum and most reliable reprocessing results are achieved using automated cleaning with subsequent steam sterilization using the fractionated vacuum steam procedure.

**Laws and standards** Adhere to national statutory regulations, national and international standards and directives, and the infection control practices that are in place at your institution for reprocessing.

**Reprocessing prior to each use** New factory products and recently repaired devices must both undergo the entire reprocessing procedure before use.

The alternating use of different reprocessing procedures can cause premature aging of the product.

Excessive doses of neutralizers and cleaning agents can damage the device and cause laser engravings to fade.

**www.a-k-i.org** Please see [www.a-k-i.org](http://www.a-k-i.org) for more detailed information about hygienically safe, value-preserving reprocessing methods that also preserve the materials.

## 4.5 Cleaning Agents

The following products can be used for **manual cleaning**:

- Prolystica 2X Concentrate Enzymatic Presoak and Cleaner (Steris)
- Cidezyme/Enzol (Johnson & Johnson)

The following detergent can be used for **automated cleaning**:

- neodisher MediClean forte 0.5 %  
(Chem. Fabrik Dr. Weigert GmbH & Co. KG)

Wherever possible, use the above cleaning agents only.

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Before use, read the manufacturer's user information carefully and follow the specifications regarding concentration, temperature, usage time, water quantities, and contact time.

## 4.6 Manual Cleaning with Immersion

### 4.6.1 Overview

| Stage | Work Step     | Temperature (°C/°F)                 | Time (mins.) | Water Quality     | Cleaning Solution           |
|-------|---------------|-------------------------------------|--------------|-------------------|-----------------------------|
| I     | Cleaning      | per the manufacturer's instructions | 2-5          | Tap water         | Enzymatic cleaning solution |
| II    | Rinsing 2x    | < 45/113                            | 2x ≥1        | Tap water         | ---                         |
| III   | Final rinsing | < 45/113                            | ≥1           | High purity water | ---                         |
| IV    | Drying        | ---                                 | ---          | ---               | ---                         |

Table 4-1: Overview of manual cleaning with immersion.

## 4.6.2 Performance

**Complete stage I: Cleaning**

*Prepare the cleaning solution per the manufacturer's instructions.*

*Proceed as follows:*

**NOTICE!** Damage resulting from rough handling. Handle the device with appropriate care. Do not knock or bend. Set down carefully.

**NOTICE!** Scratch-sensitive surfaces. Risk of corrosion. Do not use metal brushes, metal objects, or abrasive cleaners. To remove soiling on optical surfaces, use pads impregnated with cleaning solution only.

- Soak for 2-5 mins
1. Fully immerse the endoscope in an enzymatic cleaning solution bath prepared per the manufacturer's instructions for **two (2) to five (5) minutes**.

All accessible surfaces must remain immersed in the cleaning solution bath throughout the entire cleaning time.

- Clean for at least 1 minute
2. After soaking in the cleaning solution, clean all external surfaces while the endoscope is fully immersed in the cleaning solution using a soft, clean, lint-free wipe or a soft bristled brush for a minimum of **one (1) minute**.

**Complete stage II: Rinsing**

- 2 rinses in tap water of at least 1 minute
1. Fully immerse the endoscope in a tap water bath (less than 45 °C / 113 °F) and rinse all accessible surfaces twice for a minimum of **one (1) minute** each.

Use fresh water for each rinse.

2. Allow water to drip off for a sufficient length of time.

**Complete stage III: Final rinsing**

- For at least 1 min. in high purity water
1. Fully immerse the endoscope in a high purity water bath (less than 45 °C / 113 °F) and rinse all accessible surfaces for a minimum of **one (1) minute** thoroughly.
  2. Allow water to drip off for a sufficient length of time.

**Complete stage IV: Drying**

1. Thoroughly dry the endoscope with a clean, lint-free wipe or clean, lint-free cotton weave operating room towels. If available, dry all accessible surfaces with medical-quality filtered compressed air (pmax = 0.5 bar).
  2. Visually inspect the endoscope in a well-lit area; it should be completely clean and dry.  
Use a magnifying glass. Repeat the cleaning process if necessary.
- This concludes the manual cleaning process.

## 4.7 Automated Cleaning

### 4.7.1 Overview

| Stage | Step                 | Temperature (°C/°F)                             | Time (mins.) | Cleaning Solution / Water Quality |
|-------|----------------------|---|--------------|-----------------------------------|
| I     | Pre-rinsing          | Cold  | 2            | Tap water                         |
| II    | Cleaning             | 60/140  | 10           | Alkaline detergent                |
| III   | Intermediate rinsing | per the equipment manufacturer's standard cycle | 1            | Tap water                         |
| IV    | Intermediate rinsing | per the equipment manufacturer's standard cycle | 1            | Tap water                         |
| V     | Final rinsing        | 90/194  | 5            | Purified water                    |
| VI    | Dry                  | 98.8/210  | 30           | ---                               |

Table 4-2: Automated cleaning overview.

**Note:** Automated reprocessing can cause color-anodized or plastic components (e.g. standard rings and the eyepiece cup) to fade.

### 4.7.2 Performance

#### Performing automated cleaning

Use an automatic washer that has been validated as effective and meets the performance requirements of ISO 15883-1 (or the respective country-specific version thereof).

Dry removal of the products is recommended for automated cleaning.

In the case of wet removal, do not use foaming detergents and rinse the products thoroughly prior to the automated cleaning.

Proceed as follows:

- Place the endoscopes in the automatic washer.  
Avoid areas that cannot be reached during rinsing.
- Start the cleaning cycle described in [section 4.7.1](#) in accordance with the manufacturer's instructions and instruction for use for the automatic washer.

**CAUTION!** Risk of scalding when unloading the automatic washer. Be sure to wear gloves.

- Remove the endoscopes from the automatic washer.

Visual inspection

- Visually inspect the endoscopes in a well-lit area; they should be completely clean, dry, and undamaged.

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Immediately set aside any damaged components.

- This concludes the automated cleaning process.

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*Proceed as follows:*

1. Double wrap the endoscope with two single layers of FDA-cleared sterilization wrap using a sequential double wrapping technique.

**CAUTION!**

Risk of scalding when loading the sterilizer. Wear suitable gloves.

2. Load the wrapped endoscopes in the sterilizer.

Devices should be positioned to permit adequate steam circulation and penetration, good air circulation, and sufficient condensate drainage. Do not stack wrapped devices in the sterilization chamber.

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3. Start the sterilization cycle in accordance with the manufacturer's instructions and instruction for use for the sterilizer.

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Holding time: min. 4 minutes  
(effective sterilization time)

Minimum dry time: min. 30 minutes

Note: It is the responsibility of the user to establish that the sterilizer used meets the steam sterilization parameters above.

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- This concludes sterilization.

#### 4.8.2 STERRAD® Sterilization Procedure (100S, NX, 100NX)

The STERRAD® sterilization systems, manufactured by Advanced Sterilization Products (ASP), use low-temperature, hydrogen peroxide gas plasma technology for terminal sterilization of properly cleaned, rinsed, and dried reusable medical devices.

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*Proceed as follows:*

1. Double wrap the endoscope with two single layers of FDA-cleared sterilization wrap using a sequential double wrapping technique.
2. Load the wrapped endoscopes in the sterilizer.  
Devices should be positioned to permit adequate plasma circulation and penetration. Do not allow any item to touch the wall of the sterilizer. Do not stack wrapped devices in the sterilization chamber.
3. Start the sterilization cycle in accordance with the manufacturer's instructions and instructions for use for the steam sterilizer.

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- STERRAD® 100S short cycle
  - STERRAD® NX standard cycle
  - STERRAD® 100NX standard cycle
4. Remove the sterilized products from the sterilizer.  
Ensure that the products remain sterile after reprocessing.
- This concludes sterilization.

## 5 Maintenance and Repairs

### 5.1 Troubleshooting Table

| Issue                                   | Possible Causes   | Fixes  |
|---|---|--|
| Image cloudy                            | Glass surfaces soiled   | Clean manually as per <a href="#">section 4.6</a> , then reprocess, check the water quality      |
|   | Stubborn deposits on glass surfaces   | Remove deposits as per the Cleaning chapter, check the water quality                             |
|   | Leaky, defective lens system  | Send the endoscope in for repair   |
| Image too dark, too little illumination | Glass surfaces soiled   | Clean manually as per <a href="#">section 4.6</a> , then reprocess, check the water quality      |
|   | Stubborn deposits on glass surfaces   | Remove deposits as per the Cleaning chapter, check the water quality                             |
|   | Unsuitable light guide  | Use a suitable light guide   |
|   | Light guide incorrectly attached to endoscope   | Check that the light guide is seated properly  |
|   | Defective fiber optics  | Check fiber optics as per <a href="#">section 3.3</a>  |
|   | Defective light guide or light source   | Check the light guide and light source   |
| Image is yellowed                       | Soiled fiber optics   | Clean manually as per <a href="#">section 4.6</a> , then reprocess, check the water quality      |
|   | Light guide soiled or defective   | Check the light guide (e.g. by illuminating a white surface)                                     |
| Corrosion, staining, discoloration      | Inadequate cleaning (e.g. remaining protein residue)  | Clean manually as per <a href="#">section 4.6</a> , if necessary, rub thoroughly, then reprocess |
|   | Insufficient rinsing between the reprocessing stages, especially prior to sterilization                                   | Rinse sufficiently between the reprocessing stages   |
|   | Excessive chloride concentration  | Check water quality  |
|   | Excessive concentration of minerals (e.g. limescale) or organic substances  | Check the water quality, only use demineralized water if necessary                               |
|   | Heavy metal ions and/or silicates, elevated concentration of iron, copper, and manganese in the water or sterilizer steam | Check the water quality, only use demineralized water if necessary                               |

| Issue | Possible Causes   | Fixes  |
|-------|---|--|
|       | Contaminated, frequently reused cleaning solution   | The cleaning solution should be replaced regularly |
|       | Extraneous rust, e.g. resulting from rust contained in steam or damaged or non-corrosion-resistant instruments reprocessed at the same time | Avoid contact with other components                |
|       | Contact corrosion   | Avoid contact with other components                |

Table 5-1: Troubleshooting.

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#### Contaminated device

Risk of infection

> Reprocess the device properly prior to shipping.

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### 6.1 Technical Data

#### **Item no. STS-7027-110ASA**

|                                  |          |
|----------------------------------|----------|
| Direction of view                | 0°       |
| Field of view                    | standard |
| Working length                   | 110 mm   |
| Outer diameter of insertion tube | 2.7 mm   |
| Active diameter of light post    | 4.5 mm   |

#### **Item no. STS-7027-110BSA**

|                                  |          |
|----------------------------------|----------|
| Direction of view                | 30°      |
| Field of view                    | standard |
| Working length                   | 110 mm   |
| Outer diameter of insertion tube | 2.7 mm   |
| Active diameter of light post    | 4.5 mm   |

#### **Item no. STS-7027-110CSA**

|                                  |          |
|----------------------------------|----------|
| Direction of view                | 70°      |
| Field of view                    | standard |
| Working length                   | 110 mm   |
| Outer diameter of insertion tube | 2.7 mm   |
| Active diameter of light post    | 4.5 mm   |

#### **Item no. STS-7040-175AA**

|                                  |            |
|----------------------------------|------------|
| Direction of view                | 0°         |
| Field of view                    | wide angle |
| Working length                   | 175 mm     |
| Outer diameter of insertion tube | 4 mm       |
| Active diameter of light post    | 4.5 mm     |

#### **Item no. STS-7040-175BA**

|                                  |            |
|----------------------------------|------------|
| Direction of view                | 30°        |
| Field of view                    | wide angle |
| Working length                   | 175 mm     |
| Outer diameter of insertion tube | 4 mm       |
| Active diameter of light post    | 4.5 mm     |

#### **Item no. STS-7040-175CA**

|                                  |            |
|----------------------------------|------------|
| Direction of view                | 70°        |
| Field of view                    | wide angle |
| Working length                   | 175 mm     |
| Outer diameter of insertion tube | 4 mm       |
| Active diameter of light post    | 4.5 mm     |

#### **Item no. STS-7040-175BSA**

|                                  |          |
|----------------------------------|----------|
| Direction of view                | 30°      |
| Field of view                    | standard |
| Working length                   | 175 mm   |
| Outer diameter of insertion tube | 4 mm     |
| Active diameter of light post    | 4.5 mm   |

#### **Item no. STS-7040-175FA**

|                                  |            |
|----------------------------------|------------|
| Direction of view                | 45°        |
| Field of view                    | wide angle |
| Working length                   | 175 mm     |
| Outer diameter of insertion tube | 4 mm       |
| Active diameter of light post    | 4.5 mm     |

|                                  |                        |
|----------------------------------|------------------------|
| <b>Item no.</b>                  | <b>STS-7027-187ASA</b> |
| Direction of view                | 0°                     |
| Field of view                    | standard               |
| Working length                   | 187 mm                 |
| Outer diameter of insertion tube | 2.7 mm                 |
| Active diameter of light post    | 4.5 mm                 |

|                                  |                         |
|----------------------------------|-------------------------|
| <b>Item no.</b>                  | <b>STS-7027-187-BSA</b> |
| Direction of view                | 30°                     |
| Field of view                    | standard                |
| Working length                   | 187 mm                  |
| Outer diameter of insertion tube | 2.7 mm                  |
| Active diameter of light post    | 4.5 mm                  |

|                                  |                         |
|----------------------------------|-------------------------|
| <b>Item no.</b>                  | <b>STS-7027-187-CSA</b> |
| Direction of view                | 70°                     |
| Field of view                    | standard                |
| Working length                   | 187 mm                  |
| Outer diameter of insertion tube | 2.7 mm                  |
| Active diameter of light post    | 4.5 mm                  |

|                                  |                       |
|----------------------------------|-----------------------|
| <b>Item no.</b>                  | <b>STS-7020-110AA</b> |
| Direction of view                | 0°                    |
| Field of view                    | wide angle            |
| Working length                   | 110 mm                |
| Outer diameter of insertion tube | 2 mm                  |
| Active diameter of light post    | 4.5 mm                |

|                                  |                       |
|----------------------------------|-----------------------|
| <b>Item no.</b>                  | <b>STS-7020-110BA</b> |
| Direction of view                | 30°                   |
| Field of view                    | wide angle            |
| Working length                   | 110 mm                |
| Outer diameter of insertion tube | 2 mm                  |
| Active diameter of light post    | 4.5 mm                |

## 6.2 Ambient Conditions

### Transport and storage conditions

|                      |                   |
|----------------------|-------------------|
| Temperature          | -20 °C to +70 °C  |
| Rel. air humidity    | 5 % to 95 %       |
| Atmospheric pressure | 70 kPa to 106 kPa |

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Web: [www.straussurgical.com](http://www.straussurgical.com)

## **Instructions For Use**

Sinuscopes  
Endoscopes for Sinuscopy  
(Not valid for USA)

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Sinuscopios

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Sinuscópios

TPA765-025-34\_A

Version: A

2019 - Dec - 18

|          | <b>Content</b>   | <b>Page</b> |
|----------|--|-------------|
| <b>0</b> | <b>Important Information about this Document</b>       | <b>10</b>   |
| 0.1      | Scope of Validity, Identification, Purpose             | 10          |
| 0.2      | Target Group   | 10          |
| 0.3      | Using and Storing this Document                        | 10          |
| <b>1</b> | <b>General Information about the Device</b>            | <b>11</b>   |
| 1.1      | Scope of Delivery                                      | 11          |
| 1.2      | Product Description                                    | 11          |
| 1.2.1    | Performance Characteristics and Function               | 11          |
| 1.2.2    | Visual Overview  | 11          |
| 1.2.3    | Recommended Light Guides                               | 12          |
| 1.3      | Usage  | 12          |
| 1.3.1    | Intended Use   | 12          |
| 1.3.2    | Indications  | 12          |
| 1.3.3    | Contraindications                                      | 12          |
| 1.4      | Conformity   | 12          |
| 1.4.1    | Standards and Directives                               | 12          |
| 1.4.2    | Medical Device Classification                          | 13          |
| 1.5      | Marking  | 13          |
| 1.5.1    | Pictograms and Information on the Device and Packaging | 13          |
| 1.5.2    | Pictograms in this Document                            | 14          |
| 1.5.3    | Direction of View Color Coding                         | 14          |
| 1.6      | Service Department Contact Details                     | 14          |
| <b>2</b> | <b>General Safety Information</b>                      | <b>15</b>   |
| 2.1      | Safety Messages in this Document                       | 15          |
| 2.1.1    | Safety Messages at the Start of a Chapter              | 15          |
| 2.1.2    | Safety Messages in the Body of the Text                | 15          |
| 2.2      | Product Safety   | 16          |
| 2.2.1    | Basic Safety Notices                                   | 16          |
| 2.2.2    | Staff Qualifications                                   | 16          |
| <b>3</b> | <b>Operation</b>                                       | <b>17</b>   |
| 3.1      | Safety Instructions                                    | 17          |
| 3.2      | Staff Qualifications                                   | 18          |
| 3.3      | Technical Inspection Prior to Use                      | 19          |
| 3.3.1    | Visual Inspection                                      | 19          |
| 3.3.2    | Function Check   | 20          |
| 3.4      | Connection of a Light Guide                            | 20          |
| 3.5      | Post-Use and Pre-Cleaning Activities                   | 20          |
| <b>4</b> | <b>Reprocessing</b>                                    | <b>21</b>   |
| 4.1      | Safety Notices   | 21          |
| 4.2      | Staff Qualifications                                   | 22          |
| 4.3      | Validated Procedures                                   | 22          |
| 4.4      | The Process Flow                                       | 22          |
| 4.5      | Cleaning and Disinfecting Agents                       | 23          |
| 4.6      | Manual Cleaning and Disinfection                       | 24          |
| 4.6.1    | Overview   | 24          |
| 4.6.2    | Performance  | 24          |
| 4.7      | Automated Cleaning and Thermal Disinfection            | 26          |
| 4.7.1    | Overview   | 26          |
| 4.7.2    | Performance  | 27          |
| 4.8      | Sterilization  | 28          |
| 4.8.1    | Steam Sterilization                                    | 28          |
| 4.8.2    | STERRAD® Sterilization Systems (100S, NX, 100NX)       | 29          |
| <b>5</b> | <b>Maintenance and Repairs</b>                         | <b>30</b>   |
| 5.1      | Troubleshooting Table                                  | 30          |
| 5.2      | Repair   | 31          |

|          |                                   |           |
|----------|-----------------------------------|-----------|
| <b>6</b> | <b>Product Data .....</b>         | <b>32</b> |
| 6.1      | Technical Data.....               | 32        |
| 6.2      | Ambient Conditions .....          | 34        |
| 6.3      | Spare Parts and Accessories ..... | 35        |
| <b>7</b> | <b>Disposal.....</b>              | <b>36</b> |

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**7 Descarte ..... 92**

## 0 Important Information about this Document

### 0.1 Scope of Validity, Identification, Purpose

These instructions for use apply to the rigid endoscopes listed below:

#### Item numbers

STS-7027-110ASA  
STS-7027-110BSA  
STS-7027-110CSA  
STS-7040-175AA  
STS-7040-175BA  
STS-7040-175CA  
STS-7040-175BSA  
STS-7040-175FA  
STS-7027-187ASA  
STS-7027-187BSA  
STS-7027-187CSA  
STS-7020-110AA  
STS-7020-110BA

These instructions for use are an integral component of the product and contain all the information required by users and operators for safe and proper use.

These instructions for use do not describe endoscopic procedures or techniques during surgery.

### 0.2 Target Group

These instructions for use are intended for physicians, medical assistants, and employees of central sterile supply who are entrusted with the operation, handling, and reprocessing of the device.

### 0.3 Using and Storing this Document

These instructions for use must be stored in a defined location so that they may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

## 1 General Information about the Device

### 1.1 Scope of Delivery

Check scope of delivery! Check the delivery against the delivery note for completeness and damage.

The delivery left our premises in perfect condition. If you have any cause for complaint, however, please contact our Service Department.

### 1.2 Product Description

#### 1.2.1 Performance Characteristics and Function

Rigid endoscopes permit the visualization of the inside of the body and body cavities.

Depending on the endoscopic discipline and anatomical region, the approach to the body cavity may be through a natural orifice or a surgically created opening.

Structure and function Rigid endoscopes consist of a fiber optic cable and sensitive image transmission system with eyepiece.

The fiber optic cable is used to illuminate the site inside the body. The connector for connecting the light guide to the light source is situated at the proximal end of the endoscope. The adapters required to connect the light guide are included in the scope of delivery.

The distal end of the endoscope features an objective lens that captures the image from inside the body. The image is sent through the image transmission system to the eyepiece.

The eyepiece can be used with the naked eye or connected to an endocoupler, which transmits the image to a camera. The camera's controller converts the signal so that the image can be displayed on a monitor.

Special instruments Special instruments are required for endoscopic procedures.

It is the surgeon's responsibility to make sure the instruments required for a particular endoscopic procedure are assembled and available for use.

#### 1.2.2 Visual Overview

A sinuscope design illustration is provided below.


Figure 1-1: Sinuscope illustration.

- 1** Objective
- 2** Eyepiece
- 3** Light guide connector (ACMI)
- 4** Adapter for connecting Wolf light guides
- 5** Adapter for connecting Storz light guides

### 1.2.3 Recommended Light Guides

The specifications for compatible light guides are:

- Fiber bundles of 3.5 to 4.8 mm Ø
- Length up to 300 cm

## 1.3 Usage

### 1.3.1 Intended Use

Rigid endoscopes without a working channel permit visualization of the inside of the body and body cavities.

Sinusscopes are intended for use in otorhinolaryngology, head and neck surgery. They are used to visualize and illuminate particular body cavities during diagnostic and therapeutic procedures. These include the frontal, maxillary and nasal sinuses.

The rigid endoscopes described herein are intended to be used in combination with medical electrical devices that meet at **minimum** the type **BF** classification requirements according to IEC 60601-1 for elevated protection from electric shock.


This insulation barrier can be realized by the medical electrical devices themselves or by the connection cables to the endoscope and must be assured for every connection between the endoscope and connected devices.

Devices or connection cables that meet the type BF classification requirements feature the adjacent symbol.

Devices and connection cables that meet the requirements for applied parts classified as type CF offer an even higher level of protection. They bear the adjacent heart symbol and are approved for direct conductive contact with the heart.

Defibrillation-proof devices or connection cables are also compatible. These will be marked with one of the adjacent symbols depending on the degree of protection against electric shock offered.

### 1.3.2 Indications

The indications for an endoscopic application depend on the medical condition of the patient and the individual risk/benefit analysis by the surgeon.

### 1.3.3 Contraindications

Contraindications may exist depending on the patient's general condition or specific symptoms.

The decision to perform an endoscopic procedure remains with the responsible surgeon and must be made on the basis of an individual risk/benefit analysis.

## 1.4 Conformity

### 1.4.1 Standards and Directives

Rigid endoscopes are compliant with the requirements of the following applicable standard:

- **IEC 60601-2-18** Medical electrical equipment - Part 2-18:  
Particular requirements for basic safety and essential performance of endoscopic equipment

Rigid endoscopes are compliant with the requirements of the following applicable EU Directive:

- Medical Devices Directive **93/42/EEC**

### 1.4.2 Medical Device Classification

Sinuscopes Sinuscopes are Class IIa medical devices within the meaning of Directive 93/42/EEC.

## 1.5 Marking

### 1.5.1 Pictograms and Information on the Device and Packaging

This section describes the different pictograms that feature on the product or packaging.

|   |   |
|---|---|
|    | The device is suitable for autoclave sterilization  |
|    | Item number   |
|    | Serial number   |
|    | Lot number  |
|   | Quantity of products of a particular item number  |
|  | Direction of view   |
|  | Manufacturer  |
|  | Date of manufacture   |
|  | Caution (IEC 60601-1 3rd edition) / Take note of accompanying documents (IEC 60601-1 2nd edition) |
|  | Non-sterile, sterilize prior to each use  |
|  | MR unsafe   |
|  | Permissible storage and transport temperature   |
|  | Permissible relative air humidity during storage and transport                                    |
|  | Permissible atmospheric pressure during storage and transport                                     |
|  | Caution: Federal (US) law restricts this device to sale by or on order of a physician             |

### 1.5.2 Pictograms in this Document

This section describes the pictograms used in this document.

General warning sign

Dangerous electrical voltage warning sign

Biohazard warning sign, risk of infection

Hot surfaces warning sign

Laser beam warning sign

Type BF applied part according to IEC 60601-1

Type CF applied part according to IEC 60601-1

Defibrillation-proof type BF applied part

Defibrillation-proof type CF applied part

### 1.5.3 Direction of View Color Coding

A color-coded ring situated on the light guide connector indicates the direction of view of the endoscope.

| Color  | Direction of view |
|--------|-------------------|
| Green  | 0°                |
| Red    | 30°               |
| Black  | 45°               |
| Yellow | 70°               |

Table 1-1: Color-coded rings indicating direction of view.

## 1.6 Service Department Contact Details

If you have any questions about our products, their installation, or use, or you wish to arrange servicing, please contact one of our subsidiaries.

You will find contact details on the back of this document.

## 2 General Safety Information

### 2.1 Safety Messages in this Document

#### 2.1.1 Safety Messages at the Start of a Chapter

The safety messages described in this section will be listed at the start of any chapter containing instructions that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Take time to read these safety messages carefully and bear them in mind when performing the activities concerned.

The following kind of message refers to a risk that could lead to death or serious injury:

#### **WARNING**

##### **Nature and origin of the injury risk**

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to minor or moderate injury:

#### **CAUTION**

##### **Nature and origin of the injury risk**

Potential consequences of non-observance

> Preventive action

The following kind of message refers to a risk that could lead to material damage:

#### **NOTICE**

##### **Nature and origin of the risk of material damage**

Potential consequences of non-observance

> Preventive action

#### 2.1.2 Safety Messages in the Body of the Text

The warnings described in this section will be listed within the body of the instructions directly prior to any steps that carry a particular risk.

The severity of the potential risk is expressed by the signal word in the beginning of the message.

Read through these safety messages carefully and take the preventive action indicated.

The following kind of message refers to a risk that could lead to death or serious injury:

**WARNING!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to minor or moderate injury:

**CAUTION!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

The following kind of message refers to a risk that could lead to material damage:

**NOTICE!** Nature and origin of the risk. Potential consequences of non-observance. Preventive action.

## 2.2 Product Safety

### 2.2.1 Basic Safety Notices

Our products are developed and manufactured to the highest quality standards.

Risks despite highest of quality standards

Although this product corresponds to the current state of the art, risks could still arise during initial operation, use, reprocessing or maintenance.

Therefore, it is important that you read through these instructions for use carefully. Observe the warnings indicated.

The device must be operated only in a fault-free condition in accordance with the intended use and the instructions for use. Before each use, check that the device and accessories to be used are free of damage and in full working order.

The original packaging should be retained and reused when returning the device for servicing.

Follow the instructions for use for all devices and instruments that are to be used in conjunction with this device.

**WARNING!** This device is supplied in a non-sterile condition. Risk of infection. Carry out reprocessing of the device and its accessories prior to first use.

**WARNING!** Unauthorized modifications to the device. Risk of serious injury to persons. Do not make any unauthorized modifications.

**WARNING!** Component failure during a surgical procedure. Risk to the patient. Keep spares on hand for emergency replacement.

**WARNING!** High-intensity light source. Risk of injury to eyes. Do not look directly into the open end of the light guide or light exit of the endoscope.

**WARNING!** Magnetic resonance imaging (MRI). Magnetic force, electromagnetic interactions, heating of metal components. Do not operate the device in the vicinity of MRI scanners.

**WARNING!** Using the device during operation of a defibrillator. Risk to persons. Before discharging a defibrillator, remove device from the surgical field.

**CAUTION!** Rough handling. Risk to the patient resulting from damaged device. Handle the product with care. If the device is dropped or subjected to high mechanical stress, stop using it and send it in to the manufacturer for inspection.

**WARNING!** Improper handling, maintenance, and use carries a risk to the patient and user, or can lead to premature wearing of the device.

### 2.2.2 Staff Qualifications

Observe qualification requirements

Specific qualifications are required for operation and reprocessing of the device. The qualifications required for the personnel are stated in the respective chapters of this document and must be observed.

## 3 Operation

### 3.1 Safety Instructions

#### **WARNING**

##### **Interactions between devices in simultaneous use (e.g. lasers, electrosurgery)**

Risk to the patient and user, image interference, damage to the device

- > Ensure that all the devices in use meet at minimum the requirements classified as type BF according to IEC 60601-1
- > Observe the labeling and instructions for use of the devices used
- > Avoid direct contact between the endoscope and conductive parts with active RF electrodes
- > Do not activate RF electrodes in the vicinity of flammable gases or liquids
- > Before using RF devices, ensure that all potentially explosive gas mixtures and liquids have been extracted

#### **WARNING**

##### **Displaying a recorded image instead of the live image or an altered image orientation**

Risk to the patient

- > Ensure that the monitor is displaying the live image from the endoscopic camera
- > Ensure that the live image is displayed in the correct orientation (not mirrored)

#### **WARNING**

##### **Use of non-sterile parts**

Infection risk for the patient

- > Use only properly reprocessed endoscopes and endoscopic accessories
- > Always carry out a visual inspection prior to use
- > Anchor the light guide in the surgical field and secure against slipping, allowing some slack for maneuvering

#### **WARNING**

##### **Recontamination resulting from improper handling**

Infection risk for the patient

- > Adhere to the infection control practices that are in place

**⚠ CAUTION**

**High temperatures resulting from light source usage**

Irreversible tissue damage to patient or undesirable coagulation, user injury, material damage

- > Use only light guides that are suitable for use with the endoscope
- > Avoid the prolonged use of intense light
- > You should always use the lowest possible light output setting that will allow you to illuminate the target area
- > Do not touch the light source in the vicinity of the lamp
- > Do not touch light guide connectors
- > Do not allow the distal end of the endoscope to come into contact with patient tissue or with combustible or heat-sensitive materials
- > Do not exceed the maximum permitted ambient temperature for the light source

**⚠ CAUTION**

**Coupling laser beams in the endoscope**

Damage to the eyes caused by looking directly into the eyepiece

- > Wear laser safety glasses

**NOTICE**

**Damage resulting from improper handling**

Risk of damage to the product

- > Set down carefully without bending or knocking
- > Avoid scratching the surface
- > Take care when removing instruments from the surgical field

### 3.2 Staff Qualifications

The device may be operated only by physicians and medical assistants who have received instruction in how to use the device and who have the requisite training or advanced training, knowledge, and practical experience in the endoscopic procedure to be applied as defined by the provisions in place at the site of operation.

### 3.3 Technical Inspection Prior to Use

#### 3.3.1 Visual Inspection

##### Performing a visual inspection

*A visual inspection must be performed **prior to each use**.*

*Proceed as follows:*

- |   |   |
|---|---|
| Is there any external damage?             | 1. Check the endoscope and all the components to be used for external damage.<br><br>Do not use the product if it has any sharp edges or corners, bulges, or rough surfaces that might cause injury to the patient.   |
| Are the fiber optics in working order?    | 2. Check that the fiber optics in the endoscope are in full working order.<br><br>To do this, hold the distal end toward a bright light (not a cold light source) and the light guide connection toward your eyes. Gently move the endoscope from side to side, taking note of the brightness of the fibers.<br><br>If more than 20 % of the fibers remain dark, the device must not be used. Submit the endoscope to repair. |
| Are the surfaces clean and smooth?        | 3. Check that the glass surfaces and fiber optic end faces of the endoscope are clean and smooth.<br><br>Do not use the endoscope if the surfaces are soiled or scratched. Set aside the endoscope and follow the steps described in <u>section 5.1</u> .   |
| Is the image clear, bright, and complete? | 4. Make sure you can see a clear, bright, complete image.<br><br>Look through the eyepiece and assess the image quality.<br><br>Do not use the endoscope if the image is yellowed, dark, speckled, or cropped. Set aside the endoscope and follow the steps described in <u>section 5.1</u> .   |
| Is the OR equipment mutually compatible?  | 5. Check that all the equipment to be used for the procedure is mutually compatible.<br><br>Only use the OR equipment if the isolation barrier for protecting against electric shock meets <b>at minimum</b> the type <b>BF</b> classification requirements according to IEC 60601-1 (see <u>Intended Use</u> ).  |
|   | <input checked="" type="checkbox"/> This concludes the inspection.  |

### 3.3.2 Function Check

#### Performing a function check

*A function check must be performed **prior to each medical procedure.***

*Proceed as follows:*

1. Connect a camera.
  2. Switch on the system  
Switch on all the system components that you intend using for the procedure.
  3. Focus the image  
Direct the camera head at an object in the room and focus the image.  
Check that you are able to achieve a sharp, bright, good quality image.  
Do not use the system if the image contains vertical lines or color variations, flickers, or you are unable to achieve a sharp, bright, good quality image.
- This concludes the function check.

### 3.4 Connection of a Light Guide

#### Connecting the light guide

*Different light guides require different adapters to connect to the endoscope.*

*Proceed as follows:*

1. Screw the suitable adapter to the light guide connector on the endoscope and, if applicable, to the light guide.
2. Connect the light guide to the light guide connector.

The light guide is connected.

### 3.5 Post-Use and Pre-Cleaning Activities

#### Stopping usage and completing pre-cleaning

*Complete the pre-cleaning process **immediately after use.***

*Proceed as follows:*

1. Remove and discard the sterile drape, if present.
2. Uncouple the endoscope from the endocoupler.
3. Detach the light guide and all detachable parts.

**NOTICE!** The eyepiece cap is non-detachable.

4. Remove visible surgical residue as completely as possible using a moist, lint-free wipe.
5. Dry the product using a soft, lint-free wipe.
6. Place the product into a container for dry removal and seal.
7. Arrange for reprocessing, ensuring that the product is reprocessed within 6 hours.

This concludes use and pre-cleaning.

## 4 Reprocessing

### 4.1 Safety Notices

#### **WARNING**

**If the product is potentially contaminated with Creutzfeldt-Jakob pathogens, it is not suitable for reprocessing**

Risk of cross contamination during use and reprocessing

- > Do not reprocess products that are potentially contaminated
- > Products that are potentially contaminated must be disposed of

#### **CAUTION**

**Improper cleaning and disinfection**

Risk to patient as a result of inadequate cleaning and disinfection, damage to the device

- > Use an automatic washer/disinfector that meets the performance requirements of ISO 15883-1
- > Load the washer/disinfector with care: all contents must be rinsed and cleaned fully (no rinsing 'blind spots')
- > Connect instruments with lumens and channels directly to the dedicated rinsing ports on the device
- > Disassemble all valves on instruments
- > Ensure that the washer/disinfector is properly maintained
- > Only use cleaning and disinfecting agents that are approved for the product
- > The devices must be pre-cleaned immediately after use and reprocessed within 6 hours
- > When pre-cleaning, do not use fixing temperatures in excess of 45 °C
- > When pre-cleaning, do not use fixing cleaning and disinfecting agents (active ingredients: Aldehyde, alcohol)

#### **NOTICE**

**Contact with chloride-based solutions**

Corrosion and destruction of the device

- > Avoid contact with chloride-based solutions, which can be found, for example, in surgical residue, lotions, drugs, saline solutions, and cleaning and disinfecting agents
- > If the products should come into contact with chloride-based solutions, rinse thoroughly with DI water and dry completely

#### **NOTICE**

**Use of ultrasonic bath**

Damage to the endoscope

- > Do not expose endoscope to ultrasonic bath

**NOTICE****Unsuitable cleaning agents and process chemicals**

Corrosion damage, premature aging, and visible material deterioration

- > Only use cleaning and disinfecting agents that are approved for the product
- > Only use process chemicals that are compatible with the device's materials according to the chemical manufacturer's recommendations
- > Follow all the chemical manufacturer's application specifications regarding temperature, concentration, and contact time
- > Do not use process chemicals that can cause stress cracking or brittleness of plastic materials

**4.2 Staff Qualifications**

In many countries, the qualifications required by personnel responsible for the reprocessing of medical devices are regulated by law.

In any event, the reprocessing of medical devices must always fall under the responsibility of qualified personnel who have the necessary knowledge and expertise.

This knowledge and expertise can be acquired by completing further training in this field or as a result of completing a dedicated qualification with practical experience, supplemented by appropriate further training measures, where required.

**4.3 Validated Procedures**

The effectiveness of the following procedures:

- Manual cleaning and disinfection
- Automated cleaning and thermal disinfection
- Sterilization

as described in this document, has been fully validated.

**Operator responsibility** It is the operator's responsibility to introduce, document, implement, and maintain a validated reprocessing procedure.

Make sure that the equipment used for reprocessing is properly maintained.

**4.4 The Process Flow**

**Pre-cleaning** The product has to be prepared for reprocessing immediately after use by pre-cleaning.

Refer to [section 3.5](#) for pre-cleaning instructions.

**Work steps** The reprocessing procedure described in this document consists of the following steps:

- Pre-cleaning immediately after use
- Cleaning and disinfection (manual or automated) and
- Sterilization

**Protective clothing** During reprocessing, staff should wear personal protective clothing.

**Careful cleaning** The product must be cleaned thoroughly at the beginning of reprocessing. It is essential that the sterilization medium reaches all parts of the product.

|  |   |
|--|---|
| Automated cleaning and steam sterilization | The optimum and most reliable reprocessing results are achieved using automated cleaning and disinfection with subsequent steam sterilization using the pre-vacuum steam procedure.                       |
| Legislation and standards                  | Adhere to national statutory regulations, national and international standards and directives, and the infection control practices that are in place at your institution for reprocessing.                |
| Reprocessing prior to each use             | New and returned products following repair have to undergo complete reprocessing before first use.<br><br>The alternating use of different reprocessing methods can cause premature aging of the product. |
| www.a-k-i.org                              | Please see <a href="http://www.a-k-i.org">www.a-k-i.org</a> for more detailed information about hygienically safe, value-preserving reprocessing methods that also protect the materials.                 |

#### 4.5 Cleaning and Disinfecting Agents

The **manual cleaning** process was validated using:

- Cidezyme/Enzol (Johnson & Johnson)

The **manual disinfection** process was validated using:

- Cidex OPA (Johnson & Johnson)

The **automated cleaning** process was validated using:

- neodisher MediClean forte 0.5 %  
(Chem. Fabrik Dr. Weigert GmbH & Co. KG)

Wherever possible, use the above cleaning and disinfecting agents only.

Material compatibility releases also exist for

- Gigasept FF (Schülke & Mayr GmbH)  
for manual cleaning and disinfection
- Lysetol FF (Schülke & Mayr GmbH)  
for manual cleaning and disinfection
- Helipur HplusN (B. Braun Medical AG)  
for manual disinfection

Follow the manufacturer's instructions! Before use, read the manufacturer's user information carefully and follow the specifications regarding concentration, temperature, usage time, water quantities, and contact time.

## 4.6 Manual Cleaning and Disinfection

### 4.6.1 Overview

| Stage | Work step            | Temperature (°C/°F)     | Time (min) | Water quality  | Cleaning solution                 |
|-------|----------------------|-------------------------|------------|----------------|-----------------------------------|
| I     | Cleaning             | 35-45 / 95-113          | 5          | Drinking water | Enzymatic cleaning solution       |
| II    | Intermediate rinsing | Room temperature (cold) | 3 x 1      | Drinking water | ---                               |
| III   | Disinfection         | 20-25 / 68-77           | 12         | ---            | Disinfectant solution (undiluted) |
| IV    | Final rinsing        | Room temperature (cold) | 3 x 2      | DI water*      | ---                               |
| V     | Drying               | Room temperature        | ---        | ---            | ---                               |

Table 4-1: Overview of manual cleaning and disinfection.

\*DI water = Deionized water (demineralized, low-germ, max. 10 germs/ml, as well as low endotoxin contaminated, max. 0.25 endotoxin units/ml)

### 4.6.2 Performance

#### Complete stage I: Cleaning

Proceed as follows:

**NOTICE!** Damage resulting from rough handling. Handle the product with care. Do not knock or bend. Set down carefully.

**NOTICE!** Scratch-sensitive surface. Risk of corrosion. Do not use metal brushes, metal objects, or abrasive cleaners. To remove soiling on optical surfaces, use pads soaked with cleaning solution only.

1. Completely immerse the fully disassembled device in the cleaning solution.

Make sure that the cleaning solution reaches all accessible surfaces. When placing hollow components into the bath, ensure that all air escapes completely from the lumens.

2. While the product soaks in the solution, use a soft wipe or suitable brush to remove all visible residue from its exterior surfaces.

Brush all surfaces of products with hidden crevices, lumens with working channels, or complex geometries, for at least **one** minute or until no more residue can be removed.

3. These surfaces should then be rinsed at least **five** times thoroughly with the cleaning solution. Use a disposable 20 ml syringe.

**Complete stage II: Intermediate rinsing***Proceed as follows:*

1. Fully rinse all accessible surfaces three times for at least **one** minute each.

Use fresh water for each rinse.

2. Allow residual water to drip off for a sufficient length of time.

**Complete stage III: Disinfection***Proceed as follows:*

1. Completely immerse the product in the disinfectant solution bath.

Make sure that the disinfectant solution reaches all accessible surfaces.

Remove all adherent air bubbles from the component surfaces. When placing hollow components into the bath, ensure that all air escapes completely from the lumens.

2. Thoroughly rinse all surfaces of products with hidden crevices, lumens with working channels, or complex geometries at least five times. Use a disposable 20 ml syringe.

**Complete stage IV: Final rinsing***Proceed as follows:*

1. Rinse the device in three rinses of two minutes each. Use fresh water for each rinse.
2. Thoroughly rinse all surfaces of products with hidden crevices, lumens with working channels, or complex geometries at least **five** times.

Use fresh water for each rinse.

3. Allow residual water to drip off for a sufficient length of time.

**Complete stage V: Drying***Proceed as follows:*

1. Dry the device using a soft, lint-free wipe.

Use medical-quality filtered compressed air ( $p_{max} = 0.5$  bar) to dry the areas you cannot reach with the wipe.

- Visual inspection
2. After drying, visually inspect the devices in a well-lit area (see [section 3.3.1](#)) and make sure they are undamaged and free of visible residue.

Use a magnifying glass. Repeat the cleaning process if necessary.

Immediately set aside any damaged products.

- This concludes the manual cleaning and disinfection process.

## 4.7 Automated Cleaning and Thermal Disinfection

### 4.7.1 Overview

| Stage | Work step            | Temperature (°C/°F) | Time (min) | Water quality  | Cleaning solution / notes  |
|-------|----------------------|---------------------|------------|----------------|--|
| I     | Pre-rinsing          | < 25/77             | 2          | Drinking water | ---  |
| II    | Cleaning             | 55/131              | 10         | Drinking water | Alkaline detergent, e.g. 0.5 % neodisher MediClean forte (5 ml/L), ph > 10 |
| III   | Rinse I              | > 10/50             | 1          | Drinking water | ---  |
| IV    | Rinse II             | > 10/50             | 1          | DI water*      | ---  |
| V     | Thermal disinfection | > 90/194            | 5          | DI water*      | A <sub>0</sub> value > 3000  |
| VI    | Drying               | ---                 | ---        | ---            | ---  |

Table 4-2: Overview of automated cleaning and disinfection.

\*DI water = Deionized water (demineralized, low-germ, max. 10 germs/ml, as well as low endotoxin contaminated, max. 0.25 endotoxin units/ml)

**Note:** Automated reprocessing can cause color-anodized or plastic components (e.g. serial rings and the eyepiece cap) to fade.

## 4.7 Automated Cleaning and Thermal Disinfection

## 4.7.2 Performance

**Performing automated cleaning and thermal disinfection****AUTOCLAVE**

*Thermal disinfection should be performed only for products that are labeled as autoclavable.*

*Use a washer/disinfector that has been validated as effective and meets the performance requirements of ISO 15883-1 (or the respective country-specific version thereof).*

*Dry removal of the products is recommended prior to automated cleaning.*

*In the case of wet removal, do not use foaming detergents and rinse the products thoroughly prior to the automated cleaning.*

*Use deionized water for thermal disinfection.*

*Proceed as follows:*

1. Load the washer/disinfector as per the validated loading plan. Follow the manufacturer's instructions and instructions for use for the device.

When loading, take care to avoid creating rinsing blind spots.

2. Start the cleaning/disinfection cycle in accordance with the manufacturer's instructions and instructions for use for the device.

**CAUTION!** Risk of scalding when unloading the washer/disinfector. Be sure to wear gloves.

3. Take the reprocessing tray and the endoscopes out of the washer/disinfector.

- Visual inspection
4. After drying, visually inspect the products in a well-lit area and make sure they are dry, undamaged and free of visible residue.

Use a magnifying glass. Repeat the cleaning and disinfection process if necessary.

Immediately set aside any damaged products.

- This concludes the automated cleaning and disinfection process.

**4.8 Sterilization**

**4.8.1 Steam Sterilization**

**Performing steam sterilization**

**AUTOCLAVE**

*Steam sterilization should be performed only for products that are labeled as autoclavable.*

*Use demineralized feed water that meets the requirements of European Standard EN 285.*

*Proceed as follows:*

1. Make sure that manual or automated cleaning and disinfection is complete.
2. Place all parts in a reprocessing tray.
3. Double wrap the tray with two single layers of polypropylene wrap.

**CAUTION!** Risk of scalding when loading the sterilizer. Wear suitable gloves.

4. Load the sterilizer as per the validated loading plan.

Follow the manufacturer's instructions!

5. Start the sterilization cycle in accordance with the manufacturer's instructions and instruction for use for the sterilizer.

Use the following validated parameters for steam sterilization:

|              |  |
|--------------|--|
| Cycle type   | Fractionated vacuum processes                    |
| Pulses       | 4  |
| Temperature  | 134 °C (273 °F)                                  |
| Holding time | at least 3 min<br>(effective sterilization time) |

**CAUTION!** Risk of scalding when unloading the sterilizer. Wear suitable gloves.

**NOTICE!** Products can be shock sensitive when hot. Avoid knocking or shaking.

**NOTICE!** Damage resulting from a sudden change in temperature. Allow products to cool to room temperature; do not use additional cooling measures.

6. Remove the sterilized products from the sterilizer.  
Ensure that the products remain sterile after reprocessing.

This concludes sterilization.

#### 4.8.2 STERRAD® Sterilization Systems (100S, NX, 100NX)

The STERRAD® Sterilization Systems, manufactured by Advanced Sterilization Products (ASP), use low-temperature, hydrogen peroxide gas plasma technology for terminal sterilization of properly cleaned, rinsed, and dried reusable medical devices.

Refer to the STERRAD® Sterilization Systems User Guide for detailed instructions for use of any STERRAD® unit, ASP's STERRAD® Sterility Guide (SSG) at [www.sterradsterilityguide.com](http://www.sterradsterilityguide.com) or contact ASP customer service.

##### Performing sterilization

**NOTICE!** STERRAD® sterilization may result in cosmetic device changes that do not necessarily affect the device function.

**NOTICE!** All products must be thoroughly dried before loading into the STERRAD® sterilizer. Loads containing moisture may cause a cycle cancellation.

*Proceed as follows:*

1. Place all parts in a reprocessing tray. Place a STERRAD® indicator strip in each tray.

Double wrap the tray with two single layers of polypropylene wrap.

2. Load the wrapped tray in the sterilizer.

The tray should be positioned such that the plasma can fully surround it. Do not allow any item to touch the wall of the sterilizer.

Follow the manufacturer's instructions!

3. Start the sterilization cycle in accordance with the manufacturer's instructions and instructions for use for the sterilizer.

STERRAD® sterilization was validated for the following cycles:

- STERRAD® 100S Short Cycle
- STERRAD® NX Standard Cycle
- STERRAD® 100NX Standard Cycle

4. Remove the sterilized products from the sterilizer.

Ensure that the products remain sterile after reprocessing.

- This concludes sterilization.

## 5 Maintenance and Repairs

### 5.1 Troubleshooting Table

| Issue                                   | Possible causes   | Fixes  |
|---|---|--|
| Image cloudy                            | Glass surfaces soiled   | Manual cleaning as per <a href="#">section 4.6</a> , then re-process, check the water quality      |
|   | Leaky, defective lens system  | Send the endoscope in for repair   |
| Image too dark, too little illumination | Glass surfaces soiled   | Manual cleaning as per <a href="#">section 4.6</a> , then re-process, check the water quality      |
|   | Unsuitable light guide  | Use a suitable light guide   |
|   | Light guide incorrectly attached to endoscope   | Check that the light guide is seated properly  |
|   | Defective fiber optics  | Check fiber optics as per <a href="#">section 3.3</a>  |
|   | Defective light guide or light source   | Check the light guide and light source   |
| Image is yellowed                       | Soiled fiber optics   | Manual cleaning as per <a href="#">section 4.6</a> , then re-process, check the water quality      |
|   | Light guide soiled or defective   | Check the light guide (e.g. by illuminating a white surface)                                       |
| Corrosion, staining, discoloration      | Inadequate cleaning (e.g. remaining protein residue)  | Manual cleaning as per <a href="#">section 4.6</a> , if necessary, rub thoroughly, then re-process |
|   | Insufficient rinsing between the reprocessing stages, especially prior to sterilization                                   | Rinse sufficiently between the reprocessing stages   |
|   | Excessive chloride concentration  | Check water quality  |
|   | Excessive concentration of minerals (e.g. limescale) or organic substances  | Check the water quality, only use demineralized water if necessary                                 |
|   | Heavy metal ions and/or silicates, elevated concentration of iron, copper, and manganese in the water or sterilizer steam | Check the water quality, only use demineralized water if necessary                                 |
|   | Contaminated, too frequently reused cleaning and disinfecting solution  | The cleaning and disinfecting solution should be replaced regularly                                |

| Issue | Possible causes   | Fixes  |
|-------|---|--|
|       | Extraneous rust, e.g. resulting from rust contained in steam or damaged or non-corrosion-resistant instruments reprocessed at the same time | Check the supply systems, if multiple items are being reprocessed together, pay attention to material compatibility and signs of prior damage; prevent the different items from touching one another |
|       | Contact corrosion   | Avoid contact with other devices   |

Table 5-1: Troubleshooting.

## 5.2 Repair

Should you need to arrange a repair for the device, please contact our Service Department.

When you return the device, please also enclose a detailed **description of the fault**.

### WARNING

#### Contaminated device

Risk of infection

- > Reprocess the product properly prior to shipping

Only ship devices that have been properly reprocessed.

The device should be returned in its original packaging, if possible.

Affix a label to the outer packaging indicating the reprocessing status.

We reserve the right to refuse acceptance of unlabeled packages and to return them to sender.

## 6 Product Data

### 6.1 Technical Data

|                                     |                        |
|-------------------------------------|------------------------|
| <b>Item no.</b>                     | <b>STS-7027-110ASA</b> |
| Direction of view                   | 0°                     |
| Field of view                       | Standard               |
| Working length                      | 110 mm                 |
| Outer diameter of insertion portion | 2.7 mm                 |

|                                     |                        |
|-------------------------------------|------------------------|
| <b>Item no.</b>                     | <b>STS-7027-110BSA</b> |
| Direction of view                   | 30°                    |
| Field of view                       | Standard               |
| Working length                      | 110 mm                 |
| Outer diameter of insertion portion | 2.7 mm                 |

|                                     |                        |
|-------------------------------------|------------------------|
| <b>Item no.</b>                     | <b>STS-7027-110CSA</b> |
| Direction of view                   | 70°                    |
| Field of view                       | Standard               |
| Working length                      | 110 mm                 |
| Outer diameter of insertion portion | 2.7 mm                 |

|                                     |                       |
|-------------------------------------|-----------------------|
| <b>Item no.</b>                     | <b>STS-7040-175AA</b> |
| Direction of view                   | 0°                    |
| Field of view                       | Wide angle            |
| Working length                      | 175 mm                |
| Outer diameter of insertion portion | 4 mm                  |

|                                     |                       |
|-------------------------------------|-----------------------|
| <b>Item no.</b>                     | <b>STS-7040-175BA</b> |
| Direction of view                   | 30°                   |
| Field of view                       | Wide angle            |
| Working length                      | 175 mm                |
| Outer diameter of insertion portion | 4 mm                  |

|                                     |                       |
|-------------------------------------|-----------------------|
| <b>Item no.</b>                     | <b>STS-7040-175CA</b> |
| Direction of view                   | 70°                   |
| Field of view                       | Wide angle            |
| Working length                      | 175 mm                |
| Outer diameter of insertion portion | 4 mm                  |

|                                     |                        |
|-------------------------------------|------------------------|
| <b>Item no.</b>                     | <b>STS-7040-175BSA</b> |
| Direction of view                   | 30°                    |
| Field of view                       | Standard               |
| Working length                      | 175 mm                 |
| Outer diameter of insertion portion | 4 mm                   |

|                                     |                       |
|-------------------------------------|-----------------------|
| <b>Item no.</b>                     | <b>STS-7040-175FA</b> |
| Direction of view                   | 45°                   |
| Field of view                       | Wide angle            |
| Working length                      | 175 mm                |
| Outer diameter of insertion portion | 4 mm                  |

|                                     |                        |
|-------------------------------------|------------------------|
| <b>Item no.</b>                     | <b>STS-7027-187ASA</b> |
| Direction of view                   | 0°                     |
| Field of view                       | Standard               |
| Working length                      | 187 mm                 |
| Outer diameter of insertion portion | 2.7 mm                 |

|                                     |                         |
|-------------------------------------|-------------------------|
| <b>Item no.</b>                     | <b>STS-7027-187-BSA</b> |
| Direction of view                   | 30°                     |
| Field of view                       | Standard                |
| Working length                      | 187 mm                  |
| Outer diameter of insertion portion | 2.7 mm                  |

|                                     |                         |
|-------------------------------------|-------------------------|
| <b>Item no.</b>                     | <b>STS-7027-187-CSA</b> |
| Direction of view                   | 70°                     |
| Field of view                       | Standard                |
| Working length                      | 187 mm                  |
| Outer diameter of insertion portion | 2.7 mm                  |

|                                     |                       |
|-------------------------------------|-----------------------|
| <b>Item no.</b>                     | <b>STS-7020-110AA</b> |
| Direction of view                   | 0°                    |
| Field of view                       | Wide angle            |
| Working length                      | 110 mm                |
| Outer diameter of insertion portion | 2 mm                  |

|                                     |                       |
|-------------------------------------|-----------------------|
| <b>Item no.</b>                     | <b>STS-7020-110BA</b> |
| Direction of view                   | 30°                   |
| Field of view                       | Wide angle            |
| Working length                      | 110 mm                |
| Outer diameter of insertion portion | 2 mm                  |

## 6.2 Ambient Conditions

### Transport and storage conditions

|                       |                   |
|-----------------------|-------------------|
| Temperature           | -20 °C to +70 °C  |
| Relative air humidity | 5 % to 95 %       |
| Atmospheric pressure  | 70 kPa to 106 kPa |

Reprocessed devices must be protected against recontamination, in a dark, dry, well-ventilated, dust-free and light-controlled location with a consistent temperature.

Direct sunlight, high temperatures, high humidity or radiation can damage the device or present a risk of infection.

When placing into storage, make sure that the device cannot be damaged by other instruments. It is therefore best to store the device individually or use containers in which it can be secured.

### Operating conditions

|                       |                   |
|-----------------------|-------------------|
| Temperature           | +15 °C to +28 °C  |
| Relative air humidity | 5 % to 95 %       |
| Atmospheric pressure  | 70 kPa to 106 kPa |

### 6.3 Spare Parts and Accessories

Use original spare parts and accessories only.

| Image   | Designation                              | Item number               |
|---|--|---------------------------|
|  | Light guide adapter<br>for Storz systems | STS-ACC-0000-LGUA         |
|  | Light guide adapter<br>for Wolf systems  | STS-ACC-0000-LGUA-<br>RWF |

Table 6-1: Spare parts and accessories.

## 7 Disposal

### **WARNING**

#### **Contaminated device**

Risk of infection

- > The device must be reprocessed prior to disposal

When disposing of or recycling the device and its components, you must adhere to the applicable national regulations governing waste disposal and recycling.

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