

See the difference



Bipolar Resection

Autoclavable Instrument Product Guide

Bipolar Resection Set



Recommended Set

Part Number	Description	Size
STS-2040-302BA	HD 30° Cystoscope	4mm x 302mm
STS-2040-302ESA	HD 12° Cystoscope	4mm x 302mm
STS-URO-340-20-TL-SET	Resectoscope Sheath Set includes: Outer sheath, Rotatable Inner Sheath With Ceramic Tip and Obturator	26Fr.
STS-URO-230-24*	Visual Obturator With 1 Instrument Channel for 4mm Cystoscope	24/26 Fr.
Working Element	Passive Working Element compatible with Karl Storz or Olympus	26Fr.
STS-URO-400-10 or 10QL	Sheath Connection With Locking For Ellik Evacuator; in Bayonett-Lock or Quick-Lock	-
STS-URO-940-02	Luer-Lock Tubing Connector With Stopcock	-
LIGHT CABLE	Your choice of a Clear Case Fiber Optic	10Ft.
STS-STB-H52-100P*	Sterilization Tray for Cystoscopy and Resectoscope kits	-

* Visual obturator has 5 Fr. working channel which can be used for instrument insertion or a laser bridge.

Light Cables

Part Number	Description	Size
STS-LC-KK-10CLR	Clear Case Fiber Optic Cable Storz-Storz Connection	10Ft.
STS-LC-KO-10CLR	Clear Case Fiber Optic Cable, For Olympus Light Source and Storz/Olympus Endoscope	10Ft.
STS-LC-STR-10CLR	Clear Case Fiber Optic Cable Stryker-Storz Connection	10Ft.



3020 NW 82nd Avenue
Miami, Florida USA 33122
sales@endoscopia.com | Tel: +1 (305) 468-1455
ISO 13485:2016 – Certified Quality Management System