





# Bipolar System



## Strauss Bipolar

This bipolar system delivers precise energy concentration at the tip of the electrode for the least amount of thermal spread to the tissue. This leads to a more precise cut which allows the surgeon a greater amount of control during procedures. Strauss Surgical offers the finest single stem bipolar classic style electrode in the US market.



This bipolar system provides a regulated current flow to the tip of the electrode. This gives you precise energy for the cut and coagulation. The Strauss bipolar system also alleviates unnecessary energy being in contact with patient tissue since it does not use the sheath as the return like hybrid bipolar systems.

## Generator

- The generator provides instant plasma ignition, which provides a clean consistent cut that leads to minimal to no tissue charring.
- This bipolar system does not require fluid warmers for optimal performance unlike other systems.





## Generator

Part Number	Description
MCB	Electrical Surgical unit

## Passive Working Element

Part Number	Description
VP1S	Passive working element, compatible with Karl Storz
VP3S	Passive working element, compatible with Olympus

## Bipolar Electrodes

Part Number	Description
MV1S	Thin bipolar cutting loop
MV2S	Thick bipolar cutting loop
MVVS	Vaporization bipolar electrode
MVKS	Knife bipolar electrode

**strauss**<sup>®</sup>  
S U R G I C A L

+1 (305) 436-0599  
[sales@straussurgical.com](mailto:sales@straussurgical.com)  
[www.straussurgical.com](http://www.straussurgical.com)